# PERMAPLAST LH HYBRID

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Stand: 300817

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### GB Instructions

**Description and indication:** M+W Permaplast LH hybrid is a light-curing, radiopaque composite developed for restoration of the front and side tooth regions (Classes I-V). The high total filler content in M+W Permaplast LH hybrid ensures high strength, excellent colour stability, excellent abrasion properties and high polishability. M+W Permaplast LH hybrid is based on Bis-GMA resin and inorganic filler particles between 0.07 and 1.5 µm.

**Composition:** Contents: Barium glass (silan.), pyrogenic silicic acid (silan.), silicon dioxide, Bis-GMA, catalysts and pigments. Total filler content: 83 % (77 vol. %). Total inorganic filler content: 57 % (77 vol. %).

**Contraindications:** The application of synthetic restoration materials is contraindicated:  
- if they cannot be applied under dry conditions or if it is not possible to employ the methods for use specified below,  
- for direct application to dentine in close proximity to the pulp or on exposed pulp,  
- in the case of proven allergies to components of the synthetic restoration materials.

**Side effects:** No systemic side effects have been reported. Isolated cases of type IV allergic reactions have been observed.

**Safety precautions for use:** Suitable protection for the pulp/dentine is required to avoid pulp reactions. The working area is to be dried.

Avoid contact between the material and skin and eyes. In case of accidental contact, wash thoroughly with water.

**Interactions with other substances:** Do not use in conjunction with eugenol/clove oil-containing materials, as these may interfere with the curing of the material and cause staining. Avoid water and oil-containing air, as these may prevent polymerisation in the region of contact. Contact with cationic mouthwashes, plaque indicators and chlorhexidine may cause staining.

**Use:**  
**Isolation:** Use of a cofferdam is recommended.

**Cleansing:** Before preparation, cleanse the tooth with a rubber cup and a fluoride-free and oil-free polishing paste. The residue is carefully removed using a water jet and the surface of the tooth is then dried with oil-free air.

**Cavity preparation:** Prepare the cavity in the usual way, minimising loss of tooth material and retaining occlusal regions as far as possible. The production of undercuts is not required unless the restoration is exposed to high pressure and tensile loads. To obtain durable filling edges, it is recommended that the edges be chamfered to the enamel edges, especially in the occlusal region. Rinse the cavity thoroughly with water and dry with oil-free air. The cavity should be dry and free of contaminants. Avoid overthickening the dentine. For Class II cavities, use a combination of interproximal wedges and

preformed matrix strips to obtain secure interproximal contact.

**Pulp protection:** To avoid any potential pulp reactions, the exposed dentine over the pulp chamber in deep cavities is covered with a thin layer of a suitable cavity lining material (such as spot-application of calcium hydroxide-containing cavity lining materials).

**Etching and bonding:** Please follow the application instructions of the product in question.

**Cavity lining application:** When using a cavity lining material please follow the relevant instructions for use of the product in question.

**Application:**  
**Cavifill:** place Cavifill in a dispensing gun and apply material into cavity. Please note that the Cavifill presentation form is an item for single use. A Cavifill may only be used for one patient.

**Rotating syringe:** dispense required amount onto a mixing block or glass plate. Then apply into cavity using a suitable instrument. Please close syringe immediately after use.

Apply composite in cavity with layers 1.5 to 2 mm thick and cure using a suitable light polymerisation device. Cure the first layer from buccal and lingual/palatal, or via a light prism with a halogen lamp. When using metal matrices, the layers should have a maximum thickness of 1 mm as curing can only be performed from occlusal. Repeat the appli-

cation technique until the cavity is filled. Polymerise each of the individual layers for 30 to 40 seconds. Do not loosen a matrix, if used, during the application and curing steps.

The light-curing material is not suited for curing with a UV light.

**Shaping of cured filling:** The finishing process is carried out with a diamond or carbide finishing bur or finishing strips. High-shine surfaces can be achieved by a final finishing procedure using fine strips or silicone polish.

**Classification:** Meets the requirements of ISO 4049:1988.

**Storage:** The paste material has a shelf life of 3 years when stored at a maximum of 25 °C.

Do not use the product after the expiry date.

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Müller & Weygandt GmbH • Reichardsweide 40 • D-63654 Büdingen • Tel.: +49 (0)6042 – 88 00 88