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
GB Instructions of use

Read special instructions for use carefully.

Ingredients:

Paste 1: Dipropylene glycol diglycidyl ether, diethylamine, diethylenetriamine, barium sulphate, hydroxylapatite, 3-Aminomethyl-3,5,5-trimethylcyclohexylamine, polyhexamethylene biguanide hydrochloride, silicon dioxide, aqua

Paste 2: Bis(4,4'-glycidylphenoxy)propane, C12-C14 glycidyl ether, barium sulphate, tri-calcium phosphate

Properties:

M+W Permaroot is a permanent root filling material based on familiar epoxide-amine chemistry. The material is self-sterilizing to avoid reinfection of the root canal. It is x-ray-opaque, dimensionally stable. The tissue will not react with the material. Easy to place the ready-to-use paste in the root canal with familiar root canal treatment techniques. The material fulfils the requirements of EN ISO 6876. Strong adhesion to the dentin. Recommended for all orthograde filling techniques.

Application:

After preparing the root canal using standard methods, and if necessary using a temporary filling material the root canal is filled with M+W Permaroot and guttapercha points. Before use, the cap of the double-chamber syringe must be removed by turning counter-clockwise. A mixing cannula supplied with the product is attached by inserting the spring of the mixing cannula into the groove in the syringe. The cannula is tightened by turning in a clockwise direction. When you press the piston, both pastes are mixed in the mixing cannula in a ratio of 1:1 to form a ready-to-use paste. The paste can be applied and worked from a mixing block using a lentulo or guttapercha point. Alternatively, the material can be applied directly to the canal openings and coronal part of the canal using Endo tips. After treatment is completed, the cap must be screwed back on to the syringe the mixing cannula discarded.

Important notes!

After treatment is completed, the cap must be screwed back on to the syringe and the mixing cannula discarded. Otherwise premixed material in the cannula may flow back into the syringe and cure the epoxy component thus to destroy the contents in the syringe. Apply only into dry root canals. Do not use if the syringe is damaged.

Cleaning of instruments:

Clean all instruments and mixing slabs immediately after use with a 50% dilution of ethanol in water and traces of detergent.

Working time / Setting time:

The mixed root filling paste is workable for 15 minutes approx. Thermic techniques will shorten the workable time. Full setting is reached at 37 °C body temperature within 24 hours.


Storage:

Store M+W Permaroot at room temperature and in a dry place (do not refrigerate). Do not expose to direct sunlight. Seal properly after use. Once opened, use within 6 months. Do not use after the expiry date. Keep out of the reach of children. Improper storage will reduce the shelf life and may impair the product's effectiveness.

Warning:

Do not use the material on patients who are allergic to epoxy resins, amines or other components of the product. Do not allow unpolymerised pastes to come into contact with the skin or eyes. If eye contact occurs, rinse thoroughly with water immediately. Avoid contact with mucous membranes. If the root filling material comes into contact with mucous membranes remove it with suitable materials (e.g. cotton pellets). Rinse out the oral mucosa with plenty of water. In the event of continuing sensitization, seek medical advice. If accidentally swallowed, drink plenty of water. In the event of nausea or feeling unwell, seek medical advice.

Indicate batch number (LOT) to identify product in case of any occurrence. Dispose product according to official regulations.

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