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## CB Instruction sheet

One for everything! Our universal phosphate-bonded speed investment covers the entire spectrum in the dental field. The perfectly coordinated universal liquid allows the expansion of the various applications to be controlled in a targeted manner.

One investment material for all applications, only one liquid, highly precise, can be used conventionally or in the speed method.

The corresponding expansion is controlled by means of the expansion table. (see table)

Mixing ratio:	100 g powder: 24 ml liquid 150 g powder: 36 ml liquid
Processing time:	approx. 5 min.
Processing temperature:	19 °C

### Processing:

First pour the required amount of liquid into the beaker. Then add the investment material. Mix by hand for 20 – 30 sec, then stir in the vacuum mixer for 60 sec.

Setting time:	approx. 20 min.
Furnace temperature: :	max. 850 °C for rapid casting
Compressive strength:	6 – 7 MPa

### Preheating / burnout

#### Conventional (with holding stages)

- Fill the investment ring with investment material and allow to harden until it has cooled down to room temperature.
- Roughen the surface of the investment ring and place it in the cold furnace.
- Leave the investment ring at 250°C for 20 – 30 min., depending on size.
- Heat oven evenly to final temperature.
- Allow the investment ring to stand for another 30 – 60 min. (for final temperature, see recommendation of the of the alloy or pressed ceramic manufacturer).

#### Quick heating

- Fill casting ring and allow to harden for at least 20 min.
- Roughen the surface of the investment ring and place it in the furnace heated to final temperature.
- Allow muffle to stand at final temperature for 25 – 70 min. Please observe the specifications of the alloy or pressed ceramic manufacturer. (25 min. 1-cavity muffle; 30 min. 3-cavity muffle; 50 min. 6-cavity muffle; 70 min. 9-cavity muffle)

Always store the investment material in a cool, dry place.

Please use only the M+W SELECT CHROMAVEST LIQUID.

The investment is suitable for metal muffle rings as well as ring-free systems.

We recommend the use of specially adapted products in combination with the investment material.

### Warranty

All recommendations regarding the application, are based on our own experience / tests and can therefore only be seen as guide values. Our products are subject to continuous further development. We therefore reserve the right to make changes in design and composition. Our information and recommendations are based on the currently known state of science / technology and are to be regarded as correct according to our knowledge and experience at the present time. The above version supersedes all previous information.

This investment material contains quartz and cristobalite.  
Inhalation of dust must therefore be avoided!

### M+W select Chromavest Liquid Press ceramic

Inlays / Partial crowns	45 – 55%
Crowns / Bridges	55 – 65%

### M+W select Chromavest Liquid

#### Casting technique

##### EM

Inlays / Partial crowns	35 – 45%
Crowns / bridges / primary parts	45 – 55%
Secondary parts	50 – 60%

##### NEM

Crowns / bridges / primary parts	70 – 80%
Abutments	80 – 90%

### Recommendations for working with 3D printed objects

Object	EM	NEM
Crown 3D	30 – 40%	50 – 60%
Bridge 3D	30 – 40%	50 – 60%

### Conversion table for 100 g powder

Quantity in ml (total = 24 ml)

Quantity in %	Liquid	H <sub>2</sub> O
30	7	17
35	8	16
40	10	14
45	11	13
50	12	12
55	13	11
60	14	10
65	15	8
70	16	7
75	17	6
80	18	5
85	20	4
90	22	2

**ATTENTION!**  
3D objects must be cast with steel ring  
and be cast with a double fleece insert.