# 3M ESPE

## CoJet® Sand

- en **Blast-Coating Agent 30 µm**
- de **Beschichtungsstrahlmittel 30 µm**
- fr **Produit pour poncer les couches 30 µm**
- it **Materiale di rivestimento a 30 µm**
- es **Material arenado para recubrimiento 30 µm**
- pt **Produto de jorro especial para revestimento de 30 µm**
- nl **Straalmiddel voor het aanbrengen van een coating 30 µm**
- el **Υλικό αμολογίας για επίστρωση 30 µm**
- de **Blästringsmaterial für coating 30 µm**
- fi **Pinnovaitea erikoisohiekkaupahuus materiaali 30 µm**
- da **Specielt strålemiddel til overfladecoating 30 µm**
- no **Materiale for pålegging av sjikt ved sandblåsing 30 µm**

Instructions for Use  
Gebrauchsinformation  
Mode d'emploi  
Informazioni per l'uso  
Informação de uso  
Instrucciones de uso

Gebrauchsinformation  
Οδηγίες χρήσης  
Bruksanvisning  
Käyttöohje  
Bruksanvisning  
Bruksanvisning

3M Deutschland GmbH  
Dental Products  
Carl-Schurz-Str. 1  
41453 Neuss - Germany

3M ESPE  
Dental Products  
St. Paul, MN 55144-1000 U.S.A.  
3M ESPE Customer Care/MSDS Information  
U.S.A.: 1-800-834-2249 and Canada 1-866-363-3685.  
3M, ESPE, CoJet and Sinfony are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada.  
© 2012, 3M. All rights reserved.

## ENGLISH

### Product Description

CoJet Sand is a blast-coating agent for silicatization (ceramic-like coating) of metal, ceramic and composite materials. CoJet Sand is applied with micro-blasting units, e.g. CoJet® Prep. Due to its small grain size of 30 µm, CoJet Sand is particularly well suited for gentle intraoral applications. Silicatization followed by silanization with ESPE® Sil Silane solution are essential steps to form a chemical bond to composite materials. Materials treated with CoJet Sand and ESPE Sil show substantially stronger bonding as compared to materials processed with rotary instruments. CoJet Sand and ESPE Sil are components of the CoJet® System.

\* For details on all mentioned products please refer to the corresponding Instructions for Use.

\*\* Instructions for Use should not be discarded for the duration of product use.

### Indications

- Silicatization as a preparatory step in the repair of ceramic or composite veneers with or without exposed metal, all-ceramic or all-composite restorations or composite fillings
- Silicatization of ceramic, composite or metal restorations as a preparatory step of adhesive cementation, e.g. of inlays, onlays, facings, etc.
- Mechanical increase of the retentive surface of metal restorations prior to conventional cementation: e.g. crowns, bridges, inlays, onlays, endodontic pins, pin bulk-ups, etc.

### Precautionary Measures

3M MSDS can be obtained from [www.3m.com](http://www.3m.com) or contact your local subsidiary.

### Preparation

- In intraoral applications place a rubber dam.
- Use e.g. a matrix to cover dental areas that do not to be coated.
- Dental personnel and patient must wear suitable protective goggles; the dental personnel must also wear respiratory protection.
- Fill the CoJet Sand into the micro-blasting unit.
- Set the air pressure to 2-3 bar (30-45 psi).
  - If the air pressure is less than 2 bar, the level of impact energy will be too low and no coating with silicates will occur. If the air pressure exceeds 3 bar, the workpiece will be coated with silicates, but there will also be an unnecessary abrasion.
- Test the operation of the blasting unit on a work sample (see Application).
- Gently smooth the surface to be repaired using rotary instruments to render the surface evenly accessible to blasting. Then clean with a water spray and dry with air.

### Processing

- Blast the work sample from a distance of 2-10 mm and perpendicular to the surface of the work sample while providing for aspiration.
- Objects of the size of a veneering unit should be evenly coated: for 15 sec, where smaller or larger surfaces should be treated for proportionally shorter or longer times, respectively.
- Metal surfaces become evenly dark, provided the pressure is appropriate and the entire surface is coated.
  - Working with ceramic or composite parts, the success of coating cannot be ascertained by inspecting the color of the dental work. Thus, it is essential to strictly comply with the specified coating times!
- Remove any residual blast-coating agent with a stream of oil-free air.
- Avoid contamination by fluids, e.g. saliva, to prevent detrimental effects on coating.

### Postprocessing

- Directly after coating, dose ESPE Sil in a Dappen dish and use a clean brush to wet the entire surface with the solution.
- The use of ESPE Sil is not required in CoJet Sand applications, which only aim at roughening the surface of metal or metal-ceramic restorations for subsequent conventional cementation.
- Allow the silane solution to dry for 30 sec in intraoral applications or for 5 min if applied extraorally.
- Avoid contamination by fluids, e.g. saliva, to prevent detrimental effects on coating.
- In repairs, coat any exposed metal parts with Sinfony® Opaqueur.
- Use a disposable brush to apply a thin layer of Visio Bond to ceramic or composite materials, and then cure with light. In cementations, do not use any bonding agent.
- Proceed in accordance with the filling or fixation material's Instructions for Use.

### Storage and Shelf-life

Store the product in a dry place. Do not use after the expiration date.

### Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

### Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

### Limitation of Liability

Where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty contract, negligence or strict liability. Date of the information: January 2012

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