# 3M | CapMix™

ESPE

- (en)** Capsule Mixing Device
- (de)** Kapselmischgerät
- (fr)** Appareil à malaxer les capsules
- (it)** Miscelatore di capsule
- (es)** Aparato mezclador de cápsulas
- (pt)** Misturador para cápsulas
- (nl)** Mengapparaat voor capsules
- (el)** Συσκευή ανάμιξης για κάψουλες
- (sv)** Kapselblandare
- (fi)** Kapselinsekoitin
- (da)** Blandeapparat til dentale materialer i kapsler
- (no)** Blandeapparat for kapsler

Operating Instructions  
Betriebsanleitung  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones de servicio  
Instruções de serviço

Bedieningshandleiding  
Οδηγίες χρήσεως  
Driftanvisning  
Käyttöohjeet  
Brugsanvisning  
Driftsveiledning

**en** ENGLISH

# CapMix™

## Capsule Mixing Device

Table of Contents	Page
1. Safety	1
2. Product Description	2
3. Fields of Application	2
4. Technical Data	2
5. Installation	2
5.1 Setting the Mixing Time	3
5.2 Changing the Basic Settings	3
6. Operation	3
6.1 Insertion of the Capsule	3
6.2 Mixing	3
6.3 Removal of the Capsule	3
7. Notes	3
8. Operational Errors	3
9. Maintenance and Care	3
9.1 Fuse Replacement	3
9.2 Care	4
10. Return of Old Electric and Electronic Equipment for Disposal	4
10.1 Collection	4
10.2 Return and Collection Systems	4
11. Customer Information	4
11.1 Warranty	4
11.2 Limitation of Liability	4
Symbol Glossary	5

## 1. Safety

### CAUTION!

Please read these instructions carefully before you set up and start the device. As with all technical equipment, this unit will function properly and be safe to operate only if both the general safety regulations and the specific safety instructions laid out in these Operating Instructions are followed.

- The unit must be operated in accordance with the following instructions. The manufacturer accepts no liability for any damage arising from use of the unit for any other purpose.
- Prior to setting-up the device, ensure that the voltage stated on the rating plate is identical to the available main voltage. Use of the device at different voltages may lead to the destruction of the unit.
- The unit must be operated with the cover intact and closed, otherwise, there is a risk of injury.
- The cover must not be opened until the capsule holder has come to a complete stop.
- Only the capsule to be mixed but no other objects may be present in the space below the cover when the device is operated.
- Placement or operation of the unit on non-level surfaces is not permissible.
- When transferring the unit from a cold to a warmer ambient temperature, there may be a hazard due to condensation of water on the device. Hence, the unit should only be operated after fully equilibration to ambient temperature.
- For operation of the unit, connect the main plug to an adequate electrical socket with grounded lead. If an extension cable is used, ensure that the grounded lead is not interrupted.
- Before replacing fuses, the device must be turned off and the main plug must be disconnected from the socket.
- Only the appropriate fuses may be used for fuse replacement. The use of recycled fuses or short-circuiting of the fuse holder is prohibited.
- In order to prevent electrical shock hazards do not insert any objects into the unit with the exception of replacement parts properly installed in accordance with the present instructions.
- If there is any reason to believe that the safety of the device may be compromised, the unit must be put out of operation and labeled appropriately to ensure the device is not inadvertently used by a third party. The safety may be compromised, for example, when the unit fails to function properly or shows visible signs of damage.
- The device should be opened and repairs be carried out only by qualified experts.

14. Only appropriately labeled original 3M replacement parts may be used in the repair and maintenance of this unit. Any damage resulting from the use of third-party replacement parts voids the liability of 3M.
15. Do not attempt to operate this unit by bypassing the protective mechanisms (e.g. by turning off the automatic shut-off switch on the cover of the unit). Violation of the safety regulations may cause serious injury and voids all warranties and guarantees.
16. Do not expose the device to solvents, flammable liquids, and strong heat sources as these may damage the plastic housing.
17. Cleaning the unit, ensure that no cleaning agent enters the device since this may lead to electric short-circuiting or hazardous malfunctions.
18. Pay attention to capsule leak tightness when mixing. Capsules that are not leak tight contaminate the device with mercury. Contaminated devices cannot be repaired and must be disposed of as contaminated waste according to local laws.  
See the safety data sheet for the amalgam used for first aid, measures to take in case of unintended release, handling, storage, exposure limitation, and personal protective clothing.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

## 2. Product Description

The CapMix unit is a mixing device for dental filling materials and cements in capsules. Desired mixing times may be programmed into the mixing programs "P1" - "P4". Moreover, two pre-set mixing times, "10 s" and "15 s", can be selected.

These Instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly readable. For details on all additionally mentioned products please refer to the corresponding Instructions for Use.

### Intended Purpose

Intended users: educated dental professionals, e.g. general dentists, dental assistants, dental hygienists, who have theoretical and practical knowledge on usage of dental products.

## 3. Fields of Application

Dental materials packaged in capsules with a max. capsule weight of 7 g.

Capsules with a weight of more than 7 g must not be used!

## 4. Technical Data

Drive unit:	Ball-bearing motor with integrated self-resetting overheat cut-off
Operating mode:	Short-term operation, 1 min On/7 min Pause
Controller:	Microprocessor-controlled, high-precision mixing time control by quartz timing mechanism
Capsule holder:	Self-locking
Capsules:	Diameter max. 19 mm, total length max. 58 mm, length without tubes max. 38 mm, weight max. 7 g
Housing:	Housing and protective cover made of high impact-resistant plastics
Mixing speed:	4650 rpm for 7 g-capsules
Operating temperature:	18°C...40°C/64°F...104°F
Relative humidity:	max. 80% at 31°C/88°F, max. 50% at 40°C/104°F
Main voltage:	100 V, 120 V, 230 V, 50/60 Hz, for preset voltage refer to the rating plate
Power requirements:	330 VA
Fuses:	230 V: 2 A time-lag fuse 120 V: 3.2 A time-lag fuse
Dimensions:	(W x H x D) 195 x 160 x 195 mm
Weight:	3750 g
Continuous sound pressure level:	>70 dB(A)
Year of manufacture:	see rating plate

## 5. Installation

### Prior to the initial start-up, remove the transport safety guard.

Connect the CapMix unit to the power supply and switch on the unit with the power switch on the back of the device. The display shows the mixing time of program 1.

The following basic settings are preset at the factory:

Program 1 (key "P1"):	5 sec mixing
Program 2 (key "P2"):	7 sec mixing
Program 3 (key "P3"):	9 sec mixing
Program 4 (key "P4"):	12 sec mixing
Fixed time (key "10 s"):	10 sec mixing
Fixed time (key "15 s"):	15 sec mixing

### 5.1 Setting the mixing time

Mixing times in the range of 1 sec to 39 sec may be selected. By pressing the keys "+" or "-" the device increases or decreases, respectively, the mixing time by 1 sec. Depressing and holding the "+" or "-" keys, the device continues to count up or down, respectively.


### 5.2 Changing the basic settings

Set the mixing time as described above, then depress one of the program buttons, "P1" - "P4", until two short signals are issued. This confirms that the desired settings have been stored in the selected program.


## 6. Operation

### 6.1 Insertion of the capsule

Open and fold back the cover of the CapMix device. Press down on the release button to release the capsule holder.

Keep the release button depressed and push to the left to open the capsule holder.

Place the activated capsule between the two cups of the holding mechanism and release the release button. The device clamps the capsule automatically. Insert capsules with tips such that the tip points towards the left side.

Do not operate the device without a capsule inserted!

### 6.2 Mixing

Close the protective cover and press the "Start/Stop"-button. This starts-up the mixing process. The countdown of the mixing time is shown on the time display.

After the mixing time has elapsed, the latest time setting entered is re-displayed so that there is no need to re-enter any value if the same time is to be used.

The selected program can be interrupted at any time by pressing the "Start/Stop"-button, upon which the program returns to its original setting.

**CAUTION:** For safety reasons, the unit cannot be started up unless the cover is closed. Attempting to operate the unit with the cover being open causes the time display to flash and an error message to be issued. Opening of the protective cover during mixing is associated with a high risk of injury. Hence, the mixing process is aborted. This may render the capsule useless! Therefore, never open the cover before the mixing process is completed and the capsule holder has come to a complete stop!

### 6.3 Removal of the capsule

Open the lid of the CapMix unit.

Open the capsule holder by pressing the release button down and towards the left, and remove the capsule.

## 7. Notes

Capsule materials of third-party manufacturers (including amalgam capsules) may be mixed with the CapMix device provided they are suitable for the device in terms of weight and dimensions (please refer to the section, Technical Data). The optimal mixing time for mixing with the CapMix device should be inquired from the capsule manufacturer.

Pay attention to capsule leak tightness when mixing. Capsules that are not leak tight contaminate the device with mercury.

## 8. Operating Errors

Error	Possible cause	Remedy
Time display is off, even though the unit is switched on	Power failure	Check both the power cable and main voltage
Motor fails to run when the unit is switched on	Cover of the unit is open	Close protective cover
	Motor is overheated	Allow the unit to cool

**NOTE:** The CapMix device should be repaired by 3M or an 3M-authorized repair facility only.


## 9. Maintenance and Care

The device requires no maintenance.

### 9.1 Fuse Replacement


Shut the unit off and disconnect the power cable. The fuses are located on the back of the unit next to the power switch.

With a screw driver press down on the small bracket of the fuse unit, squeezing the unit out a small distance.

Remove the fuse unit from the device.

Press the catch lightly in an upward direction so that the fuse holder can be pulled out. Remove the fuse from the fuse holder and replace with a new fuse.

Attention: Use only fuses with identical ratings.

Insert the fuse holder in the fuse unit until the catch audibly locks in place, then replace the fuse unit in the device. Make sure that the bracket of the fuse unit points towards the power switch of the device.

## 9.2 Care

Clean all components with a soft cloth and, if necessary, a mild cleaning agent (e. g. dish-washing detergent). Solvents or abrasive cleaners may not be used in any case, as they can damage the plastic components!

- Cleaning agents must not enter the unit!

To disinfect all components, spray the disinfectant on a towel and use it to disinfect the unit. Do not spray the disinfectant directly on the device.

- Disinfectant agents must not enter the unit!
- Dry residual disinfectants on the device with a soft and fluff-free cloth, as they damage the plastic components.

If necessary, ask the manufacturer of the disinfectant if its constant use will damage plastic surfaces.

## 10. Return of Old Electric and Electronic Equipment for Disposal

### 10.1 Collection

Users of electric and electronic equipment are required to collect their old equipment separately from other waste in accordance with the regulations of the specific country. Old electric and electronic equipment must not be disposed of with unsorted household waste. This separate collection is a prerequisite for recycling and reprocessing as an important method for preserving environmental resources.

### 10.2 Return and Collection Systems

When your CapMix is no longer usable, do not dispose of the device with household waste. 3M has set up special disposal facilities to handle the equipment. Details about the procedure for the specific country can be obtained from the pertinent 3M subsidiary.

## 11. Customer Information

No person is authorized to provide any information that deviates from the information provided in this instruction sheet.


### 11.1 Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

### 11.2 Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

## Symbol Glossar

Reference Number and Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.1.1 Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
ISO 15223-1 5.1.3 Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.5 Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6 Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.1.7 Serial number		Indicates the manufacturer's serial number so that a specific medical device can be identified.
ISO 15223-1 5.4.4 Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
CE Mark		Indicates conformity to European Union Medical Device Regulation or Directive.
Medical Device		Indicates the item is a medical device.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a dental professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1).
PAP 20		Indicates paper material is recyclable.
Recycle electronic equipment		DO NOT throw this unit into a municipal trash bin when this unit has reached the end of its lifetime. Please recycle.
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

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3M, ESPE and CapMix are trademarks of 3M or 3M Deutschland GmbH.
Used under license in Canada. © 2019, 3M. All rights reserved.
3M ESPE Customer Care/MSDS Information: U.S.A. 1-800-634-2249.
3M, ESPE et CapMix sont des marques de commerce de 3M ou 3M Deutschland GmbH.