Cercon® xt ML

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Instructions for use
Notice d'utilisation
Instrucciones de uso
Istruzioni per l'uso

THE DENTAL
SOLUTIONS
COMPANY™

 Dentsply
Sirona

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Instructions for use Cercon® xt ML

Product description:

Cercon® xt ML blanks are made of yttrium oxide- (yttria-) stabilized zirconium oxide (zirconia) (Y-TZP). They are build out of different shade layers which create a natural gradient of the restoration after sintering. They are used in fabricating frameworks for fixed prosthetic restorations.

Depending on the framework design, Cercon® xt ML frameworks can be ceramically veneered or delivered as fully contoured restorations. Which blank is selected will depend on the tooth shade to be reproduced and the space available for the veneer. With fully contoured restorations, no space is required for the ceramic veneer, which may allow the dentist to preserve more of the tooth substance during preparation.

Framework material	Zirconium oxide (Y-TZP)
Definitive cementation	<ul style="list-style-type: none"> • Adhesive cementing • Conventional cementing

Objects are individually fabricated to your digital design specifications such as anatomic contour, framework and wall thickness, connector diameter and cementing gap.

Technical specifications Cercon® xt ML:

- Type II, class 4 (pursuant to DIN EN ISO 6872:2015)
- CTE: 10.1 µm/m·K (25–500 °C)
- Modulus of elasticity: 210 GPa
- Flexural Strength: approx. 750 MPa (three-point flexural testing)

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Cercon® xt ML is indicated for all ceramic restorations for anterior and posterior locations:

- Crowns
- 3-unit bridges (up to the 2nd premolar)

Contraindications:

- This medical product may not be used in patient hypersensitivity to zirconia (Y-TZP) or one of the other ingredients
- Bruxism or recalcitrant parafunctional habits (for ceramically veneered frameworks)
- Insufficient available space
- Endodontic posts
- Endosseous implants
- Inlay bridges
- 3-unit bridges in the molar region

Warnings:

Possible cross-reactions or interactions of this product with other products or material already present in the oral environment must be taken into consideration by the dentist when selecting this product.

Precautions:

Please note:

- Keep product dust away from eyes.
- Avoid any contact with mucosa.
- After use, wash your hands and apply a hand cream.
- Do not smoke, eat or drink while handling the product.
- Do not swallow the product.
- Do not inhale dust particles during grinding.
- Use local vacuum suction and suitable mouth/face protection during manual machining at the workplace.

The safety and warning notes listed here describe how to use our product in a safe and risk-free manner. Notify the dentist in charge of all factors described above if you use this product for a custom design and make sure to comply with the pertinent Material Safety Data Sheets (MSDS).

Adverse effects:

If properly processed and used, adverse effects of this product are highly unlikely. However, reactions of the immune system (such as allergies) to substances contained in the material or localized paraesthesia (such as taste disturbances or irritation of the oral mucosa) cannot be completely ruled out as a matter of principle. Should you hear or be informed of any adverse effects – even when doubtful – we would like to request notification.

Any serious incident involving the product shall be reported to the manufacturer and the competent authority in accordance with local requirements. For a summary of the safety and clinical performance of this product, see <https://ec.europa.eu/tools/eudamed>

Special notes on framework design

Framework design for fully contoured restorations:

Indicated especially where limited occlusal space is available, for frameworks that are not to be veneered or for frameworks with body stains.

The surface of fully contoured frameworks can carefully be optimized with rotary tools such as fine cutters prior to sintering.

Make sure that the occlusal surface is not changed by subsequent deepening of fissures, as notching action may compromise the strength of the material. Please note that flat occlusal reliefs can extend the life expectancy of fully contoured restorations. When making manual adjustments, make sure never to separate the interdental spaces of the frameworks with cutting discs or other rotary instruments, doing so might damage the framework and compromise the strength of the material!

Important note:

Please ensure that the minimum wall thickness of the framework in the area of the occlusal surface is respected even after occlusal adjustments.

Framework design for veneered restorations:

Frameworks to be ceramically veneered are designed to reduced anatomical contour to provide maximum support for the veneer.

The frameworks can be veneered using the press-on technique or the build-up technique.

Wall and border thickness:	Cercon® xt ML
Wall thickness, single crowns	0.7 mm
Margin thickness, single crowns	-
Wall thickness, bridges	0.7 mm
Margin thickness, bridges	-
Additional dimensional requirements for the anterior region:	Cercon® xt ML
Number of pontics	1
Connector cross-section	12 mm ²
Additional dimensional requirements for the posterior region:	Cercon® xt ML
Number of pontics	1
Connector cross-section	16 mm ²

Nesting

Please consider the 4 different layers for an optimal nesting:

Disc Height	14 mm	100 %	18 mm	100 %	25 mm	100 %
Incisal	3.6 mm	25.70 %	3.6 mm	20 %	3.6 mm	14.40 %
Transition Layer 2	1.8 mm	12.85 %	1.8 mm	10 %	1.8 mm	7.20 %
Transition Layer 1	1.8 mm	12.85 %	1.8 mm	10 %	1.8 mm	7.20 %
Dentin	6.8 mm	48.60 %	10.8 mm	60 %	17.8 mm	71.20 %

Positioning as an example in an 18 mm disc, unsintered.

Entering the shrinkage factor for the milling unit:

In CAM software allowing the entry of 3 dimensions, enter X, Y and Z values.

In CAM software allowing the entry of 2 dimensions, enter X or Y and Z values.

In CAM software allowing the entry of 1 dimension, enter the X value.

Separating

Notes on separating the objects:

Separate the objects from the disk by sandblasting with aluminium oxide (50 µm, max. 1.5 bar).

Sintering

Sintering in the Cercon® heat plus P8:

- 1500 °C in the Cercon® heat plus P8
 - Program #4 for bridges up to 3 units, $T_{max} = 1500$ °C

Sintering in the heat DUO/Multimat2Sinter:

- Program #6: Speed program for bridges up to 3 units, $T_{max} = 1540$ °C
- Program #7: Standard sintering for bridge frameworks up to 3 units, $T_{max} = 1520$ °C

Sintering in the inLab Profire:

- #1 Cercon base_ht_xt Speed: Speed program for bridges up to 6 units (Cercon® base / Cercon® ht); for bridges up to 3 units (Cercon® xt, Cercon® xt ML),
- #2 Cercon base_ht_xt Standard: Standard sintering for bridge frameworks up to 8 units (Cercon® base / Cercon® ht); for bridges up to 3 units (Cercon® xt, Cercon® xt ML)

Sintering in third-party furnaces:

Sintering results may be adversely affected by e.g.:

- Incorrect sintering temperatures
- Insufficient heating power
- Incorrect temperature curves
- Incorrect object placement
- Insufficient heat-storing capacity of the furnace over the time of the sintering cycle
- Manufacturer-related or age-related variations in heating performance
- Object contamination by oxidation products emitted by non-enclosed heating elements

Any of these factors by itself or in combination may reduce the maximum strength of our above-mentioned zirconium dioxide materials and compromise the life expectancy of the frameworks.

For these reasons, we cannot grant general approval for the use of third-party furnaces for sintering crown and bridge frameworks made of Cercon® xt ML. We will, however, open the system technically for the use of furnaces of third-party furnaces only on condition that the following requirements are met:

Additional notes:

The programming you use for your sintering furnace should be analogous to the Dentsply Sirona sintering programs.

General Sinter programs Cercon® xt ML all shades

Material:	Start-Temp. °C	Ramp time min	Temp. 1 °C	Holding time min	Ramp time min	Temp. 2 °C	Holding time min	Cooling
Sinter program for bridge frameworks up to 3 units	RT ¹⁾	40	900	0	55	1500	145	With closed furnace cooling down to 200 °C
Speed-Sinter program for bridge frameworks up to 3 units	RT ¹⁾	90	1540 ²⁾	35	20	1150	0	Gradual opening of the furnace within 35 Min down to 200 °C

Material:	Start-Temp. °C	Ramp time °C/min	Temp. 1 °C	Holding time min	Ramp time °C/min	Temp. 2 °C	Holding time min	Cooling
Sinter program for bridge frameworks up to 3 units	RT ¹⁾	22	900	0	11	1500	145	With closed furnace cooling down to 200 °C
Speed-Sinter program for bridge frameworks up to 3 units	RT ¹⁾	17	1540 ²⁾	35	18	1150	0	Gradual opening of the furnace within 35 Min down to 200 °C

Material:	Start-Temp. °C	Ramp time °C/h	Temp. 1 °C	Holding time h:min	Ramp time °C/h	Temp. 2 °C	Holding time h:min	Cooling
Sinter program for bridge frameworks up to 3 units	RT ¹⁾	1320	900	00:00	660	1500	02:15	With closed furnace cooling down to 200 °C
Speed-Sinter program for bridge frameworks up to 3 units	RT ¹⁾	1020	1540 ²⁾	00:35	1080	1150	00:00	Gradual opening of the furnace within 35 Min down to 200 °C

1) Room temperature

2) valid for closed sinter bowls, otherwise 1520°C

Sintering programs, Multimatt2Sinter / heat DUO / inFire HTC-speed sintering furnace:

Cercon® xt ML speed sintering:

Step	Heating rate °C/min	Temperature °C	Holding time min
S4	99	750	0
S3	7	1520	35
S2	14	1350	0
S1	23	1100	0

Cercon® xt ML standard sintering:

Step	Heating rate °C/min	Temperature °C	Holding time min
S4	25	0	0
S3	99	300	0
S2	11	1500	135
S1	22	880	0

Sintering temperatures are recommendations. If necessary, carry out a trial sintering cycle and adapt the sintering temperatures or times as needed.

Manual finishing after sintering:

- Sandblast the inside and outside of the framework with aluminium oxide (110-125 µm, max. 2-3 bar, 45 ° angle).
- Eliminate premature contacts one by one until the framework has reached its final position on the die(s).
- During the try-on and fitting of the framework, keep the dies on the cast and try on the framework as a whole.
- Once the try-on and fitting has been completed, do not perform any additional adjustments such as finishing the entire framework.

Note: Sintered zirconia should be finished using diamond instruments under proper irrigation only. Keep the pressure on the framework material to a minimum and work only in one direction.

- Sandblast the finished areas once again with aluminium oxide (110-125 µm, max. 2-3 bar, 45 ° angle).
- Finally, clean the framework using a steam cleaner.

Staining technique

For individualization of fully contoured Cercon® xt ML restorations we recommend Dentsply Sirona Universal Stain & Glaze.

Layering technique

We recommend Cercon® ceram Kiss / Celtra® Ceram ceramic veneering materials for veneering zirconia frameworks. Make sure to follow the applicable Instructions for Use.

Annealing:

Based on the results of our scientific testing of zirconia frameworks, we consider a separate annealing (“healing”) step both unnecessary and inappropriate.

Surface treatment in the dental lab:

Unveneered Cercon® xt ML frameworks should be highly polished or glazed with a high-gloss ceramic glaze. This also makes the restoration more amenable to proper oral hygiene procedures.

Polishing in the dental office:

Extensive studies have shown that the abrasive action of Cercon® xt ML on antagonists is less than that of conventional ceramic veneers and not more than standard lithium disilicate ceramics even after finishing and polishing. Important: After making minor final occlusal adjustments at chairside, we recommend prophylactic polishing of the adjusted sites to a high lustre or adding glaze with a smooth surface prior to temporary or definitive cementing to protect the antagonists from possible abrasion. This also makes the restoration more amenable to proper oral hygiene procedures. Tooth surfaces whose shade was obtained by staining may present brighter spots where occlusal adjustments were made.

Transport and storage: No particular requirements.

Shelf Life: Shelf life of Cercon® xt ML is 7 years from date of manufacture.

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DeguDent GmbH