BIOSTAR®

Bedienungsanleitung
Operation Manual
Mode d'emploi
Manual de Instrucciones
Istruzioni per l'uso

SCHEU
SCHEUGROUP

Abb./Fig. 1



Abb./Fig.2/230V

Abb./Fig.3/100,115V

Abb./Fig. 4

Technische Daten / Technical data

Spannung / AC 230 V **7** **8**, 115 V, 100 V **7** **18** / 50/60 Hz

Leistung / Power 750 VA

B x H x T / W x H x D 520 x 260 x 440 mm

Gewicht / Weight 14 kg

Druck / Pressure **9** max. 6 bar / 87 psi

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1. Safety Instructions

Carefully read the operation and safety instructions before commissioning and connecting the device, ensuring that all requirements are met.

General information

- Inspect the machine for damage before commissioning. If any damage is found, it shall not be operated.
- Always supervise the machine when it is switched on.
- Only use blanks of a diameter of 125 mm!
- Keep the device out of reach of children!
- Avoid any contact with liquids in order to protect the device against electric shock.
- Infrared heater and heater arm reach high temperatures during operation.
- Do not operate the device close to highly flammable materials.
- Ensure to operate the device only using its handles to avoid the risk of burns.
- The device shall only be transported when it has cooled down.
- Make sure that the heater is not exposed to any load during transport.
- Safe operation of the device can only be ensured using the original SCHEU-DENTAL spare parts.
- Any kind of service operation shall always be performed by service technicians instructed by SCHEU-DENTAL or at the SCHEU-DENTAL facilities. Any conversion or modification of the device compromises product safety and leads to the expiration of warranty and guarantee claims.
- Never open the device without authorization. Caution: Risk of injury!
- Never carry out repairs yourself!
- The device is designed for electricity networks with alternating voltage of 100 V, 115 V or 230 V (please check the declared voltage on the serial number on the back of the unit).

Installation of on-site power supply

The installation of power supply for the connection of our equipment shall be carried out in accordance with international regulations and the regulations derived thereof.

Our equipment is built according to VDE (German Association for Electrical Engineering, Electronics and Information Technology), protection class I and shall be connected to a socket with protective conductor.

Environment

The BIOSTAR® is made from reusable or recyclable materials. In accordance with the 2012/19/ EU Waste Electrical and Electronic Equipment Directive, SCHEU-DENTAL will take back the unit free of charge if it is returned by free delivery. Alternatively, the machine can be sent to a designated collection point for recycling of electrical and electronic equipment.

2. General Commissioning

⚠ For optimal use of the device, ensure to understand and follow the below instructions. When setting up the BIOSTAR® in its working environment, make sure that the infrared heater ①, the pressure chamber ② and the locking shaft ③ are freely movable. Avoid setting up the machine in an environment close to steam jet and sand blasting units or places with high plaster and dust exposure.

⚠ Fill the stainless steel pellets provided with the machine in the model cup ④ and the remaining quantity in the pellet drawer ⑤; when using the model platform ⑥ ensure its surface in the model cup is free from pellets. Ensure that the surfaces above and beneath the model cup are always free from pellets during the thermoforming process.

⚠ Please note that it is not recommended to use any other material to embed the model than the provided stainless steel pellets (e.g. rice or sand). Those materials may lead to damages of the pneumatic parts and exhaust air openings by swirling during pressure built-up and depressurization. Moreover, their weight may not be appropriate for precise covering and embedding the models. Ensure to use solely the original power and air connections provided with the delivery. Make also sure that the compressed air supply is completely free from oil residues and humidity and has a minimum feed pressure of 6.0 bar.

3. System settings

You can customize the settings of language, pressure, temperature and screen contrast to your personal needs.

Switching on the device

Press the “CE+home” button five times, once the current device version appears on the display.

Selecting the language

Select digit – confirm with “code/scan”

Selecting the pressure

Select digit – confirm with “code/scan”

Selecting the temperature

Select digit – confirm with “code/scan”

Selecting the screen contrast

Select digit – confirm with “code/scan”

4. Info menu

The following information is given in the info menu:

Software: Display of the installed firmware version.

Safety instructions: Here you will find important instructions on the safe operation of your pressure moulding device.

Indications: Discover the variety of applications in the digitized version of the application booklet "pressure moulding technique".

Troubleshooting measures: Find quick and easy solutions to problems in your working process.

Range of pressure moulding material: Get an overview on the wide range of consumables with specific information on each kind of material.

Videos: Fresh up your knowledge by watching our short videos on the fabrication process of various kinds of pressure moulded appliances.

After starting the device, you can access the info menu by pressing the button "info" **20**. The device must be in the home position (see 7).

The keypad is equipped with various keys, some of them will light up to indicate that specific information is available.

You can navigate through the system pressing the following keys:

Info button 20: Accessing the info menu

Note!

If you have already scanned a material, only specific information on this material will be available. You can not access the info menu during the various steps of the thermoforming process (heating, pressure moulding, cooling); the info menu is hidden. Information and instructions on the thermoforming process have priority and are displayed, once infrared heater, pressure chamber or locking shaft are moved from the home position.

← and → buttons: The horizontal navigation keys are only available once you accessed "indications", enabling you to navigate through "indications".

OK button 22: This key confirms the current selection.

↑ and ↓ buttons: The vertical navigation keys lead you through the current menu.

↶ button 21: By pressing the **return button**, you can return to your previous selection.

CE+home button: Pressing this button, you can return to the menu item "material selection according to indication".

5. Network connection

Connect the BIOSTAR® device to the power supply (please also refer to installation of on-site power supply as described in the safety instructions) by plugging the power cable provided with the machine into the power socket **7** on the back of the device. For the device version in 230 V, two additional device fuses are built in the main socket protecting the electrical parts from overvoltage. In order to exchange the fuses, press the small clip underneath the power socket and remove the fuse box **8**. For the device versions in 110 V/115 V these fuses can be found in two separate holders **13** that can be unscrewed. For the device version in 230 V the fuse TAA (REF 3100.940), for the device version in 115 V the fuse T10A (REF 3101.341) and for the device version in 100 V the fuse T10A (REF 3100.941) are used.

6. Compressed air supply

A 20 bar air hose **9** with filter element and hose nozzle (DN 7.2 / 6 mm inner tube diameter) is permanently connected to the device and shall not be replaced. Connect the BIOSTAR® device and the air hose provided with the machine to your compressed air system, using optionally the quick coupling (DN 7.2 / 6 mm inner tube diameter) or a screw connection. The operating pressure is set at the factory to 6.0 bar/87 psi. This pressure guarantees best thermoforming results and is being continuously indicated on the display; it shall not be modified. In case a reduction should be required, you can adjust the pressure using the pressure regulator **10** on the back of the device. The pressure regulator is unlocked and turned by pulling (↶↷). After adjustment, the pressure regulator shall be locked again by pressing it in. In case pressure is increased exceeding 6.5 bar/94 psi (which does not correspond to our recommendation), the safety valve starts to operate, which is indicated by a loud venting noise.

7. Home position

In the home position, the infrared heater **1** is swivelled backwards, the pressure chamber **2** is being opened to the left, and the locking shaft handle **3** is being swivelled backwards to the stop / 12 o'clock position (see fig. 1).

8. Switching on the device

The device is switched on using the on/off switch **13**. A test and setting mode is run for approx. 10 sec, before the standard dialogue appears that is always the start dialogue for any thermoforming process. Programming is done using the keyboard **15** or the scanner **16** and is indicated on the display **14**.

9. Fixing the blanks

You can process round blanks with a diameter of 125 mm in thicknesses from 0.1 mm to 5.0 mm in the BIOSTAR® machine. Open the pressure chamber **2** by swivelling it to the left and remove the locking ring **11** by turning its handle to the 5 o'clock position. Place the blank on the sealing ring of the balancing ring **12** that is visible in the pressure chamber and lock the pressure chamber by refitting the locking ring and turning its handle to the 6 o'clock position.

10. Preparing the model

When using hard or hard-elastic material, the model is placed in the model cup ④ that has been filled with pellets to its upper sealing rim. When using soft-elastic material, the model is placed on the model platform ⑥. Ensure in both cases that the front teeth are facing to the right in the direction of the locking shaft ③ or to the left in the direction of the pressure chamber ② in order to avoid a different material thickness between the quadrants.

11. Programming using the scanner or by entering the code or heating time

A three-digit code is assigned to all our pressure moulding blanks, containing all required information on temperature, heating and cooling time. This code is printed in plain text and as bar code on all our material except for material with a thickness of 0.1 mm or 0.15 mm. Once you have selected the desired material, press the button **"code/scan"**. The scanner ⑬ is now activated. Now present the bar code of the selected material under the scanner until you hear an acoustic signal. You also have the option to enter the code manually: Enter the three digits and confirm by pressing **"code/scan"** again. Another option is to enter the heating time manually. Doing so, press the key with the clock symbol ⑰, enter the required heating time and confirm by pressing key ⑰ again.

12. Heating

Once you have scanned, entered the code or the heating time, swivel the heater ① over the fixed blank. The temperature will be automatically adjusted by a thermo element, which is visible by a change of the heater colour.

13. Thermoforming

Acoustic signals indicate the end of the heating process (+5 sec to 0 sec = beeps / 0 sec to -5 sec = continuous tone). Once 0 sec are reached, the heater shall be swivelled backwards, upon reaching -5 sec, the heater is switched off automatically for safety reasons. Swivel the pressure chamber ② by a 180° turn to the right (**pay attention to clean sealing surfaces**) and turn the locking shaft ③ 180° to the front. The pressure chamber is filled with air and the cooling time starts automatically which is indicated by a blue LED.

14. Depressurizing

Optical signals (blue flashing LED) and acoustic signals (4 beeps) indicate the end of the cooling phase. Pressing the flashing button **"air"** depressurizes the pressure chamber, 4 further beeps after approx. 3 sec indicate the end of depressurizing. Turn now the locking shaft ③ 180° backwards, push the locking ring ⑩ ← to the left and open the pressure chamber ② ↻. Once the model has been removed, the device can be used for the next pressure moulding process.

15. Special programs for the use of the heating support

You have the option to program time intervals of 60/120/180/240 sec. For this, the pressure chamber shall be closed and the locking shaft shall be opened. These time intervals can be used for heating and bonding the occlusal surfaces using the heating support (REF 3452). Once you have entered the corresponding heating time, the heater with the attached heating support is swivelled forwards in the heating position. This operation can be carried out several times depending on the material thickness.

16. Service/maintenance

The BIOSTAR® device is generally maintenance-free. Ensure however to control the air tube with the filter element for any kind of debris or moisture. The filter element can be exchanged if necessary. The model cup ④ and the exhaust air holes as well as the model platform ⑥ can be cleaned by means of air nozzle or steam jet. Ensure the reflective surfaces in the pressure chamber are always clean for optimum heating results. The aluminium balancing ring ⑫ in the pressure chamber ② can be removed for cleaning by tipping. Doing so, clean and grease the sealing surfaces, the lateral sealing ring, the 4 rivets and springs. Make sure to refit the balancing ring on the springs and rivets with the sealing ring facing upwards. The whole unit shall be freely movable (see also Troubleshooting Measures).


17. Service updates

The following applies to devices with serial numbers beginning with 70858 for the 230 V version, with 70527 for the 115 V version and with 70103 for the 100 V version:

You have received a registration card with your device. Once you are registered, you will be automatically informed about updates that can be downloaded from our homepage <https://www.scheu-dental.com/en/service/support-and-program/>. In order to access the SD card (see fig. 4) unscrew the cover on the left side of the device. Remove the SD card and erase any existing data or reformat your SD card. Extract the download file on the computer and copy the subdirectories "data" and "fw" to a blank and formatted SD card. Insert the SD card in the **service interface ⑱** on the left side of the device. Refit the cover of the service interface on the device. The update for the SD card is now complete.

The following applies to devices with serial numbers beginning with 70001-70857 for the 230 V version, with 70001-70526 for the 115 V version and 70001-70102 for the 100 V version:

In these cases, you can use a USB flash drive for the updates: Ensure to use a blank and formatted USB flash drive. Extract the download file on the computer and copy the subdirectories "data" and "fw" to the USB flash drive. A video showing the procedure is available in the support area. Start the device press the **"ok"** button 5 times. You are asked to enter the password for the info menu. Enter the digit 1 and confirm by pressing **"ok"**. By pressing the **digit "1"** again, the update starts. Once the update is complete, quit the service menu pressing the key **"CE+home"**.

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Due to a safety-related technical modification the pressure chamber of your pressure moulding machine has been equipped with a drill for decompression.

The venting hole is sealed by locking the closed pressure chamber using the locking shaft. Please make sure that the sealing surface is always slightly greased. This is the only way to ensure easy locking and low wear of the seal.

A worn seal can be simply pulled out from the sleeve holder of the BIOSTAR®.

Before refitting a new seal, ensure to also grease the outer lateral surface of the seal (see pic.4)

You can order the brass pressure piece (6) with screw (5) using the REF 3110.828. Your accessories set also contains 3 replacement seals for the BIOSTAR® in 10 mm (REF 3000.856) and grease with pellet brush (REF 3005.800).

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ⒼB Warranty terms and conditions for the following machines of SCHEU-DENTAL GmbH

The machine you bought has been carefully checked and found to be in perfect condition, meeting highest technical demands.

SCHEU-DENTAL GmbH does not provide warranties for:

- repair services effected by the buyer or unauthorised staff
- in case of damages due to inappropriate use or handling
- in case of a failure to comply with the handling, cleaning and maintaining prescriptions
- in case of inappropriate connection or incorrect installation
- in case of force majeure or any other influences beyond the control of SCHEU-DENTAL
- in case of any modifications to the machines
- in case of using spare parts of other manufacturers
- in case of signs of wear due to age, not affecting the main functions of the machines.

Utilisation of warranties does not extend or restart the standard warranty period. In case of unjustified utilisation of our warranties we reserve to charge a fee for the corresponding services.

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