
- GB** CORDLESS LED CURING LIGHT
- F** LAMPE À PHOTOPOLYMERISER À LED SANS FIL
- I** LAMPADA DI FOTOPOLIMERIZZAZIONE A LED CORDLESS
- CZ** BEZŠŤŮROVÁ LED POLYMERIZAČNÍ LAMPA
- H** VEZETÉK NÉLKÜLI LED POLIMERIZÁCIÓS LÁMPA

# LED-POLYMERISATIONS-LAMPE

Daily dental excellence.

|                                |        |       |
|--------------------------------|--------|-------|
| <b>D</b> BEDIENUNGSANLEITUNG   | Seite  | 2-13  |
| <b>GB</b> INSTRUCTIONS FOR USE | Page   | 14-25 |
| <b>F</b> MODE D'EMPLOI         | Page   | 26-37 |
| <b>I</b> ISTRUZIONI PER L'USO  | Pagina | 38-49 |
| <b>CZ</b> NÁVOD K POUŽITÍ      | Strana | 50-61 |
| <b>H</b> HASZNÁLATI ÚTMUTATÓ   | Oldal  | 62-73 |

Stand 020821

 BG LIGHT LTD · 155, Vasil Aprilov blvd. · BGR · 4027 Plovdiv · Tel.: +359 32 644089

 Müller & Weygandt GmbH · Reichardsweide 40 · D-63654 Büdingen, Germany · Tel.: +49 (0) 60 42 - 88 00 88

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### ACHTUNG!

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## INSTRUCTION FOR USE

### ATTENTION!

Before installation and starting work with the unit, carefully read this guide!

### CONTENTS:

- I. UNIT DESCRIPTION AND FUNCTION
- II. TECHNICAL DATA
- III. ACCESSORIES AND SPARES
- IV. SAFETY PRECAUTIONS
- V. PREPARATION AND SEQUENCE OF OPERATION
- VI. ROUTINE CARE AND MAINTENANCE
- VII. WARRANTY CONDITIONS
- VIII. SYMBOLS
- IX. FAQ
- X. INFORMATION ON A.M. Edelingh LED-curing light
- XI. SERVICE DATA

### I. UNIT DESCRIPTION AND FUNCTION

A.M. Edelingh LED curing light cordless is a light source for intraoral polymerization of dental materials. The unit uses LED (Light Emitting Diode) module as the light source. The emitted light is in the wavelength range 410-490 nm.


The appliance should only be used by a dental practitioner and in a dental practice. The device consists of Handpiece, Stand and a Power Adapter.

A.M. Edelingh LED curing light is produced in conformity with MDD 93/42/EEC and 2007/47/EC.

### II. TECHNICAL DATA

1. Operating voltage 3.7V Li Ion battery (2600mAh)
2. Power consumption - mains adaptor: 110-240VAC/50-60Hz, 0.2-0.1A.  
Output: 5VDC/ 1.2A micro USB socket.
3. Dimensions:
  - Emitting handpiece head -14 x 16 mm,
  - max diameter of handpiece, - 30 mm length -215 mm
4. Weight of Emitting handpiece -215 g
5. Light intensity:
  - In RAMP mode - up to 1500 mW/cm<sup>2</sup> - constant intensity, regardless of battery charge level.
  - In HYPER/ORTHO mode - up to 3500 mW/cm<sup>2</sup> - constant intensity, regardless of battery charge level.
6. Emitted time:
  - In RAMP mode - 10 / 20 sec./±10 %
  - In HYPER mode 3/3+3\* sec./±10 %
  - in HYPER ORTHO mode -3/ 36\* sec./±10 % (\*with 1 sec pause)Start/Stop button, can be stopped in every moment by second press of the button.
7. Emitted light - blue, 410-490 nm.
8. Light source - 4 powerful LED crystals in module with reflector.
9. Charge time - 6-8 hours for full charge and 3-4 hours to 80% of full charge.

10. Working ability after full charged battery - up to 4000 sec. (400 cycles x 10 sec.)
11. visual indication - color RING.
12. low battery warning. Fully discharged battery indication.
13. working mode (1 min./ 10 min. off)
14. broadband LEDs covering all photoinitiators in 410-490 nm wavelength range.


The manufacturer of this unit declares to provide on request all additional necessary technical documentation / information which will help the user's technical staff to service the parts of the unit which the manufacturer has claimed to be a subject for repair.

### III. ACCESSORIES AND SPARES

1. Charge adaptor with micro USB output 5V DC / 1.2A - 1 pc.
2. Cordless handpiece - 1 pc.
3. Battery (inside the unit) - 1 pc.
4. Holder - 1 pc.
5. Barrier sleeves - 1 pack (50pcs)
6. Operations guide - 1 pc.

### IV. SAFETY PRECAUTIONS

As with all technical devices, the proper function and safe operation of this unit depend on the user's compliance with the standard safety procedures as well as the specific safety recommendations described in this Operations guide.

1. The unit must be used in strict accordance with the operations guide. The manufacturer accepts no liability for any damage resulting from the use of this unit for any other purposes.
2. Before start-up of the unit make sure that the operating voltage indicated on the adapter label corresponds to the mains voltage. Operation of the unit at a different voltage may damage the unit. Electrical safety is ensured by Protection class II against electrical current damages according to EN 60601-1. Use of "Shuko" plug or unit's grounding is not necessary.
3. A.M. Edelingh LED curing light should be operated in indoor premises, in the following conditions:
  - temperature +10° - +40° C;
  - relative air humidity - 30 - 75%;
  - atm. pressure - 700 - 1060 hPa;
  - absence of chemically active or flammable substances.

4. Irradiation of the eyes bears certain risk for their damage. Therefore, the light must not be directed toward the eyes! Exposure must be restricted to the working area. Use of protection glasses for operator and patient that blocks effectively the blue spectral range is highly recommended!

This curing unit produces high-output curing energy! A significant increase in curing energy is possible compared with equipment previously used. It is important to observe the following precautions and procedures:

- Do not place light directly on or towards unprotected gingiva or skin.
- Adjust curing techniques in accordance with the increase in curing energy.

ATTENTION! The emitted light is in the wavelength range 410-490 nanometers. Use protection goggles during work!

5. The unit can be used only with doctor's advice on persons with heart pacemakers implant; persons with photobiological reactions; persons currently taking photosensitive medication; persons with cataract surgery; persons with retinal diseases, etc.

6. Solvents, flammable liquids and powerful heat sources must not be kept near the unit as exposure to one of these factors may damage the plastic parts of the unit.

7. Cleaning agents should not enter the unit as this may short-circuit the unit or cause potentially dangerous malfunction.

8. Only factory-authorized technicians are authorized to open and repair the unit.

9. Only original A.M. Edelingh LED curing light spare parts should be used to replace defective parts. The unit warranty does not cover any damage resulting from the use on non-original replacement parts. The unit or any of its parts should not be disassembled while it is connected to the mains!

10. A.M. Edelingh LED curing light must be used with barrier sleeves to avoid cross contamination and liquid damage. Dispose the barrier sleeve after each patient procedure.

11. **FRAGILE!** Be extremely careful during transportation, storage and handling of the unit!

12. Pursuant to Directive 2012/19/EU, this symbol shows that the product must not be disposed of as urban waste at the end of its operating life. The product must be taken to centers specialized in separate collection of electric and electronic equipment. Correct disposal of the equipment, which is no longer used, prevents negative consequences for the environment and human health!

## V. PREPARATION AND SEQUENCE OF OPERATION

A.M. Edelingh LED curing light is designed for use as a battery powered unit

1. Take the curing light out from the packing.  
Attention! The unit should be placed at least 1 meter away from any heat source.
2. Put a new barrier sleeve on the front part of the handpiece.
3. Connect the power adapter to charge battery if needed.


**IMPORTANT!** Connect power charger in proper way!

A.M. Edelingh LED curing light has 2 curing modes. They are indicated on the color ring on handpiece.

### When ring is lit

RED - HYPER mode;  
GREEN - RAMP mode;  
BLUE - change of mode, sound signal.  
YELLOW - thermal protection activated.

First select desired curing mode by pressing and holding the Start/Stop button.

The device has soft start (gradual increase in LED module intensity) in standard curing mode, which helps prolong LED module's life and at the same time avoid stress and shrinking in photocomposite.

4. The unit is activated by pressing the Start/Stop button of the handpiece. By second pressing light can be stopped at any moment.


### Timer selection:

#### - In RAMP mode:

- press once (1x) for 10 sec;
- press twice (2x) for 20 sec;

#### - In HYPER mode:

- press once (1x) for 3 sec;
- press twice (2x) for 3+3 sec;

Blue light appears at the end of handpiece head.


In HYPER mode, after each 3 / 3+3 sec. cycle, there is a 3 sec. pause, when indication ring blinks, and the device cannot be activated until it stops blinking. This pause is created for dental pulp protection from overheating.


**ORTHO version** has difference in HYPER mode - HYPER ORTHO mode - double (2x) pressing of Start button = 12 cycles x 3 sec. Curing position must be changed every 3sec to avoid overheating and damage of tissues!

5. A.M. Edelingh LED curing light has RAMP and HYPER mode of light curing that guarantee full and fast hardening of the restorative material:

**RAMP mode (10 / 20 sec.)** - During the first 3 sec. light intensity slowly increases to 1500 mW/cm<sup>2</sup> - "soft start" to avoid stress in the restorative material.

**HYPER mode (3 / 3+3 sec.)** - In this mode the light intensity is at maximum level - 3500 mW/cm<sup>2</sup> that ensures fast polymerisation and is suitable for work with liquid composites.

6. Direct the emitting light of handpiece to the curing spot at the closest distance (2-7 mm) and make sure the emitting head does not touch the photocomposite during the first several seconds and after that, when surface layer is hardened enough, the head can touch the photocomposite.

7. After second pressing of the Start/Stop button or after the selected 10 / 20 or 3 / 3+3 sec. exposure the unit is deactivated.

8. Switching off of curing light in the end of the day shall be done by unplugging of adapter from the wall socket.

## Overheating protection

### 9. Overheating protection:

The unit is equipped with overheating protection, which is activated if the temperature of the handpiece rises up to 40°C. When overheating protection is activated the color ring flashes in yellow for the next 60 sec. and is then deactivated until the temperature of the handpiece reaches 35°C.

### 10. BATTERY MAINTENANCE, CHARGING AND REPLACEMENT

The cordless handpiece is equipped with electronic charging control that ensures optimal battery charge and helps avoid overcharging.

**10.1** Charge adaptor is SYS1421-0605-W2E - 100-240VAC input and 5VDC output - micro USB type. It is not allowed to use another type power supply.

#### 10.2 Low Battery indication:

-rarely blinking of lightguide shows that battery is almost discharged. Handpiece can be used for several minutes, light intensity remains stable.

**10.3** No danger for battery overcharge - inside is battery protector. When battery is fully charged, the handpiece can work for approximately 4000 sec.

**10.4** Curing light can be used during charge time.

#### 10.5 To replace the Li-Ion battery:




- Place the battery end of handpiece on a non-slip surface and unscrew the cylindrical part 3 of the handpiece. Carefully remove the coupling with the battery out.

- Install new battery in reverse order.

- To assemble handpiece, place the battery end of handpiece on a non-slip surface, then hold firmly part 1 to fit part 2 together. In the same time screw cylindrical part 3 back to front cone - part 1.

**Note: The fitting of parts 1 and 2 - tooth and groove - must match! Battery cable must not be rotated in any way! Use only original A.M. Edelingh LED curing light Li-Ion battery!**

**10.6** The battery must not be immersed in water or thrown in fire, plus- and minus pole must not be short-circuited. Utilize according to local requirements. When discarding old battery, observe environment protection rules!

## VI. ROUTINE CARE AND MAINTENANCE

### Cleaning of the unit

For disinfection of the unit, spray the disinfection agent onto a piece of soft cloth/ cotton and clean the handpiece and power adapter. Do not spray directly on handpiece. Prevent liquids from entering openings of the unit. Do not immerse units in disinfectants. Do not use abrasives or solvents as these may damage the unit!


Suitable cleaning agents: products from M+W Dental, Isopropyl alcohol-based cleaners, Ethyl alcohol-based cleaners, Lysol disinfectant, Other no-bleach and non-abrasive disinfectants or cleaners.

Do not spray directly onto handpiece or in button / charge connector!

### 3. Infection control

To prevent cross-contamination, it is obligatory to use a disposable barrier sleeve over the handpiece for each patient.

## VII. WARRANTY CONDITIONS

1. Warranty period of the A.M. Edelingh LED curing light is 24(twenty four) months from the date of purchase. If date of purchase is not written, warranty starts at date of production.

2. During the warranty period the replacement of the defective parts (LED module and battery - 6 months warranty) is done free of charge by manufacturer's service.

3. If during the guarantee period the unit is damaged, due to incorrect operation (mechanical, chemical, thermal or electrical damages), damage caused by improper use, storage or any other reasons in user's fault, the repair is paid by the user. If any damage is noticed on electric cables, the device must be brought immediately to manufacturer's service. The device must not be used if cables are damaged. If liquids, such as water and solvents, aggressive or flammable substances and their vapors, also pouring or wiping the device with the same, insects or rodents have entered the device, it must be brought immediately to manufacturer's service. In such cases warranty is lost. Damages and claims are not accepted due to electric shocks, thunder storms, non-following of measures for electric safety or insufficient protection of patients, personnel or outer persons toward light irradiation. The manufacturer does not owe compensation for lost benefits in cases of defects or imperfect functioning of the device regardless of the cause. No claims and damages are accepted due to incorrect security and attention for protection during transportation, unpacking, moving, work and preservation of device. Warranty becomes void if any of the above mentioned happens.

If a dispute arises on the application and interpretation of this Operations guide, it shall be solved according to the applicable laws in Germany, and the court of Friedberg, Germany shall be the court of jurisdiction.

4. This warranty becomes void if any repair to or servicing of the product is performed by any unauthorized party or if other than original A.M.Edelingh LED curing light spare parts are used therefore.

5. The warranty servicing is performed at the manufacturer's service at the following address:

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## VIII. SYMBOLS

|   |  |
|---|--|
|  | Manufacturer<br>BG LIGHT LTD, 155, Vasil Aprilov blvd.,<br>4027 Plovdiv, Bulgaria  |
|  | Distributed by:<br>M+W Dental GmbH, Reichardsweide 40,<br>D-63654 Büdingen, Germany  |
|  | Pursuant to Directive 2012/19/EU,<br>this symbols shows that the product must not be disposed<br>of as urban waste at the end of its operating life. |
|  | Pay attention to the instructions<br>accompanied by this symbol.   |
|  | Fragile!   |
|  | CE-marked  |
|  | by this symbol   |
|  | Applied part type B.   |
| $\Omega$  | OHM (electric resistance unit)   |
| s   | Second (time unit)   |
| W   | Watt (power unit)  |
| Hz  | Hertz (frequency unit)   |
| mm  | Millimeter (length unit)   |
| A   | Amper (electric current unit)  |
| V   | Volt (electric voltage unit)   |
| Pa  | Pascal (atmosphere pressure unit)  |
| °C  | Degrees Centigrade (temperature unit)  |
| g   | Gram (weight unit)   |
| <b>SN</b>   | Serial number  |
| <b>LOT</b>  | Lot  |

## IX. FAQ

| Problems  | How to solve  |
|---|---|
| No response by button starting  | Fully discharged battery. Need to connect to charger to begin charge. If button continues not to responding - send the device to service.   |
| The button not response right   | Problem of button is entering of composite, bond, disinfectant, liquid, when device is used without barrier sleeves. Need cleaning with wet cotton with alcohol. If button continues not to response correctly - send the device to service.  |
| Can't start in HYPER mode, color ring blinking in red.  | HYPER mode has a 3 sec safety time where can't be started to prevent damage of soft tissues from overheating.   |
| Meaning of COLOUR RING colours  | Green - RAMP MODE- 1500mW/cm <sup>2</sup><br>Red - HYPER MODE - 3500 mW/cm <sup>2</sup><br>Blinking red - safety time<br>Blue - 100% charged battery<br>Slow blinking blue - in charging mode<br>Quick blinking blue - fully discharged battery<br>Yellow - LED module overheating  |
| Slow blinking of blue light of main LED   | This is the indicator that the battery only allows approx. 5-10 cycles until it is completely discharged. Urgent recharging need to required.   |
| Weak curing effect  | Need to clear output of the lightguide with non-metal stick an clear with alcohol for stuck composite or other. If the result is weak - send the device to service. Measuring intensity requires special equipment and a simple lightmeter is not suitable to check the intensity. If the problem is in LED module or electric control part - send back to service. |
| Patient feels discomfort during the curing process, especially during a long curing time 20-30 sec. | Cure at short cycles of 10 sec. or less with short pauses of 1-3 sec. sensibility is higher if light is near dental pulp an is decreased by every layer. The final curing can't make any sensitivity.   |
| Power cord or mains adapter is damaged  | Send back to service  |
| Entering some liquid into device  | Send back to service  |
| After drop on the floor device is not working properly  | Send back to service  |

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|  |   |   |
|--|---|---|
| <br>155, Vasil Aprilov blvd.<br>4027 Plovdiv, Bulgaria<br>office@bglight.com  | Декларация за съответствие<br>Declaration of conformity   | <b>ТД 11</b>  |
|  |   | Редакция <b>01</b><br>Лист <b>15</b> от <b>18</b>                                 |
| Manufacturer: <b>BG LIGHT LTD</b>  |   |  |
| Address: 155, Vasil Aprilov blvd., 4027 Plovdiv, Bulgaria<br>Tel.: +359 32 644089, Fax: +359 32 641913<br>BULSTAT 115841960, VAT N: BG115841960  |   |   |
| Product: <b>BLUEDED XPRESS – LED curing light</b>  |   |   |
| Classification: Active device in Class I of the Medical Device Directive 93/42/EEC+AC:2007   |   |   |
| Classification is done by the manufacturer according to MDD 93/42/EEC+AC:2007, Rule 12, Annex IX.  |   |   |
| Notified body: TUV NORD Polska Sp. z o.o., ul.Mickiewicza 29, 40-085 Katowice, Poland.   |   |   |
| This Declaration of conformity is valid only in combination with our certificates of Notified body TUV NORD Polska Sp. z o.o. Certificates N: AC090 100/1971/4047/2020, AC090 MD/1971/4047/2020, TNP/MDD/0334/4047/2020. |   |   |
| The manufacturer declares under its sole responsibility that the products are developed and produced in conformity with MDD 93/42/EEC+AC:2007 and the following applicable standards:                                    |   |   |
| EN 60601-1   | Medical electrical equipment - Part 1: General requirements for safety  |   |
| EN 60601-1-2   | Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests |   |
| EN 62304   | Medical device software. Software life cycle processes.   |   |
| EN 62353   | Medical electrical equipment. Recurrent test and test after repair of medical electrical equipment  |   |
| EN 62366   | Medical devices. Application of usability engineering to medical devices.   |   |
| EN ISO 14155-1   | Clinical investigation of medical devices for human subjects. General requirements.   |   |
| EN ISO 14155-2   | Clinical investigation of medical devices for human subjects. Clinical investigation plans.   |   |
| EN ISO 14971   | Medical devices – Application of risk management to medical devices   |   |
| EN ISO 15223-1   | Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements  |   |
| EN 1041  | Information supplied by the manufacturer with medical devices   |   |
| EN ISO 10993-1   | Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process  |   |
| 2012/19/EC   | Directive on waste electrical and electronic equipment (WEEE)   |   |
| All company products are manufactured under the current Quality Management System, ISO 9001 and ISO 13485.   |   |   |
| Dipl. Eng. Plamen Karaivanov<br>Manager<br>Bulgaria<br>BG LIGHT LTD  |    | 01.09.2020<br>Plovdiv,  |

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**Müller und Weygandt GmbH**
BG LIGHT LTD · 155, Vasil Aprilov blvd. · BGR - 4027 Plovdiv · Tel.: +359 32 644089