PONTIFORM AUTOMIX 10:1

Stand: 2/2018

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ⓑ Instructions for use

M+W Pontiform automix 10:1 is a pasty, multi-functional methacrylate-based composite material in a cartridge system that can be mixed at the time of use. The product is used for fabricating short and long-term temporary crowns and bridges, inlays, onlays and veneers directly in the patient's mouth. The strong final hardness and high elasticity make the material especially well-suited for fabricating large bridges. M+W Pontiform automix 10:1 fluoresces under UV light.

M+W Pontiform automix 10:1 is free of methyl methacrylate. The curing temperature is less than 40 °C. The completed temporary restoration protects the prepared teeth from outside effects and maintains the occlusion. After curing, the temporary restoration exhibits excellent flexural strength and abrasion resistance, low polymerisation shrinkage and outstanding fitting precision. The product is also distinguished by its good esthetics, optimum polishing qualities and high colour stability and fluorescence.

Temporary restorations fabricated with M+W Pontiform automix 10:1 can be easily repaired in case of breaks because the previously bonded material chemically bonds with the freshly mixed material or light-cured composites. The high stability becomes especially apparent on the thin edges of crowns. The material is therefore very well-suited for use with extensive bridgework.

Application

1. Preparation for impressions: Before preparing the die or the planned extraction, a situational impression is taken using addition cross-linked silicone impression materials (storage-stable models) or alginates. To improve the stability of the temporary restoration, cut the interdental areas into the impression.

Note: Level out undercuts in the silicone impression and provide any drainage grooves that may be needed.

2. Preparing the cartridges: Insert the cartridge into the appropriate dispenser, remove the cap and attach the mixing tip. The cap must not be re-used. From this point on, the used mixing tip serves as the cap.

Note: Discard some of the material that is first discharged from the mixing tip (about a pea-sized amount). The mixture will then be ideal. The same applies each time you

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perform mixing again.

3. Application: The substance is applied in the situation impression by lightly pressing it out of the mixing tip. To prevent bubbles from forming, always allow the mixing tip to be submerged in the material and fill from the bottom up.

4. Moulding the temporary restoration

a. Fill the situation impression in the necessary places with M+W Pontiform automix 10:1.

b. Place the cast in the patient's mouth. Steps a and b must be completed less than 50 seconds (processing time).

c. After the cast containing M+W Pontiform automix 10:1 has been in the mouth for 1-2 minutes, remove it from the mouth. (Alternative: 3-4 minutes when curing on the model).

Note: The temporary restoration should be in an elastic state when it is removed from the mouth.

5. Hardening and handling: After removing the plastic temporary from the situational impression (or from the die, if necessary), remove any excess material. The moulding should ideally be hardened in hot water (45°C - 55°C), for example in a polymerisation pressure pot. If this step is skipped, setting at room temperature requires approx. 6 minutes before the plastic temporary restoration can be processed and burnished with rotary instruments. **Do not inhale dust from grinding. Wear a face mask or use an exhaust system!**

Note: Before working on the temporary restorations, the inhibition layer on the surface which is the product of oxygen must be removed with a solvent (e.g. ethyl alcohol).

6. Attaching the temporary restoration: Temporary restorations should be fixed in place, preferably using eugenol-free temporary cements. If eugenol-containing temporary cements are used, it should be noted that using composite fixing cements at a later point can result in problems during the curing process.

7. Repairing the temporary restoration: M+W Pontiform automix 10:1 is characterised by its high mechanical stability. If a temporary restoration does break, however, we

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recommend the following procedure:

a) Temporary restoration breaks shortly after fabrication:

Reconnect the breakage by applying new M+W Pontiform automix 10:1 from the cartridge.

b) If a temporary restoration breaks while being worn by the patient:

The breakage point is roughened with a miller or sand blaster and an undercut is provided.

The breakage point prepared in this way is bonded using freshly mixed M+W Pontiform automix 10:1. To speed up the complete polymerisation, the repaired temporary restoration can be placed in water heated to 50°C for several minutes.

8. Special information: Material that has not set can be removed with cloth soaked in alcohol or similar solvents.

9. Warning: M+W Pontiform automix 10:1 is methyl methacrylate-free but does contain other acrylates, however. It is not possible to rule out sensitive patients experiencing additional sensitivity when using M+W Pontiform automix 10:1. Should allergic reactions occur, cease using M+W Pontiform automix 10:1. Do not use on patients who have allergic reactions to acrylates. Avoid contact with mucous membranes and eyes. If the material comes into

contact with the skin, wash off immediately with soap and water. If it comes into contact with the eyes, immediately flush with plenty of water and consult a doctor if necessary.

10. Storage information: Do not store above 25°C. Do not use after the expiry date.

Cartridge contents: 50 ml

Warranty


P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. user is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior Dental Materials GmbH product.

Limitation of liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep out of the reach of children! For dental use only!

Timeline for processing and curing of M+W Pontiform automix 10:1.

If curing is completed at room temperature, the subsequent processing of the temporary restoration can be performed after approx. 6 minutes (from the time of mixing).

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Manufacturer: P.L. Superior Dental Materials GmbH · 22607 Hamburg · Germany