
SEAL TEST

REFERENZ / REFERENCE / RÉFÉRENCE / RIFERIMENTO / REFERENCE / REFERENCIA / REFERENS / VIITEMATERIAALI

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GB SEAL TEST for routine inspection and functional assessment according to DIN EN ISO 11607-2:2006

A perfect sealed seam can only be guaranteed if the sealing parameters are set correctly and the material to be sealed has been adapted accordingly. If one of the process parameters (sealing temperature, throughput speed, contact pressure) shows deviations, then this is made visible on the indicator strip. According to DIN EN ISO 11607-2:2006, sealed seams must demonstrate the following quality properties as part of their functional assessment (OQ):

- Intact sealing over a defined seal width
- No punctures or tears
- No delamination or separation of materials
- No channel formation of open sealed seams

Place the M+W Seal Test into the packaging material with its test field being visible, and seal with the sealing device to be checked. The selected packaging material must have a width of at least 100mm. Following a visual inspection for a correct sealed seam over the entire width, the Seal Test is released and filed for documentation purposes. Optical magnification in the area of the intact sealed seam allows easy checking for consistency. In addition, there is the option of documenting the process data required for documentation on the seal test.

The M+W Seal Test offers the following fields for complete documentation:

- Sealing temperature
- Sealing pressure
- Sealing speed
- Sealing device
- Meter reading
- Date
- Release
- Signature

Should you sealed seam not be perfect (see reference photos), then we urgently recommend you to contact your trade partner. Even the best sterilisation becomes void without perfectly sealed seams!

With the indicator strip of the M+W Seal Test one can check whether these quality properties have been fulfilled. A subjective interpretation of the sealed seam is excluded, regardless of the type of sealing device used. **In addition, the standard DIN 58953-7 requires continuous monitoring of these quality properties and thus makes daily testing mandatory. Furthermore, the standard requires a routine inspection of the sealing devices.**

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