VOCO

Solobond Plus

Instructions for use

0482

Product description:
Solobond Plus is a light-curing two-stage universal bonding agent for a durable bond without marginal gaps between dental hard substance (dentine/enamel) and composite. The Solobond Plus system consists of a primer to condition the dentine and an adhesive as bonding agent between conditioned dentine, etched enamel and composite.

Solobond Plus tolerates residual humidity and permits "wet bonding". The special composition of Solobond Plus enables self-conditioning priming (*without etching* the dentine, see method of use A), as well as use of the total-etch-technique (*including etching* the dentine, see method of use B). Thus shallow cavities can be restored by using the total-etch-technique. For deep cavities in proximity of the pulp, self-conditioning priming is recommended.

Solobond Plus can be used with both light-curing and self-curing composite materials.

Indications:
Solobond Plus is suitable for direct composite restorations of all classes as well as for adhesive luting of crowns, bridges, inlays, onlays, veneers etc. with composite cements.

Method of use:

A) Direct composite restoration with self-conditioning priming

1. Drying
The application of Solobond Plus requires dry tooth surfaces due to the subsequent application of the composite. Avoid contamination with blood or saliva. The use of rubber dam is recommended.

2. Preparation
Prepare the cavity according to the principles of the adhesive filling technique. Do not remove the smear layer! Protect areas in proximity of the pulp with a suitable calcium hydroxide based liner.

3. Etching of enamel margins
Apply VocoCID etching gel to the areas to be etched and let act for ca. 20 - 30 s. Then remove gel by suction and carefully rinse with water. Remove residual humidity (air syringe).

4. Application of the primer
Dose primer onto a mixing pallet and apply to the moist dentine surface using a disposable brush. Rub primer into the dentine for 30 s.

Remove excess primer with a brush and thoroughly dry the cavity with oil-free air. Attention: Do not rinse with water!

5. Application of the adhesive
Repeat this step a second time when using Rebilda DC.

Dispense adhesive onto a mixing pallet and apply evenly to all prepared dentine and enamel surfaces using a new disposable brush. Thoroughly brush the adhesive into the primer for at least 15 s. Afterwards spread the adhesive gently with a faint air jet and polymerise with light for 20 s. To light-cure this material conventional polymerisation devices should be used. The light output should be a minimum of 500 mW/cm² on halogen polymerisation devices and 300 mW/cm² on LED devices. Do not remove or contaminate the unpolymerised inhibition layer as it is essential for bonding to the composite.

In order to produce a perfect bond between the hydrophilic dentine and hydrophobic composite, it is essential to rub the primer into the smear layer (see item 4) and thoroughly brush the adhesive into the primer/smear layer.

6. Composite filling

Composite filling material applies and adapt it to the cavity walls. Build up composite filling layer by layer, polymerise and finish.

B) Direct composite restoration with total-etch-technique

Steps 1. and 2., see A) 1. and A) 2.

3. Etch the entire cavity with VocoCID (etching gel); start with the enamel margins (enamel for ca. 30 s, dentine max. 15 s) remove the etching gel by suction, rinse for approx. 20 s and remove residual humidity (air syringe), but do not dry.

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In Case of linings containing composite materials or pulp capping materials, apply Solobond Plus intensively only to dentine areas (materials can be partly dissolved).

Avoid contact of Solobond Plus with the oral mucosa.

In case of contact of Solobond Plus with the eyes, rinse with plenty of water and seek medical advice.

Information on storage and application:

Tightly close Solobond Plus container after use. Store at temperatures of 4 °C to 23 °C.

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Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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