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## Intended use

Bite registration material based on vinyl polysiloxane

## Indications for use

StoneBite® is an addition-curing silicone for bite registrations.

## Indication

- Bite registrations
- Key for intraoral registration

## Warning

StoneBite® should be used solely on the occlusal surface of the teeth. Do not use StoneBite® for conventional impressions. Due to the high end hardness the removal is difficult and may lead to inadvertent extraction of teeth. Prior to application, ensure that undercuts, clasps, brackets, or retainers are blocked out. Material flowing into the undercuts of orthodontic or prosthetic appliances can significantly complicate removal and lead to mechanical damage.

## Contraindication

The product must not be used in cases of known allergy to platinum compounds.

## Side effects

There are no known side effects when used and processed properly. However, immune reactions (e.g. allergies) or local disorders of sensation (e.g. irritation of the oral mucosa) cannot be ruled out in principle. If in case doubt, it is advisable to carry out an allergy test. Professional users can contact the manufacturer to inquire about allergenic ingredients in the products.

## Target patient group

Patients who need dental care.

## Intended users

Dentist


## Product description

StoneBite® is an addition-curing silicone for bite registration with a high final hardness. The high final hardness enables a problem-free milling or grinding. In case of breaks, clearly defined pieces will be obtained.

## Processing 50 ml double cartridge

StoneBite® from the double cartridge is mixed automatically at a 1:1 ratio. To do this, insert the double cartridge into the injector. Remove the closure cap and express a small amount of material from the double cartridge to ensure that both components emerge evenly. Now, attach a mixing tip to the double cartridge. Check whether the material is being expressed from the tip evenly and in a mixed state. Using application tips, StoneBite® is to be applied cleanly and without waste onto the occlusal surface. Due to its

thixotropic properties, StoneBite® does not flow into the interdental spaces. Application must be carried out without interruption, as the oral temperature leads to accelerated curing. Achieve maximum intercuspation within the 30 s processing time and maintain closure for an additional 50 s. Once the total curing time has elapsed (80 s from start of mixing), the bite registration can be removed from the patient's mouth easily and without deformation. Please find processing times in "technical data" overview.

## Cleaning / disinfecting

Bite registration – After removal from the mouth rinse the bite registration under running water. Disinfection can be carried out with commercially available, suitable disinfection solutions (e.g. based on quaternary ammonium compounds with benzalkonium chloride or glutaraldehyde), according to the manufacturer's instructions. Mixing tip and application tip – After use leave the mixing tip on the double cartridge, disinfect it by wiping the mixing tip with a disinfectant and do not exchange it before the next application. It serves as a sealing cap. The mixing tip and application tip are for single use only and must be replaced before the next application. DS 50 Injektor – DS 50 Injektor can be used several times if properly disinfected. Disinfection can be carried out using a commercially available dipping disinfectant (e. g. glutaraldehyde) according to the manufacturer's instructions.

## Precautions

- The product is intended exclusively for its specified use and shall only be used by professionally qualified and trained personnel.
- Wear protective clothing! Set material is chemically resistant. Stains on your clothing cannot be removed.
- Due to production reasons, components may initially emerge unevenly from new cartridges. For this reason, express a small amount of material before attaching the mixing tip. Also, discard the first increment expressed after attaching the mixing tip.
- Attach a new mixing tip before each use of the cartridge. Using a mixing tip that contains partially or fully cured material may cause the cartridge to rupture, rendering the remaining material unusable.

## Batch number / use-by date

This product was developed for dental use and has to be applied according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore, the user is responsible for verifying the suitability and applicability of the product for the intended purposes prior to use.

- Do not use latex gloves, hand creams, cleaning agents, etc. during processing, as the chemical reaction of the two components is inhibited and the material no longer hardens completely. For this reason, we recommend that you wear commercially available polyethylene gloves.
- The setting times can be increased by lower temperatures or reduced by higher temperatures.
- Solutions used prior to bite registration (e.g. retraction solutions) may affect the curing reaction and must therefore be removed completely.
- Remove all residues of bite registration material from the patient's mouth.
- After use, leave the mixing tip on the double cartridge as a cap.
- Use suitable personal protective equipment (respiratory protection, safety glasses, and, if necessary, a suitable exhaust ventilation) during the finishing of the cured material.

## Keep away from children!

## Disposal

Disposal of contents/container in accordance with local/regional/national/international regulations and safety data sheet.

## Hazard identification

Hazard and safety notes can be found on the product label and the corresponding safety data sheet (Download at <https://dentamid.drewe.de/service/downloads>).

## Serious incidents

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

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