# Self Cure Activator

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<b>DIRECTIONS FOR USE - ENGLISH</b>
<p>Self Cure Activator (SCA) is for use with XP BOND, Prime&amp;Bond XP, Prime&amp;Bond NT, prime&amp;bond one etch&amp;rinse or Spectrum bond in conjunction with a dual cure or self cure resin cement such as Calibra<sup>®</sup> Esthetic Resin Cement or Core•X<sup>®</sup>-flow. SCA is suitable for bonding of all indirect restorations including metal, porcelain and composite crowns, inlays, onlays, veneers, bridge retainers and endodontic posts. Use of resin cements other than Calibra<sup>®</sup> Cement or Core•X<sup>®</sup>-flow or the above mentioned adhesives with SCA is at the discretion and sole responsibility of the dental practitioner.</p>

XP BOND, Prime&Bond XP, Prime&Bond NT, prime&bond one etch&rinse and Spectrum bond are self-priming dental etch&rinse adhesives designed for bonding resin-based restorative materials to enamel and dentine.

#### COMPOSITION

Urethane dimethacrylate (UDMA); 2-Hydroxyethyl methacrylate (HEMA); Catalyst; Photoinitiators; Stabilisers; Acetone; Water

#### INDICATIONS

When used with XP BOND, Prime&Bond XP, Prime&Bond NT, prime&bond one etch&rinse or Spectrum bond, SCA is indicated for adhesive cementation of indirect restorations in conjunction with dual cure or self cure resin cements.

#### CONTRAINDICATIONS

- SCA and XP BOND, Prime&Bond XP, Prime&Bond NT, prime&bond one etch&rinse and Spectrum bond are contraindicated for:
  - direct or indirect pulp-capping,
  - use in patients known to be allergic to methacrylate resins or any of the other components.

#### WARNINGS

- SCA contains polymerisable methacrylate monomers. Avoid prolonged or repeated contact with skin (allergic contact dermatitis), oral soft tissues, and eyes. Avoid prolonged inhalation. Do not take internally.
- Eye Contact: SCA contains methacrylates which may be irritating to eyes. Before using this product wear protective glasses and cover the patient's eyes to protect from splashing material. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Skin Contact: SCA contains polymerisable monomers which can cause skin sensitisation (allergic contact dermatitis) in susceptible individuals. In case of contact with skin, immediately wipe off thoroughly with cotton and alcohol and then wash well with soap and water. If skin rash and sensitisation or other allergic reactions occur, discontinue use and seek medical attention.
- Oral mucosa contact: Avoid contact with oral soft tissues. If accidental contact occurs, flush mucosa with plenty of water and expectorate water. If sensitisation of mucosa persists, seek medical attention immediately.
- Ingestion: Do not swallow or take internally. If accidental swallowing occurs, drink lots of water. If nausea or illness develops, seek medical attention immediately.
- SCA contains acetone. Do not inhale vapor. Inhalation: Give oxygen or artificial respiration if necessary.

#### PRECAUTIONS

- This product is only intended for use as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the dental practitioner.
- SCA is for use in conjunction with XP BOND<sup>®</sup>, Prime&Bond<sup>®</sup> XP, Prime&Bond<sup>®</sup> NT<sup>™</sup>, prime&bond<sup>®</sup> one etch&rinse or Spectrum<sup>®</sup> bond and DeTrey Conditioner 36. For correct use, please also refer to the Directions for Use of the respective products.
- Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
- Dental materials containing eugenol should not be used in conjunction with this product because they may impede curing and cause softening of the polymeric components of the material.
- Insufficient data exist to support the use of desensitising agents and/or cavity cleansing agents with SCA. Recommended pre-treatment is pumice or prophylaxis paste with a rubber cup.
- Contact with saliva, blood and/or some astringent solutions during adhesive procedures may cause failure of the restoration. Use of a rubber dam or adequate isolation is recommended.
- Exercise caution if mineral-impregnated (e.g., ferric compounds) retraction cords and/or hemostatic solutions are used in conjunction with adhesive procedures. In-vitro studies have suggested compromised adhesive bond strength to enamel and dentine contaminated with haemostatic agents. Marginal seal may be adversely affected, allowing microleakage, subsurface staining and/ or restoration failure. If gingival retraction is necessary, use of plain, non-impregnated cord is recommended.
- In-vitro studies have shown a relationship between adhesive performance and air-thinning/drying/ solvent evaporation techniques. Improper over-thinning or under-evaporating may compromise bond strength, leading to microleakage, post-operative sensitivity and/or restoration failure. (See Step-by-Step Instructions)
- The use of a dual-cure adhesive system, such as XP BOND<sup>®</sup>, Prime&Bond<sup>®</sup> XP, Prime&Bond<sup>®</sup> NT<sup>™</sup>, prime&bond<sup>®</sup> one etch&rinse or Spectrum<sup>®</sup> bond, with SCA can shorten the working time of a dual cure resin cement system. This effect should be investigated in the laboratory prior to clinical use.
- The SCA bottle should be tightly closed immediately after use.
- When mixed with SCA, XP BOND<sup>®</sup>, Prime&Bond<sup>®</sup> XP, Prime&Bond<sup>®</sup> NT<sup>™</sup>, prime&bond<sup>®</sup> one etch&rinse or Spectrum<sup>®</sup> bond can still be cured with visible light. Proceed immediately once materials have been placed in mixing well or protect from ambient light.
- Must only be used in well ventilated areas.
- Flammable:** SCA contains acetone. Keep away from sources of ignition.
- STORAGE: SCA should be kept out of direct sunlight and stored in a well ventilated place at room temperature not exceeding 25°C. Refrigerated storage is not required, but is acceptable when not in use. Inadequate storage conditions will shorten the shelf life and may lead to the product malfunctioning. Allow material to reach room temperature prior to use. Protect from moisture. Do not freeze. Do not use after expiration date.

#### ADVERSE REACTIONS

- Product may irritate the eyes and skin. Eye contact: irritation and possible corneal damage. Skin contact: irritation or possible allergic response. Reddish rashes may be seen on the skin. Mucous membranes: inflammation, oedema, sloughing. (See Warnings)