Superlite Power Pen/ Color Pen

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Superlite Power Pen/ Color Pen

Instructions for use

Thank you for trusting our products!

Carefully read this Operating instructions manual before installing and operating the unit to use and maintain it successfully!

Save this Operations guide for your reference in the future.

CONTENTS

I. UNIT DESCRIPTION AND FUNCTION

II. SYMBOLS

III. SAFETY PRECAUTIONS

IV. TECHNICAL DATA

V. COMPLETE SET

VI. PREPARATION AND SEQUENCE OF OPERATION

VII. ROUTINE CARE AND MAINTENANCE

VIII. WARRANTY CONDITIONS

IX. SERVICE DATA

X. FAQ

XI. DECLARATION OF CONFORMITY

XII. M+W DATA

I. UNIT DESCRIPTION AND FUNCTION

Dental LED curing light M+W Superlite Power Pen/ Color Pen is a medical device - specialized light source for intraoral polymerization of dental materials, sensitive to the blue part of the light spectrum.

The device is intended for use only by a qualified dental practitioner and in a dental practice.

M+W Superlite Power Pen/ Colour Pen consists of a cordless handpiece and charging adapter. The Superlite Power Pen also includes a holder.

M+W Superlite Power Pen/ Color Pen is manufactured in conformity with the requirements of Regulation on medical devices MDR 2017/745 and standards ISO 13485:2016, ISO 9001:2015.

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 BG LIGHT PRO LTD · 155, Vasil Aprilov blvd. · 4027 Plovdiv/Bulgaria

II. SYMBOLS

	Caution! Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Consult instructions for use Indicates the need for the user to consult the instructions for use.
	Dangerous voltage To indicate hazards arising dangerous voltage.
	Hazardous light emission To indicate hazards arising from light radiation.
	Hazardous thermal effects To indicate hazards arising from thermal effects.
	Manufacturer Indicates the medical device manufacturer.
	Date of manufacture Indicates the date when the medical device was manufactured
	Medical Device Indicates the item is a medical device
	Unique Device Identifier Indicates a carrier that contains Unique Device Identifier information
	Catalogue number Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Serial number Indicates the manufacturer's serial number so that a specific medical device can be identified.
	Batch code Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Applied part type B according to electric safety classification
	Temperature limit Indicates the temperature limits to which the medical device can be safely exposed.
	Humidity limitation Indicates the range of humidity to which the medical device can be safely exposed.
	Waste Electrical and Electronic Equipment (WEEE) According to Directive 2012/19/EU, this symbol indicates that the product should not be disposed of as urban waste at the end of its operating life.
	Fragile Indicates a medical device that can be broken or damaged if not handled carefully.
	European Conformity Indicates conformity with local laws and regulations within the European Economic Area

III. SAFETY PRECAUTIONS

GENERAL WARNINGS

M+W Superlite Power Pen/ Color Pen is a Class I medical device and it meets the requirements of the Medical Devices Regulation - MDR (EU) 2017/745. In order to be used safely for staff and patients, the following rules must be observed:

- Do not allow unauthorized and untrained personnel to use the device to avoid risks.
- Disconnect the device from the mains after completing the procedures.
- Do not use or store the device in a dusty environment.
- Do not expose the device to direct sunlight.
- Do not spray disinfectant directly into the device - only rubbing with a swab drained of disinfectant is acceptable.
- Do not get wet or drop liquid on the device, cables, adapter to avoid electric shock or damage to the device.
- Store the device in a dry place, moisture can cause electric shock and damage.
- The device must not be used if any of its parameters are not normal (timer, light intensity, heat radiation).
- Strong electromagnetic fields in the building can cause interference and malfunction of the device. If their source cannot be determined, change the location of the device and plug it into another socket or other room, even in another building.
- Opening and repairing the appliance may only be carried out by authorized service technicians from the manufacturer.
- Only original M+W Superlite Power Pen/ Color Pen parts must be used when replacing defective parts. The device or any of its parts must not be disassembled while it is connected to the mains!

- Fragile! Use caution when transporting, using and storing the device!

- According to Directive 2012/19/EEC, this symbol indicates that the product should not be disposed of as a general waste at the end of its lifespan. The product must be taken to a specialized center for the separate collection of electrical and electronic equipment according to local regulations. Proper disposal of equipment that is no longer used prevents negative consequences for the environment and human health!

- In accordance with the requirements of MDR (EU) 2017/745, user and / or the patient must report any serious accident that have occurred during use of the device to the manufacturer and the competent authority of the Member State in which the user/patient is established.
- All packaging materials of the product must be kept away from children to avoid risks of injury / suffocation.

SAFETY MEASURES AND RISKS

The device must be used in strict accordance with the Operating Instructions Manual.

1. Electrical safety

Before starting the appliance, make sure that the voltage and the type of plug correspond to the mains supply in the country. Use only the original adapter type SYS1638-0605-W2E.

Electrical safety is ensured by class II protection against electric shock according to EN 60601-1.

M+W Superlite Power Pen/ Color Pen must only be operated indoors, under the following conditions:

- temperature from + 10 ° to + 40 ° C;
- relative humidity 30 - 75%;
- lack of dust in the room;
- atmospheric pressure 700 - 1060 hPa;
- absence of chemically active and flammable substances;

- no part of the device should be wetted or immersed in water;
- the device or any of its parts must not be disassembled while it is connected to the mains!

Protect the cables of the appliance from insulation damage and breakage from sharp objects, strong pulling, rodents, chemicals. If such damage is noticed on the electrical cables, it is necessary to take the device for repair to the company service. The device must not be used with damaged cables.

In case of thunderstorms, the procedures must be stopped and the adapter must be disconnected from the mains.

Risk: Failure to comply with these instructions may result in electric shock to users of the device.

2. Light radiation

M+W Superlite Power Pen/ Color Pen is a source of intense light in the blue range, to which the human eye has a high sensitivity. This results in serious measures to be taken for patients, medical staff and accidentally nearby people, animals and plants.

As such, use protection goggles for the operator.

Irradiation of the eyes and skin with intense light carries a risk of damage from light and heat.

The light should never be directed at the eyes! Irradiation should be limited to the workplace area. The special safety goggles that meet the requirements must be used:

- to cover the eyes and temples tightly, even if the person is wearing optical glasses.
- be made of volumetric colored impact-resistant plastic.
- do not transmit light with a wavelength of 410 to 490 nm.
- reduce the intensity of the blue spectrum by more than 100 times.
- have a stable mechanical structure, no scratches, cracks and damage to its surface.

The device can be used only after a doctor's consultation on or by persons suffering from photo-biological reactions; persons taking photosensitive drugs; persons undergoing cataract surgery, persons with retinal diseases, etc.

The risk of improper irradiation is severe eye irritation, temporary spots in the visual field, severe visual impairment in direct radiation, to loss of vision.

3. Thermal radiation

The thermal effect is due to the absorption of the energy of the blue light in the tissues, during which the energy is converted into heat. The risk is only with prolonged overdose.

Risk of pain, burning of soft tissues.

4. Fire safety

- Keep the device away from solvents, flammable liquids and powerful heat sources.
- Do not expose to direct sunlight.
- Do not allow liquids and detergents to enter the device, as this may cause a short circuit and fire or cause potentially dangerous damage.
- If the product emits an odor or smoke – disconnect from the mains, do not attempt to repair it, take it to a service center.

Risk of fire, ignition and damage.

5. Contraindications

- The device can be used only after medical consultation on or by: persons with implanted cardiac pacemaker; persons suffering from photobiological reactions; persons taking photosensitive drugs; persons undergoing cataract surgery; persons with retinal diseases; people with allergies; people who have recently undergone cosmetic surgery on the face or lips, including injections of hyaluronic acid or botox; people with very sensitive skin or dermatitis, etc. If you are taking photosensitizers or medicines, check the package leaflet.


IV. TECHNICAL DATA

1. Operating voltage - 3.7V Li Ion battery type 14500 (800mAh capacity)
2. Power consumption - mains adaptor: SYS1638-0605-W2E
 - 110-240VAC / 50-60Hz input voltage
 - Output: 5VDC / 1.2A micro USB socket
 - 0.5A max power consumption
3. Dimensions of cordless handpiece:
 - min. diameter – 18mm, max diameter - 25mm
 - length - 148 mm (without lightguide)
4. Weight of cordless handpiece (with lightguide and protection filter)
 - M+W Superlite Power Pen (plastic handpiece) - 80g
 - M+W Superlite Color Pen (metal handpiece) - 100 g
 - Holder – 276g
5. Polymerization mode – with soft start.
6. Light source – 2 band LED module with reflector optics.
7. Emitted light – blue, visible spectrum, 410-490 nm.
8. Light intensity - up to 1100 mW/sq.cm Intensity is constant and not dependent on battery charge level.
9. Emitting time: - 10 / 20 / 30 sec. / $\pm 10\%$
10. Working ability after fully charged battery - up to 2100sec. (210 cycles x10sec).
11. Charge time - 3 hours for full charge and 2 hours to 80% charge.
12. Visual indication - color indication ring and audible signal every 10 sec.
13. Low battery warning. Fully discharged battery indication.
14. Fiberoptic lightguide (rotatable 360 degrees)
15. Working mode (1 min. work / 10 min. pause)
16. Degree of protection against electric current - applied part type B.

The manufacturer of this unit declares to provide on request all additional necessary technical documentation / information which will help the user's technical staff to service the parts of the unit which the manufacturer has claimed to be a subject for repair.

V. COMPLETE SET

1. Charge adaptor with micro USB output 5V DC / 1.2A - 1 pc.
2. Cordless handpiece - 1 pc.
3. Battery pack - 1 pc.
4. Lightguide - 1 pc.
5. Protection shield - 1 pc.
6. Protection cap - 1 pc.
7. Holder - 1pc.
(M+W Superlite Power Pen)
8. Operations guide - 1 pc.

Control and indication:

1. USB charging connector
2. Start / Stop button
3. Illuminated indicator ring
4. Battery pack

VI. PREPARATION AND SEQUENCE OF OPERATION


M+W Superlite Power Pen/ Color Pen is designed for use as a battery powered unit.

1. Take the curing light out of the packing and assemble the holder by inserting the upper part into the base part (as shown on picture below). Place the holder on a horizontal surface. The holder is only included with the Superlite Power Pen. For the Superlite Color Pen, it is optionally available in the colours black, white and silver.

2. Mount lightguide and protection shield or protection cap.

See VII. ROUTINE CARE AND MAINTENANCE for instructions of initial cleaning and mounting.

For polymerization of posterior teeth
(attach as shown)
(1) Protection cap

For polymerization of anterior teeth / general application

Connect the power adapter to charge battery if needed.

The device can be operated during charging.

Unplug the charger after charging the device or at the end of the day.

3. Indication

3.1. Color indication ring

Deactivated mode - Color indication ring is not lit.

Work mode – polymerization.

Battery charge / sound signal.

Fully charged battery – device can be unplugged from adapter.

Overheating – **Thermal protection activated**, wait until device cools down.

Battery is completely discharged - plug into the charging adapter.

3.2 LED-Modul

LED module indication:

Blinking means that the battery is discharged - and the device will soon turn off.

After the initial press of the button, the device is in standby mode. Select the timer:

Press the start and stop button once: 10 seconds of light polymerisation

Press the start and stop button twice: 20 seconds of light polymerisation

Press the start and stop button 3 times: 30 seconds of light polymerisation


After the set time has elapsed, the unit will turn off the output light.

If desired, the dentist can interrupt the cycle by pressing the Start / Stop button again.

During operation, the fiberoptic lightguide must be at a distance of at least 2 mm from the obturation

Diagram of curing mode

5. Overheating protection

The unit is equipped with overheating protection, which is activated if the temperature of the LED module rises up to 50°C. When overheating protection is activated the color ring will blink in yellow. It is necessary to wait for device's cooling down.

6. Battery – maintenance, charging and replacement

The device uses Li-Ion battery type 14500 / 3,7V / 800 mAh.

Charging the battery - by connecting the charging adapter to the mains.

Attention! Connect properly to the socket!

Unplug the adapter after charging.

Charge adaptor is SYS1638-0605-W2E - 100-240VAC input and 5VDC output - micro USB type.

It is not allowed to use another type power supply.

When to replace the battery

- at reduced capacity - frequent charging is required;
- in case of battery failure.
- in case of damaged charging connector

To replace a Li-Ion battery

- Disconnect the adapter from the socket and the USB socket.
- Place the handpiece on a non-slip surface and unscrew Part 3. Release Part 2.

Carefully remove the connector (4).



- The installation of a new battery is done in the reverse order.

- When assembling the handpiece - place Part 2 (containing the battery) on a non-slip surface and press Part 1 to Part 2 at the same time. At the same time, screw the cylindrical Part 3 back to the front conical Part 1.

Note: The parts of the assembly between 1 and 2 must fit!

Do not twist the battery cables!

Use only original Li-Ion batteries by the device manufacturer (14500 / 3,7V / 800 mAh)!

The battery must not be short circuited, mechanically damaged, immersed in water or to be heated or thrown into a fire.

- Non-functioning batteries can be taken to battery collection points.
- When disposing of the battery, observe the environmental protection regulations!

VII. ROUTINE CARE AND MAINTENANCE

1. Examining of lightguide / protection cap.

The light guide is fitted by pushing it on and removed by pulling it off (can be rotated 360 degrees).

Every day the lightguide tip should be examined for:

- stuck restorative material (it should be removed very carefully with blunt non-metal object);
- mechanical damages on the lightguide.

Remove the light guide and looking through it at strong light: if there are darkened areas over 10% of the area, replace it with a new one.

Before every patient lightguide should be sterilized by autoclave at 134°C for 5 min.

The same must be done for protection cap.

The protective filter is wiped only with disinfectant.

Note: Lightguide and protection cap are delivered non-sterile.


2. Cleaning of the unit.

For disinfection of the unit, spray the disinfection agent onto a piece of soft cloth / cotton and clean the handpiece and power adapter.

Do not use abrasives or solvents as these may damage the plastic parts of the unit!

Do not spray directly onto handpiece or in button / USB charge connector!

VIII. WARRANTY CONDITIONS

1. Warranty period of M+W Superlite Power Pen/ Color Pen dental curing light is 24 (twenty four) months from the date of purchase. If date of purchase is not written, warranty starts at date of production.
2. During the warranty period the replacement of the defective parts is done free of charge by manufacturer's service. Lightguide is not covered by warranty. LED module and battery have 6 months warranty.
3. If during the guarantee period the unit is damaged, due to incorrect operation (mechanical, chemical, thermal or electrical damages), the repair is not under warranty.
If the charging cable is damaged, it must be replaced. The device must not be charged or operated with a damaged cable.
4. No damages or complaints are accepted as a result of electric shocks, thunderstorms, non-compliance with electrical safety measures or insufficient protection of patients, personnel and other persons from light radiation.
5. Claims due to improper or insufficient security and care, protection and security during transport, unpacking, relocation, operation and storage of the device are not accepted. The warranty is void for the above events.
6. The device must only be transported to manufacturer's service only in original packaging in order to avoid unwanted damage. The light guide must be packed separately from the cordless handpiece. The warranty card must be presented (this Guide, with completed Chapter XII M+W data) or a copy of the sales invoice showing the serial number of the device.
7. When carrying out repairs by unauthorized persons outside the company service or if non-original elements / parts are used, the user loses the right to warranty service.
8. Disclaimer: The manufacturer is not liable for compensation and lost profits in the event of damage or incomplete functioning of the device, regardless of the cause or in the event of problems arising from failure to comply with safety measures in the event of risks (chapter III, Safety measures and risks).
9. In case of disputes in connection with the application and interpretation of this Operating instruction manual, they will be decided by the court in Plovdiv, by virtue of the current Bulgarian legislation.
10. The warranty service is performed in the company service at the address:

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11. The device's service life is set at 5 years from the date of purchase. After the end of its life cycle, the user must discontinue use of the device.

At the user's request, the life cycle can be extended. For this purpose, the device must be returned to the manufacturer, where its characteristics must be checked and all main elements (emitters, adapter, electronic boards) tested and replaced.

IX. SERVICE DATA

X. FAQ

Problem	Solution
No response by button starting	Fully discharged battery. Need to connect to charger to begin charge. Use only the original adapter type SYS1638-0605-W2E. If button continues not to respond - send the device to service.
The button does not respond properly	You need to clean with a cotton swab soaked in alcohol. If the button still does not respond correctly - send the device to service
Weak curing effect	It is necessary to clean the light exit window of the light guide with a suitable non-metallic object and to remove any adhering composite (restorative material) or other impurities with alcohol. If the result is poor - send the device to the service.
The patient feels discomfort during curing process especially during long curing time 20-30sec	Reduce the exposure time to 10 seconds and take a short break before the next exposure. The sensitivity is increased, especially for exposures in the near-pulse range.
Power cord or mains adaptor is damaged	Replace power supply adapter.
Entering some liquid into device	Send to service.
After drop on the floor device is not working properly	Send to service.

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 TECHNICAL FILE M+W Superlite LED curing light	EU Declaration of conformity <i>Developed in conformity with MDR (EU) 2017/745</i>	TD 7.2
		Revision 03
		Page 1 of 2

Manufacturer: **BG LIGHT PRO LTD**
 SRN:
 Address: 155, Vasil Aprilov blvd., 4027 Plovdiv, Bulgaria
 Tel.: +359 32 644089, +359 888 809256, email: office@bglight.com
 BULSTAT UIC 208000520, VAT N: BG208000520

Dental LED curing lights			
N	Product code BGL	Product code M+W	Name of device
1	200-003.01MW	200027	M+W Superlite Power Pen cordless (black)
2	200-003.02MW	200025	M+W Superlite Power Pen cordless (silver)
3	200-003.03MW	200026	M+W Superlite Power Pen cordless (white)
4	200-003p.02MW	200006	M+W Superlite Color Pen (dark grey)
5	200-003p.04MW	200000	M+W Superlite Color Pen (white)
6	200-003p.05MW	200001	M+W Superlite Color Pen (orange)
7	200-003p.06MW	200002	M+W Superlite Color Pen (light green)
8	200-003p.09MW	200003	M+W Superlite Color Pen (blue)
9	200-003p.10MW	200004	M+W Superlite Color Pen (violet)
10	200-003p.12MW	200005	M+W Superlite Color Pen (pink)
11	200-004MW	100045	M+W Superlite 1300
12	200-004h MW	1052297	Handpiece for M+W Superlite 1300
13	200-009ort.01MW	200089	M+W Superlite LED (black) (XPRESS ortho)
14	200-009ort.03MW	200090	M+W Superlite LED (white) (XPRESS ortho)

Basic UDI: 38001748200000G8
 EMDN code: Q0190

Classification: Active medical device of **Class I** of the Regulation on medical devices - MDR (EU) 2017/745

Intended purpose: M+W Superlite is designed for photopolymerization of composites and materials used in dental practice (irradiation of blue light 410-490 nm).

The manufacturer declares under its own responsibility that the specified medical device complies with the applicable GENERAL SAFETY AND PERFORMANCE REQUIREMENTS, defined in Annex I of the normative act described below and normative technical documents, when used for its intended purpose and in accordance with the safety requirements.

Document	Title	Edition / date of issue
Regulation (EU) 2017/745	REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL <i>of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC</i>	05.05.2017 (last change 24.04.2020)

To achieve compliance, the requirements of the following standards are met:

EN ISO 13485:2016 +/AC:2017/ /AC:2018/ A11:2022 +/AC:2017/ /AC:2018/ A11:2022	Medical devices - Quality management systems - Requirements for regulatory purposes
EN ISO 9001:2015	Quality management systems - Requirements
EN ISO 60601-1:2006 /A1:2013/AC:2014/A2:2022	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.
EN 60601-1-2:2015/A1:2021	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests
EN 60601-1-6:2010+ /A1:2015 / /A2:2021	Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
EN 60601-1-8:2007+ /A1:2013 /A11:2017 /A2:2021	Medical electrical equipment — Part 1-8: General requirements for basic safety and essential performance — Collateral standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems

EN ISO 10650:2018	Dentistry - Powered polymerization activators
EN 62304:2006/A1:2015	Medical device software. Software life cycle processes.
EN 62353:2014	Medical electrical equipment. Recurrent test and test after repair of medical electrical equipment
EN 62366-1:2015+ AC:2016/A1:2020	Medical devices. Application of usability engineering to medical devices.
EN ISO 14155:2020	Clinical investigation of medical devices for human subjects - Good clinical practice
EN ISO 10993-1:2018	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
EN ISO 14971:2019+A11:2022	Medical devices – Application of risk management to medical devices.
CEN ISO/TR 24971:2020	Medical devices - Guidance on the application of ISO 14971 (ISO/TR 24971:2020)
EN ISO 15223-1:2021	Medical devices - Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements
EN ISO 20417:2021	Medical devices - Information to be supplied by the manufacturer
Directive 2012/19/EC	Directive on waste electrical and electronic equipment (WEEE)

Classification is done by the manufacturer according to Regulation on medical devices - MDR (EU) 2017/745, Annex VIII, Rule 13. Conformity assessment procedure according to article 52, paragraph 7 of MDR (EU) 2017/745.

The declaration of conformity is issued in implementation of Annex IV "EU Declaration of conformity" of EU Regulation 2017/745, based on the results of tests carried out and assessment of compliance with the General safety and performance requirements defined in Annex I, implemented and certified Quality Management System - ISO 9001:2015 certificate – No AC090 100/2409/6029/2025, ISO 13485:2016 – No AC090 MD/0008/6029/2025 from TUV NORD Polska Sp. z o.o. (NB 2274).

BG LIGHT PRO LTD maintains data on the provision, evaluation and maintenance of compliance of the medical device, according to the requirements of Annex II "Technical documentation" of MDR (EU) 2017/745.

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**Müller & Weygandt GmbH · Reichardsweide 40