# PERMABOND

## Gebrauchsinformation

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## GB Instructions for use

### Intended purpose

M+W Permabond is a light-curing, single-component bonding agent for dentine and enamel for use with light-curing composites. Before applying M+W Permabond, enamel and dentine surfaces must be etched.

### Indication

Light-curing bonding agent for enamel and dentine, used in conjunction with the total-etch technique for light-curing composites.

### Recommended use

**Note:** When using M+W Permabond, it is essential that both the enamel and the dentine surfaces are etched beforehand!

1. Prepare the cavity in accordance with the general rules for adhesive techniques. Bevel enamel margins and protect areas close to the pulp in deep cavities using an appropriate base lining such as calcium hydroxide lining material. It is advisable to use a rubber dam.

2. Etch dentine and enamel with 37% phosphoric acid gel (e.g.

M+W Ätzel). Observe the manufacturer's instructions for the corresponding etching gel. The dentine must not be exposed to the etching gel for more than 15 seconds. The enamel should have a chalky white appearance after the etching process.

3. Rinse the etched surface for at least 15 seconds with water, then dry the surface with oil and water-free air. Make sure that the dentine does not over-dry.

4. Apply M+W Permabond onto the etched surface and work in for 20 seconds. Briefly dry for approx 5 seconds with oil and water-free air. The treated surface should be covered by an even, glassy coating. If this is not the case, a little more material must be applied.

5. Light-cure the M+W Permabond using a suitable light-curing unit for 10 seconds.

6. Use the restorative material in according with the instructions of the manufacturer in question.

### Please note

• Light units should have a luminous intensity of 400 mW/cm<sup>2</sup>, emit at 450 nm, and their performance should be checked regularly. The light should be positioned as closely to the material as possible.

• The etching gel must be used on both enamel and dentine surfaces. DO NOT use the material without etching the dentine.

• The oxygen inhibition layer that develops during use facilitates bonding with the materials used over it, and therefore must not be removed.

• Contamination of the bottle's contents with a contaminated disposable brush, and premature polymerization due to ambient light, should be avoided.

• Where possible, a matrix should only be used after application of the bonding agent in order to prevent excess material collecting between the tooth and the matrix (pooling effect).

• Do not use chemically curing and dual-curing composites with M+W Permabond.

• Avoid skin contact with the fluid. In case of accidental contact, rinse with water and soap. In case of contact with the eyes, rinse immediately with lots of water and consult a physician.

### Interactions

M+W Permabond should not be used in conjunction with materials that contain eugenol, as these can inhibit polymerization.

### Contra-indications

The use of M+W Permabond is contraindicated if a dry working field cannot be ensured during use, where there is exposed pulp or where there is an existing allergy to any of the components of the material.

### Side effects

There are no known side effects to date. In some cases, hypersensitivity or a contact allergy has been described in the case of products with similar constituents.

### Composition

Bis-GMA-based resin matrix, water, additives, catalysts.

### Storage

Store under dry conditions at 2 - 25°C (35,6 - 77°F). Do not use after the expiration date.

**Keep out of the reach of children! For dental use only!**

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