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Ⓧ Instruction for use

Intended Purpose and Application

M+W Bloc Out is a light-curing blocking-out material. The material is used for the modelling on different materials, e. g. plaster or acrylic. The material remains stable during deep-drawing.

1. Fabricate the stone model and trim. Remove plaster bubbles.
2. After the model is completely dry, apply M+W Bloc Out. This is applied to the labial

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surfaces of the teeth to be whitened and is used later as a reservoir for the gel. Usually this involves teeth 15-25 and 35-45. It is of course possible to prepare a smaller or larger number of teeth or single teeth.

Important: Apply M+W Bloc Out in layers approx. 0,5 mm thick and maintain a gap of approx. 0,5-1 mm to the gingival margin. Do not cover incisal edges and occlusal surfaces with gel!

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3. Polymerise M+W Bloc Out. A standard halogen polymerisation lamp (approx. 20 sec. per tooth), a plasma light-curing unit (approx. 3-6 sec. 2-3 min., depending on the manufacturer) or a light oven (approx. 2-3 min., depending on the manufacturer) can be used for polymerising.

4. After curing, completely remove the inhibition layer on the surface using a cloth moistened with alcohol. Do not use acetone or similar cleaning agents.

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Müller & Weygandt GmbH · Reichardsweide 40 · 63654 Büdingen/Germany · Tel.: +49 (0) 60 42 - 88 00 88