# 3M Impregum™ Pentax™ Super Quick

- Ⓢ Polyether Impression Material
- Ⓢ Polyether Abformmasse
- Ⓢ Matériau d'empreinte polyéther
- Ⓢ Material e impronta a base di polieter
- Ⓢ Matériau de impressão de poliéster

Instructions for Use  
Gebrauchsinformation  
Mode d'emploi

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3M, ESPE, Garant, Impregum, Pentax and Pentamix are trademarks of 3M or 3M Deutschland GmbH. © 2016, 3M. All rights reserved.

## ENGLISH

**Product Description**  
Impregum™ Pentax™ Super Quick – Heavy Body and Medium Body are fast-setting polyether impression materials with a heavy-bodied and medium-bodied consistency. The mixing ratio, based on volume, is 5 parts catalyst paste : 1 part catalyst paste. The paste is mixed in Pentamix versions from Pentamix™ 2 (September 1999); on earlier Pentamix units (before September 1999) extrude the material too slowly. Impregum™ Super Quick – Light Body is a fast-setting polyether impression material with a light-bodied consistency in the Garant™ cartridge. The mixing ratio, based on volume, is 2 parts base paste : 1 part catalyst paste.

Impregum Pentax Super Quick and Impregum Super Quick materials are not available in all countries.

## Symbol Glossary

Reference Number & Symbol*	Symbol	Description of Symbol
5.1.1 Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
5.1.4 Use-by date		Indicates the date after which the medical device is not to be used.
5.1.5 Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
5.1.6 Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified.
5.3.7 Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed.
5.4.2 Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure.
5.4.4 Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for variety reasons, be presented on the medical device itself.

\* Symbols from ISO 15223-1 Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General.

Ⓢ These Instructions for Use must be kept for reference for the duration of product use. See the pertinent Instructions for Use for details on Astringent Retraction Paste, Intra-oral Syringe, the Pentamix™ mixing unit and accessories, Polyether Adhesive, Pentax™ mixing unit and accessories, Polyether Adhesive, Pentax™ Mixing Tip – Red, Garant™ Dispenser, Garant™ Mixing Tips and Intra-oral Tips. The product may only be used when the product labeling is clearly readable.

**Indications**  
• Impressions of inlay, onlay, veneer, crown, and bridge preparations  
• Fixation and implant impressions

Impregum Pentax Super Quick – Heavy Body and Medium Body are fast-setting polyether impression materials with a heavy-bodied and medium-bodied consistency. The mixing ratio, based on volume, is 5 parts catalyst paste : 1 part catalyst paste. The paste is mixed in Pentamix versions from Pentamix™ 2 (September 1999); on earlier Pentamix units (before September 1999) extrude the material too slowly. Impregum™ Super Quick – Light Body is a fast-setting polyether impression material with a light-bodied consistency in the Garant™ cartridge. The mixing ratio, based on volume, is 2 parts base paste : 1 part catalyst paste.

**Precautionary Measures**  
**For Patients**  
This product contains substances that may cause an allergic reaction by skin contact in certain individuals. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If an allergic reaction is observed in a patient, discontinue use of the product in this patient and seek medical attention as needed. Known allergens are 1-Dodecylimidazole, sulfonium salt, and mint flavors.

**For Dental Personnel**  
This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended.  
3M MSDSs can be obtained from www.mm.com, or contact your local subsidiary.

**Preparatory Steps**  
**Blocking Out**  
Block out any undercuts or interdental spaces sufficiently to aid impression removal from the teeth after the impression material has set. Removal from the mouth may be difficult or not least the extraction of natural teeth or prostheses. Blocking-out is especially important when using impression materials with high Shore hardness.

**Prepare Impression Tray**  
All impression trays generally used for precision impressions are suitable for use with Impregum Pentax Super Quick.  
For sufficient adhesion, apply a thin layer of Polyether Adhesive to the tray and allow to dry completely according to the Polyether Adhesive Instructions for Use.

**Prepare Pentamix/Pentax Cartridge/Poly Bag**  
Place Impregum Pentax Super Quick Heavy Body or Medium Body poly bags in the Pentamix™ Pentax cartridge.  
Place the cartridge in the Pentamix mixing unit and attach a new Pentax Mixing Tip – Red.  
If a new mixing tip is installed when the cartridge is inserted, check whether the drive shaft is engaged before the start mixing.  
With newly filled cartridges, start mixing and extrude and discard the first unevenly mixed paste prior to the first use for impression taking. Do not use the paste to take an impression until the color of the paste is homogeneous.

**Prepare Garant Dispenser/Cartridge**  
At the start of a new cartridge and before every new mixing process, first extrude a quantity of paste without using a mixing tip until both components emerge evenly; remove any paste plugs that have formed in the cartridge openings.  
Then attach a Garant Mixing Tip Purple and a Garant Intra-oral Tip White. Alternatively, 3M Intra-oral Syringe – Purple can be used for direct application into the mouth.

**Prepare Intra-oral Syringe – Purple**  
Intra-oral Syringe – Purple can be used when applying the light body material. The Intra-oral Syringes are disposable syringes that are discarded after use on a patient.  
Please follow the Instructions for Use for the Intra-oral Syringe – Purple.

**Retraction**  
Retraction cords as well as sulfates and pastes based on aluminum chloride, aluminum sulfide, or iron sulfates are suitable retraction agents.  
Recommendation: 3M ESPE Astringent Retraction Paste. The paste containing aluminum chloride is directly applied from the retraction capsule into the sulcus. The paste can be used as an alternative to, or in combination with, retraction cords or other retraction methods, depending on the clinical situation and the patient's requirements.  
Areas from which impressions are to be taken should be kept dry.

**Working Time**  
Impregum Pentax Super Quick Heavy Body  
0:45 2:00  
Impregum Super Quick Medium Body  
0:45 2:00  
Impregum Super Quick Light Body  
0:45 2:00

**Do not use**  
Do not use Impregum Pentax Super Quick Heavy Body or Medium Body for impressions of the maxilla and mandible. Do not use Impregum Super Quick Light Body for impressions of the maxilla and mandible.

**Preparation**  
Remove the cap and mix the paste with the catalyst paste. The mixing ratio, based on volume, is 5 parts catalyst paste : 1 part catalyst paste.

**Application**  
Apply the paste to the prepared tray. The paste should be applied to the tray and allowed to dry completely according to the Polyether Adhesive Instructions for Use.

**Impression**  
Press the tray into the impression. The tray should be pressed into the impression until the color of the paste is homogeneous.

**Removal**  
Remove the tray from the impression. The tray should be removed from the impression until the color of the paste is homogeneous.

**Storage**  
Store the product at 18–25 °C/64–77 °F. Do not refrigerate! Do not use after the expiration date. Store impressions dry and below 30 °C/86 °F in the dark.

**Customer Information**  
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

**Warranty**  
3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for use. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

**Limitation of Liability**  
Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.  
Information valid as of August 2018

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► To avoid introducing bubbles into the model, briefly pre-rinse with water and dry with air before pouring. Do not use surfactants that impair the quality of polyether impressions and are not necessary.

**Cleaning and Disinfection of the Garant Cartridge**  
**Caution**  
Make sure that the cleaning and disinfectant agents you have chosen do not contain any of the following materials:  
• Oxidation agents (e.g., hydrogen peroxides)  
• Oils  
• Glutaraldehydes  
Use only validated methods for cleaning and disinfection. Use solely cleaning and disinfectant agents tested for efficacy and compatibility with the medical device (e.g., FDA certification). Always observe all applicable legal and hygiene regulations for dental offices and/or hospitals.  
**Limitations on Reprocessing**  
The Garant cartridge may be reused up to 18 times, provided that it is undamaged and clean.

**General**  
This cartridge is non-sterile when delivered and must be cleaned and disinfected before initial use and prior to every further use.

**Preparation for Decontamination**  
The cartridge should be reprocessed as soon as possible after every use.

**Reprocessing**  
The reprocessing procedure must be conducted before the initial use and after every further use. The cartridge must first be cleaned manually and, after cleaning, be manually disinfected in accordance with the following instructions.

**Manual Disinfection by Wiping (Intermediate Level Disinfection)**  
The cartridge is disinfected with ready-to-use disinfectant wipes (CaviWipes). The wipes contain the disinfectant CaviCide (alcohol-quaternary ammonium disinfectant). Follow the cleaning agent manufacturer's directions for use, observing in particular the contact times.

**Manual Disinfection by Wiping (Intermediate Level Disinfection)**  
The cartridge is disinfected with ready-to-use disinfectant wipes (CaviWipes). The wipes contain the disinfectant CaviCide (alcohol-quaternary ammonium disinfectant). Follow the cleaning agent manufacturer's directions for use, observing in particular the contact times.

**Preparation for Decontamination**  
The cartridge should be reprocessed as soon as possible after every use.

**Reprocessing**  
The reprocessing procedure must be conducted before the initial use and after every further use. The cartridge must first be cleaned manually and, after cleaning, be manually disinfected in accordance with the following instructions.

**Manual Disinfection by Wiping (Intermediate Level Disinfection)**  
The cartridge is disinfected with ready-to-use disinfectant wipes (CaviWipes). The wipes contain the disinfectant CaviCide (alcohol-quaternary ammonium disinfectant). Follow the cleaning agent manufacturer's directions for use, observing in particular the contact times.

**Manual Disinfection by Wiping (Intermediate Level Disinfection)**  
The cartridge is disinfected with ready-to-use disinfectant wipes (CaviWipes). The wipes contain the disinfectant CaviCide (alcohol-quaternary ammonium disinfectant). Follow the cleaning agent manufacturer's directions for use, observing in particular the contact times.

**Preparation for Decontamination**  
The cartridge should be reprocessed as soon as possible after every use.

**Reprocessing**  
The reprocessing procedure must be conducted before the initial use and after every further use. The cartridge must first be cleaned manually and, after cleaning, be manually disinfected in accordance with the following instructions.

**Manual Disinfection by Wiping (Intermediate Level Disinfection)**  
The cartridge is disinfected with ready-to-use disinfectant wipes (CaviWipes). The wipes contain the disinfectant CaviCide (alcohol-quaternary ammonium disinfectant). Follow the cleaning agent manufacturer's directions for use, observing in particular the contact times.

**Manual Disinfection by Wiping (Intermediate Level Disinfection)**  
The cartridge is disinfected with ready-to-use disinfectant wipes (CaviWipes). The wipes contain the disinfectant CaviCide (alcohol-quaternary ammonium disinfectant). Follow the cleaning agent manufacturer's directions for use, observing in particular the contact times.

**Preparation for Decontamination**  
The cartridge should be reprocessed as soon as possible after every use.

**Reprocessing**  
The reprocessing procedure must be conducted before the initial use and after every further use. The cartridge must first be cleaned manually and, after cleaning, be manually disinfected in accordance with the following instructions.

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