Attachment-Kleber Zirconia

Stand: 02/2026
Format: 500136_V260302

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GB Processing instructions

1 Intended purpose

Fixing composites/root pin cements are used for the adhesive, permanent bonding of restorative materials to restorative materials, to hard tooth structure or to a denture base, or to secure root pins in the root canal.

2 Description of product and user

2.1 Product specific description

Attachmentkleber Zirconia fluorescent white is a fluorescent white, purely chemical curing attachment composite for adhesive bonding of zirconia constructions to restorations made of zirconia, stainless steel, non-stainless steel and ceramic. The guidelines and specifications of DIN EN ISO 4049 apply.

2.2 Users

For use by laboratory technicians in a dental laboratory as well as for use in the dental practice by dental professionals.

3 Composition

Cristobalite sand, glass powder, di-urethane dimethacrylate, monofunctional aliphatic urethane acrylate, pyrogenic silica, dibenzoyl peroxide

Total filling materials: 62 wt.% inorganic filling materials: glass filler (average particle size 5 µm), Cristobalite sand (average particle size 3 µm), pyrogenic silica (average particle size 0,005µm-0,05µm)

4 Indications

- permanent attachment of indirect restorations made of metal (stainless steel, non-stainless steel, titanium), metal oxide ceramic (e.g. zirconia, aluminium oxide) or vitreous ceramic to zirconia implant abutments; zirconia constructions to metal (e.g. titanium) implant abutments
- Luting of titanium adhesive bases with individual zirconia abutments

5 Contraindications

- If a patient has allergies or hypersensitivities to a component of this product, it should not be used or used only under the strict supervision of the attending physician/dentist
- Attachment of restorations to natural hard tooth tissue
- If at least one of the components to be luted is not made of zirconia
- If safe cleaning and drying, or the specified application technique, is not possible
- Restorations that are to be attached to natural stumps and implant posts at the same time

6 Warnings

Contains: di-urethane dimethacrylate, monofunctional aliphatic urethane acrylate, dibenzoyl peroxide.

Attention. May cause an allergic skin reaction.

7 Safety instructions

Avoid breathing vapours/spray. Wear protective gloves.

8 Interactions with other agents

- Phenolic substances (such as Eugenol) inhibit polymerisation. Therefore, do not use any lining materials (such as zinc oxide-eugenol cements) which contain substances of this type.
- Oxidative disinfectants such as hydrogen peroxide may interact with the curing system and impair the curing process. The Automix syringe should therefore not be treated with oxidative disinfectants.
- The dentist should consider known interactions and cross-reactions of the medical device with other materials already in the patient's mouth before using the product.

9 Application

9.1 Luting of permanent attachments of indirect restorations made of metal (stainless steel, non-stainless steel, titanium), metal oxide ceramic (e.g. zirconia, aluminium oxide) or vitreous ceramic to zirconia implant abutments; zirconia constructions to metal (e.g. titanium) implant abutments

9.1.1 Removal of temporary

Remove the temporary or temporary abutment and clean objects to be luted thoroughly with alcohol (e.g. isopropanol) and then dry with oil-free air. Avoid repeated ingress of moisture or contamination.

9.1.2 Insertion and monitoring of final abutment

Only the original central screw should be used to permanent attach the suprastructure in the patient's mouth. The manufacturer's instructions are to be followed! A retraction cord can be applied to increase control over restoration fit and subsequently remove any surplus attachment composite.

9.1.3 Trial of restoration and allowing to dry

Check that the restoration is occluded and fits exactly. Exercise extreme caution when checking occlusion of brittle ceramic workpieces, as there is a risk of them cracking when in an unluted state.

Adjustments may be made using fine diamond-tip drills at moderate speed and light pressure if required. Milled surfaces should be re-polished.

Relative drying of the work area is essential when composites are used for adhesive fixing, preferably using wadding and dry angles. Absolute drying, e.g. by using a rubber dam, is recommended.

9.1.4 Surface preparation of abutment

- Seal the screw channel on the abutment using, e.g. a temporary filler, to prevent ingress of attachment composite.
- Clean and dry the abutment (see 9.1.1.)
- Apply a coating of a suitable luting agent in accordance with the manufacturer's instructions.

9.1.5 Surface preparation of restoration

The framework manufacturer's instructions will generally apply. In special cases, restorations should be prepared as follows:

Restorations made from metal, zirconia or aluminium oxide ceramic

- Clean and dry the restoration (see 9.1.1.)
- Sand blast the inside surface of the restoration with 50 µm aluminium oxide at max. 2 bar
- Dry restoration with oil-free air

Note: In order to achieve optimum bonding, do not clean with phosphoric acid.

Restorations made from lithium disilicate vitreous ceramic (e.g. IPS e.max Press, IPS e.max CAD)

- Etch with a commercially available hydrofluoric acid gel in accordance with the manufacturer's instructions.
- Rinse off the restoration thoroughly with water and dry with oil-free air.
- Apply a coating of a suitable luting agent in accordance with the manufacturer's instructions.

9.1.6 Application of Attachmentkleber Zirconia fluorescent white

Attach a disposable Automix tip to the dual chamber syringe. The initial quantities of mix from the syringe should not be used for luting. Squeeze out the Attachmentkleber Zirconia fluorescent white from the Automix syringe and apply the desired quantity directly into the restoration and spread evenly. Since the attachment material hardens in the mixing tip, the tip can be used to seal the contents of the syringe until next used.

9.1.7 Insertion of restoration and removal of surplus

Apply light pressure to bring restoration in situ and fix/hold. Remove surplus attachment composite immediately using a disposable brush, foam pellet, dental floss or an implant scaler. Pay particular attention to ensuring that the surplus is properly removed from difficult to reach areas (approximal, gingival margins). Formation of an oxygen inhibition layer can be prevented by coating restoration margins immediately after removal of surplus with glycerine gel / Airblock and rinsing off with water after curing is complete.

9.1.8 Finishing restoration

Remove the retraction cord if applied. Check occlusion and functional movement and adjust as required. Polish restoration margins with polishers or disks.

9.2 Luting of titanium adhesive bases with individual zirconia abutments

9.2.1 Surface preparation of titanium adhesive base

- Flaring of the titanium base can be reduced to the implant diameter in the luting area (zirconia/titanium join) from a subgingival, anatomic viewpoint. The minimum wall thickness of the zirconia portion should be taken into account here.
- The hex connection on the titanium base should be protected against damage and should therefore be screwed to a model analogue for luting.

Important note: screw with the laboratory screw only! The scan screw is not suitable for luting purposes as it cannot be removed if adhesive enters the microgap.

- Block out the screw access on the titanium base with wax.
- Clean the titanium bases thoroughly with alcohol (e.g. isopropanol) and then dry with oil-free air. Avoid repeated ingress of moisture or contamination.

- Carefully blast the titanium adhesive base at max. 2.5 bar and 50 µm aluminium oxide.
- Apply a coating of a suitable luting agent in accordance with the manufacturer's instructions.

9.2.2 Surface preparation of zirconia abutment

- Clean the abutments thoroughly with alcohol (e.g. isopropanol) and then dry with oil-free air. Avoid repeated ingress of moisture or contamination.
- Carefully blast the inside of the abutment at max. 1.0 bar and 50 µm aluminium oxide.
- Apply a coating of a suitable luting agent in accordance with the manufacturer's instructions.

9.2.3 Luting of abutment and removal of surplus

Attach a disposable Automix tip to the dual chamber syringe. The initial quantities of mix from the syringe should not be used for luting. Express the Attachmentkleber Zirconia fluorescent white from the Automix syringe and apply the desired quantity directly onto the components to be luted and spread evenly. Since the attachment material hardens in the mixing tip, the tip can be used to seal the contents of the syringe until next used. Any surplus adhesive in the screw channel should be removed carefully with the aid of a microscope. If the screw seating inside the titanium inserts subsequently needs to be smoothed, a suitable reamer should be used.

9.2.4 Finishing of the luted abutment

The conical zirconia abutments are reworked in the milling machine with water cooling after luting until a matte silk shine is achieved. This is followed by manual high-gloss polishing.

10 Notes on Processing

10.1 Processing time

Processing and setting times depend on the ambient temperature. The Attachmentkleber Zirconia fluorescent white should be processed quickly after removal from the Automix syringe and the restoration inserted rapidly. The following times apply once Attachmentkleber Zirconia fluorescent white is squeezed out of the Automix syringe:

* 11 Troubleshooting / FAQ Liste

Fault	Cause	Corrective action
Attachment does not grip after cement is applied	Bonding surface contaminated	Ensure work areas are free of contamination after sand blasting and application of primer
	Components of the cement are not evenly mixed (air bubbles)	Use a mixing tip and discard the first 0.25 g or so of the material squeezed out of the tip
	Application takes longer than processing time	Keep within processing time and only mix small quantities
	Objects move during curing time for cement	Keep within curing time and do not touch objects, fix in place if necessary
Cement does not cure	Use of self-etch bonding systems	Do not use non-system surface primers. Use original primers only.
	Contamination of bonding surfaces with products containing Eugenol	Use products that are free of Eugenol
	Cement stored incorrectly (too warm)	Purchase new cement and store in refrigerator
	Unequal quantities of the two components mixed in	Use original mixing tips for the syringe and discard the initial cement squeezed out

	Room temperature ca. 21 °C (ca. 69 °F)	Intraoral temperature 37 °C (98,6 °F)
Processing time	2-3 min	approx. 1 min
Curing time, incl. processing time	approx. 5 min	approx. 2-3 min

10.2 Mixing ratio

Use of the Automix syringe including mixing tip ensures that Attachmentkleber Zirconia fluorescent white is always mixed at the optimum ratio of 1:1.

11 Troubleshooting / FAQ List*

12 Storage and handling:

The Attachmentkleber Zirconia fluorescent white must be stored at 3°C-9°C. Leave the used mixing tip in place to seal the Automix syringe.

The paste should be at room temperature during processing as it may be difficult to squeeze out and mix at refrigerator temperatures.

13 Shelf life

The maximum shelf life is printed on the label of each syringe. Do not use after the expiry date.

14 Warnings on side effects

With proper preparation and use of this medical device, adverse effects are extremely rare. However, immune reactions (such as allergies) or local discomfort cannot in principle be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case.

15 Instructions for disposal

Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

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