DENTA-DON

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DIRECTIONS FOR USE

Heat-curing acrylic for the production of dentures by injection or packing procedure.

If there is anything in this patient information leaflet that you do not understand, please contact our customer service department before using the product.

Properties:

- color-stable
- excellent physical properties
- easy handling

Mixing ratio:

2.5 parts powder : 1 part liquid (by weight)

Processing:

- Insulate all plaster surfaces with the alginate separating agent.
- Add M+W Denta-Don powder to the liquid in a mixing bowl and stir well with a spatula for about 30 sec. to obtain a homogenous mass.
- The material can be processed after a dwelling time of 10 min.
- If a longer processing time is desired, we recommend to change the mixing ratio to 2,0 parts powder and 1,0 part liquid (by weight).
- If necessary, mixed acrylic may be stored in a refrigerator for several hours. The processing time at room temperature is:
mixing ratio 2,5 : 1 = ca. 1 hour
mixing ratio 2,0 : 1 = ca. 2 hours
- Subsequent to the dwelling time, the material must be kneaded for ca. 30 sec.
- The ductile M+W Denta-Don material is filled into the prepared flask and pressed.
- The material must be kept under pressure for 15 min.
- Before polymerization, a clamping bracket is put on the flask.
- In order to keep the acrylic teeth from breaking off, roughen their surfaces and prepare them with undercuts. Wet the teeth with a monomer or adhesive agent.

Polymerization:

Polymerize M+W Denta-Don in a closed flask held tight by a clamping bracket for 20 min. in boiling water. Let the denture cool down to room temperature inside the flask with the clamping bracket. Next, remove the denture from the flask using the traditional method. M+W Denta-Don obtains excellent properties also with a normal polymerization method or with a long-term polymerization method.

Processing in the UNIPRESS-System:

Process the flasks and plaster surfaces as described in the instructions for this system. M+W Denta-Don is mixed as described for the injection molding technique, placed in the cylinder and injected. The polymerization cycle is the same as that used for injection molding.

Repairs:

Repairs can be carried out with M+W Denta-Press, M+W Denta-Fast NF. Before applying the mixed acrylic, roughen and clean the fracture surfaces and coat them with monomer. M+W Denta-Don is polymerized for 15 min. in a pressure polymerizing vessel at a pressure of 2-5 bar with a water temperature of 45 °C.

Storage:

At 10-25 °C. Please observe the shelf-life.

Important:

If the temperature deviates from the room temperature: Higher temperatures decrease the processing times, lower temperatures increase the processing times.

Shades:

pink transparent

Presentation:

Powder in packages of 500 g or 1000 g.
Liquid in bottles of 500 ml or 1000 ml.

Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions:

If a patient has known allergies against or hypersensitivities towards a component of this product we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

Note:

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Warning:

Avoid skin contact with unpolymerized material and do not inhale monomer fumes.

Disposal:

Powder can be disposed of with household waste. Liquid must be disposed of as special waste.

Troubleshooting

Problem	Cause	Remedy
Acrylic does not fully cure	Incorrect mixing ratio	– Adhere to the instructions for use – Adhere to the mixing ratio
Acrylic is whitish on the fitting surface	– incorrect or poor application of separating agent – model not soaked	– Apply ample separating agent – Check the separating agent if necessary – Soak the model for approx. 10 min. in lukewarm water
Acrylic forms bubbles on the fitting surface	Model not soaked	– Soak the model for approx. 10 min. in lukewarm water
Acrylic surface white/porous	Recommended time for injecting acrylic exceeded or delay in inserting into the pressure vessel	Adhere to the processing times
Color differences with repairs	Incorrect mixing ratio	Adhere to the mixing ratio
Acrylic breaks when lifted	Poorly separated	Check the separating agent
Acrylic / denture has whitish streaks	Incorrectly mixed / too much powder used	Thoroughly spatulate the acrylic dough / adhere to the mixing ratio

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