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GB Instructions for use

Intended purpose

Light-curing material for coating/isolating the gingiva from the teeth to be treated during in-office bleaching

Description of product and user

M+W Gingiva Protector is a light-curing gingiva guard which is used similar to a rubber dam for isolating the gingiva from the teeth requiring treatment.

Patient target group

All patients who are undergoing tooth whitening with highly concentrated whitening gel in a dental practice. The M+W Gingiva Protector is used as a protective barrier and is applied prior to tooth whitening.

Users

For use in the dental practice by dental professionals

Composition

Aliphatic urethane dimethacrylate, tetramethylene dimethacrylate, glass powder sil., light initiators, stabilisers, pigment

Indications

For coating gingiva during in-office bleaching procedures

Contraindications

If a patient has allergies or hypersensitivities to a component of this product, it should not be used or used only under the strict supervision of the attending physician/dentist.

Warning

Contains aliphatic urethane dimethacrylate, tetramethylene dimethacrylate. May cause an allergic skin reaction.

Safety instructions

Avoid breathing vapours/spray. Wear protective gloves. If skin irritation or rash occurs: get medical advice/ medical attention.

Interactions with other medicinal products

The dentist should consider known interactions and cross-reactions of the medical product with other materials already in the patient's mouth before using the product.

Application / Preparation

Preparatory work

Use suitable measures (such as cheek/lip retractor) to ensure that tongue and lip are fixed during treatment (and during the subsequent bleaching procedure) and do not come into contact with the material. Remove the cap on the syringe and replace it with the top of the syringe. Check that the material flows freely.

Application

Dry the gingiva. Spread a layer of M+W Gingiva Protector approximately 4 - 6 mm wide and 2 mm deep along the area of application. (Photo 1)
Press material into the interdental spaces. Cure with a standard handheld polymerisation lamp. Polymerise for 20 to 40 sec, ensuring that the lamp does not stay in one place for too long. Carefully remove any surplus material left on the teeth after curing with a scalpel. The bleaching gel can now be applied. (Photo 2)

Completion

Remove M+W Gingiva Protector when the bleaching procedure is

completed. Use an instrument (e.g. a probe) to lift the material up at one end and draw it off. Remove the top of the syringe and replace it with the cap. (Photo 3)

Storage and handling

Close the device tightly after use. Store at 3°C - 25°C.

Shelf life

The maximum shelf life is printed on the label of each syringe. Do not use after the expiry date.

Side effects

With proper preparation and use of this medical device, adverse effects are extremely rare. However, immune reactions (such as allergies) or local discomfort cannot in principle be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case.

Instructions for disposal

Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

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