implantlink® semi Classic

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implantlink® semi Classic

GB
Instructions for use

Resin-based temporary luting cement for implant-retained restorations

Dual-curing semi-permanent implant cement with very low displacement resistance and thus extremely low film thickness. Due to highly cross-linked, ductile polymer structure, high compression strength and hermetic sealing are achieved. Excess material can easily be removed in large pieces (no crumbling) during elastic phase. Luted restorations are effortlessly removed without damage with standard devices (Crown Butler, Corona Flex). Compatible with all conventional material combinations, i.e. zirconium and metal frameworks, ceramics, titanium. Antibacterial, eugenol-free, odourless and tasteless; natural opaque colour for screening metallic abutments.


Fig. 1

1. Preparation and cleaning of implant abutments and inner surfaces of restorations

Carefully check fit, occlusion and proximal contact points. Complete all finishing procedures. Degrease, clean (e.g. alcoholic solutions) and dry prepared abutments. Make sure the interior surfaces of the superstructure are clean and dry. (Fig. 1).


Fig. 2

2. Mixing and applying

Bleed the syringe before installation of the automix cannula until material emerges uniformly from both openings (Fig. 2). Insert the plunger in the body of the cartridge. Remove the cap from the cartridge by turning anti-clockwise and attach mixing cannula. Pay attention to the guides on cannula and cartridge. Lock in place by turning clockwise (Fig. 3). Again, bleed small amount to ensure uniform and even material flow. Due to finely tuned dispensing pressure, dosing can be performed individually and precisely. Dispense the material with uniform pressure. After use, leave mixing cannula on cartridge until next application, thus sealing the cartridge.


Fig. 3

3. Inserting the restoration

Apply a thin layer of **implantlink® semi Classic** directly into the superstructure (Fig. 4). Carefully position the restoration on the implant abutment with slight pressure (Fig. 5); then press down firmly. Remove excess material during the gel phase, after around 2 – 3 minutes, with a suitable instrument (Fig. 6).

Fig. 4

4. Working time / light curing

The working time at room temperature is around 60 seconds. The cement is cured at body temperature in about 5 – 6 minutes. Direct exposure to light accelerates setting time. Light curing (marginal region, excess material) is achieved with approx. 20 seconds per surface (Fig. 7). When **implantlink® semi Classic** is used with translucent materials (e.g. thin layers of ceramic provisional crowns and bridges), increase exposure to light through the material to 30 - 40 seconds per surface.


Fig. 5

Precautions

- Do not rough or sand-blast abutment surface.
- Do not damage the underlying titanium surface when removing excess material. If necessary, use plastic instruments (implant scalers).
- Leave no cement residues in the mouth or gingival sulcus.
- Use no adhesion-reducing additives. They could uncontrollably render the retentive qualities of **implantlink® semi Classic**.
- Do not manually remove or mix cement components.
- Avoid contact of unmixed components with skin or mucosa.
- Observe the safety data sheet.
- Do not use **implantlink® semi Classic** for the luting of hybrid restorations with natural teeth and implant abutments.
- Only for intended use by qualified personnel.
- Please follow the instructions on the safety data sheet!



Fig. 6

Autopolymerisation:

Light-curing:


Fig. 7

Indication for Use:

- Temporary luting of implant-retained restorations

Contra-indication:

- In very rare cases **implantlink® semi Classic** may cause sensitization. If such reactions should occur, discontinue using **implantlink® semi Classic**.

Technical data:

- **Mixed volume:** 5 ml (mini-mix)
- **Mixing ratio:** 4:1
- **Colour code:** Base: white-opaque
Catalyst: semi-transparent
- **Mixing time:** None: automix (mini-mix)
- **Working time:** approx. 60 sec.
- **Gel phase:** 2 - 3 min.* (for removal of excess material during placement)
- **Setting time in the mouth:** 5 - 6 min.*
- **Setting time with light-curing:** approx. 20 sec. per surface; 30-40 sec. through translucent materials
- **Film thickness:** < 10 µm
- **Compression strength:** approx. 85 MPa
- **Application temperature:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity
- **Storage:**

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity. Higher temperatures accelerate, lower temperatures retard a. m. times.

Ordering information:

implantlink® semi Classic
Standard packing **03092**
cartridge of 5 ml mini-mix 4:1
10 mixing cannulas, brown

implantlink® semi Forte
Standard packing **03371**
cartridge of 5 ml mini-mix 4:1
10 mixing cannulas, brown

implantlink® semi Xray
Standard packing **02195**
cartridge of 5 ml mini-mix 4:1
10 mixing cannulas, brown

Mixing cannulas 4:1 **02591**
brown, 25 pcs.

monoprint supra
Transfer silicone
Standard packing **02525**
2 x 50 ml cartridges
6 mixing cannulas, green

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