3. Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e., by preserving as much of the tooth structure as possible. Do not prepare sharp internal edges and angles or additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. Bevel enamel margins of anterior teeth. In the posterior region, only the sharp enamel edges should be rounded (finishing diamonds, 25 – 40 µm). Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes using rubber cups or rotary brushes. Subsequently, remove all residue in the cavity with water spray and dry with water- and oil-free air.

4. Pulp protection / Base

Do not apply a base material when using an enamel/dentin bonding agent. Selectively coat only very deep areas close to the pulp with a calcium hydroxide material (e.g. ApexCal®) and subsequently cover with a pressure-resistant cement (e.g. a glass-ionomer cement, such as Vivaglass® Liner). Do not cover the other cavity walls so that they can be used to support the bond with an enamel/dentin adhesive.

5. Apply matrix / interdental wedge

Use a wrap around matrix for cavities affecting the proximal area or a sectional matrix and wedge it.

6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the Instructions for Use of the product in use. Ivoclar Vivadent recommends using the universal bonding agent Adhese® Universal, a single-component, light-cured adhesive for direct and indirect bonding procedures that features compatibility with all etching techniques.

7. Application of Tetric Prime

- For an optimum result, apply Tetric Prime layers of max. 2 mm (or 1.5 mm Dentin shades) and adapt with a suitable instrument (e.g. OptraSculpt).
 - Sufficient exposure to the curing light prevents incomplete polymerization.
 - For the recommendations regarding exposure time per increment and light intensity see table 1.
- The instructions for use of the curing light must be observed.**
- When using a metal matrix, additionally polymerize the composite material from the buccal or the lingual/palatal aspect after removing the matrix.
 - A flowable composite (e.g. Tetric EvoFlow® or Tetric PowerFlow) can optionally be used. This optional layer must be cured separately (please refer to the respective Instructions for Use).

8. Finishing / Checking the occlusion / Polishing

Remove excess material with diamond or tungsten carbide finishers after polymerization. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restorations. Use silicone polishers (e.g. OptraPol®) as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Additional notes

1. Tetric Prime can be used in combination with Tetric EvoFlow, Tetric PowerFill and Tetric PowerFlow.
2. In the case of repairs, additional Tetric Prime can be directly applied to the polymerized material. If the Tetric Prime restoration has already been polished, it must first be roughened and wetted with Adhese Universal before fresh Tetric Prime is applied.
3. Tetric PowerFill should have ambient temperature when applied. Cool temperatures render the material difficult to extrude.
4.  For single use only. If Tetric Prime is applied from the Cavifil directly into the oral cavity of the patient, the Cavifil or the application cannula of the syringe should only be used for one patient due to hygienic reasons (prevention of cross-contamination between patients).
5. Syringes or Cavifils should not be disinfected with oxidizing disinfection agents.

Warning

Avoid contact of unpolymerized Tetric Prime with skin, mucous membrane and eyes. Unpolymerized Tetric Prime may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage

- Storage temperature 2 – 28 °C (36–82 °F)
- Close syringes/Cavifils immediately after usage. Exposure to light leads to premature polymerization.
- Do not use Tetric Prime after the indicated date of expiration.
- Date of expiration: see information on Cavifils, syringes and packages.

Keep out of the reach of children!

For use in dentistry only!

The material has been developed solely for use in dentistry and should be used strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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