# VisCalor

## Dispenser

**USER MANUAL  
BENUTZERHANDBUCH  
MANUEL D'UTILISATION  
MANUAL DEL USUARIO  
MANUALE D'USO**

Deutsch

Français

Español

Italiano

# VOCO

# TABLE OF CONTENTS

Introduction.....	3
Benutzerhandbuch.....	36
Manuel D'Utilisation.....	69
Manual Del Usuario .....	101
Manuale D'Uso .....	135

## **Important!**

***PLEASE NOTE!*** Prior to installation and start-up of the device, please read these instructions carefully. As with all technical devices, the proper function and safe operation of this device depend on the user's compliance with the standard safety procedures as well as the specific safety recommendations presented in these Operating Instructions.

For use by qualified professionals only!

## **Introduction**

Congratulations on your decision to incorporate the VisCalor® Dispenser heated composite delivery system into your practice. This patent-pending device is the first of its kind, using NIR (near infrared) technology to rapidly warm highly filled composite compules. With the touch of a button, VisCalor® Dispenser is able to heat composite material to 65°C (150°F) in seconds, and maintain that temperature throughout the procedure.

Warming composite significantly lowers the viscosity of the material resulting in better adaptation, reduced voids and microleakage. Materials remain highly sculptable, non-sticky, and easily shaped during manipulation.

The VisCalor® Dispenser is designed and manufactured by InterMed, Inc. and distributed exclusively by VOCO.

**NOTE: ALL COMPONENTS ARE NON-STERILE**

### **Intended Use / Indications for Use**

Heating and dispensing of dental composite materials.

### **Contents of the VisCalor® Dispenser Kit**

- (1) Handpiece
- (3) Changeable Tops
- (1) Stand
- (1) Power Supply
- (1) Device Maintenance Kit
- (100) Protective Sleeves

**NOTE: ALL COMPONENTS ARE NON-STERILE**

### **Intended Use / Indications for Use**

Heating and dispensing of dental composite materials.

### ***Unpacking the Container***

No special assistance is required to unpack and assemble the VisCalor® Dispenser. Packaging should be inspected upon arrival for evidence of shipping damage. Damaged packaging may indicate the presence of an unsafe product and the product should not be used until carefully inspected.

### Introduction

Please reference the image below to familiarize yourself with the VisCalor® Dispenser.

### *Start Heating Cycle*

#### *LED Indicators*

The VisCalor® Dispenser is designed to display the battery state via the battery level indicator:

The LED light surrounding the ON/OFF button indicates current charge level

(< 20% charged - **Red**, > 20% charged - **Green**).

It is advised to charge the battery when the indication lights are **Red**.

### *Charging the Battery*

The dispenser battery is partially charged when shipped. To charge the battery, simply connect the micro-USB charger to the USB charging port on the device handle and plug the power adapter into a standard outlet.

During battery charging, the device activity indicator light will pulse.

Once charging is complete, the device activity indicator will remain blue.

When unplugged from the charger, the device activity indicator turns off and the device will enter a sleep mode.

### *Inserting Composite*

Snap a composite compule into the end of the device. Position the compule so that the tip is directly under the LED lens.

**The capsule must always be placed at the back position**

Note: The device will not operate in the absence of a composite compule.

### *Settings*

Settings 1	VisCalor® bulk
Settings 2	All other VOCO composite caps

**Please use only composite caps from VOCO with the VisCalor® Dispenser!**

**CAUTION:** Parts of the composite compule may become hot. Use caution and avoid touching the top of a hot composite compule.

### Selecting On Setting

Press & hold to enter programming mode. This is indicated by flashing yellow lights.

Short press to toggle between settings.

Press & hold to exit. Device will beep 3 times and shut off.

Your setting is now memorized as default.

Short press to start/stop heating cycle. The blue lights will flash while heating.

**Setting 1:** 30 seconds

**Setting 2:** 70 seconds

The blue lights will turn solid when heating cycle is complete.

The device will automatically shut off after 3 minutes of run time.

## Cleaning & Disinfection

The device may be wiped down with a standard disinfectant wipe. Never spray disinfectant directly onto the device.

**CAUTION:** Always use a protective sleeve.

Black changeable top is removable and autoclavable.

Between uses, check the lens to ensure it is clean. A dirty lens will decrease performance.

The VisCalor® Dispenser is provided non-sterile. There are no special accessories needed to disinfect the VisCalor® Dispenser heated composite delivery system.

Quaternary ammonium compound products are recommended (containing 20% alcohol or less). Wipe, do not spray, solution onto the unit. Prevent liquids from entering openings on the VisCalor® Dispenser unit.

***DO NOT AUTOCLAVE HANDPIECE.***

***CHANGEABLE TOPS CAN WITHSTAND AUTOCLAVE - (134°C / 0.22mPA for 5 minutes).***

### **CAUTION:**

- ***DO NOT*** immerse the unit or unit parts in solutions. Use of solutions other than those recommended may damage plastic parts and will void product warranty.
- ***DO NOT*** use abrasive material such as scouring powder, organic solvents, or solvent-based cleaning fluids. In case of severe contamination, gently clean the device by using diluted alcohol.
- Store the device in the box if it is not to be used for an extended period of time.

### *Infection Control Measures*

**CAUTION:** To prevent cross-contamination, a disposable plastic sleeve must be used over the VisCalor® Dispenser with each use. A low-density polyethylene plastic disposable sleeve covers the changeable top of unit and provides a hermetically sealed barrier between the handpiece and patients. The disposable sleeve limits patient-to-patient contamination. Discard used protective sleeves after each patient.

### *Routine Maintenance*

Between uses, check the lens to ensure it is clean.  
A dirty lens will decrease performance.

1. Insert the hex key provided into the trigger bolt.

2. Rotate counterclockwise and remove the trigger bolt once loose.

3. Slide the trigger and spring assembly apart from the device.

4. Tilt the device forward to allow the plunger to slide forward and be removed from the device.

5. With the changeable top removed, clean the device with the provided brush.


6. Thoroughly clean the plunger space and LED lens to remove any composite residue.

**NOTE:**

- Alcohol may be used as a cleaning solution.
- **DO NOT** insert the brush past the end of the bristles.

### Trigger Assembly

After the item is cleaned, replace the trigger.


7. Move the **plunger** in so that the slotted end is visible.


8. Place spring in trigger with short bent leg towards top.

**9.** When inserting the trigger, the spring leg must go into this **hole**.

**10.** Place the trigger all the way in, pull and release until properly seated on the slots of the plunger. When seated properly the plunger will move forward and backward.

**11.** Move the silver spring wire toward the nozzle with the hex tool and insert the screw. Tighten the bolt so it bottoms out and is snug. Do not overtighten.

**Hex Tool**

## Battery Removal & Replacement

To remove battery compartment door, press the release button and separate from device by pulling.

Remove the old battery and replace with the new battery, ensuring that +/- indicators are properly aligned as indicated on the battery door.

## Troubleshooting Guide

If the suggested solutions do not rectify the problem, please contact VOCO.


Problem	Possible Solution
<b>VisCalor® Dispenser will not turn on</b>	<ol style="list-style-type: none"><li>1. Check the unit's battery indication light. If red or no light, charge the unit.</li><li>2. Remove the battery and check orientation. Reinsert into the device in proper orientation.</li></ol>
<b>VisCalor® Dispenser is not heating</b>	<ol style="list-style-type: none"><li>1. Verify the unit is charged.</li><li>2. Verify that an appropriate temperature is selected.</li><li>3. With the handpiece turned off, and battery removed, inspect the lens for residual composites.</li></ol>



Problem	Possible Solution
<p><b>Battery will not charge</b></p>	<ol style="list-style-type: none"> <li>1. Inspect the USB charging port of any foreign objects that could interfere with a proper connection.</li> <li>2. Remove the battery and check orientation. Reinsert into the device in proper orientation.</li> <li>3. Make sure the battery is properly inserted in the handpiece.</li> <li>4. Make sure the charger is plugged in and verify the outlet is receiving power.</li> </ol>
<p><b>Composite compules are melting</b></p>	<ol style="list-style-type: none"> <li>1. Ensure that you are using a composite cap from VOOCO.</li> <li>2. Ensure the device is in the correct setting.</li> </ol>

## Technical Data

Technical Information	VisCalor® Dispenser
<p><b>Charger</b></p>	<p><b>Input:</b> 100-240 VAC, 50-60 Hz  <b>Nominal Consumption:</b> 6W max  <b>Manufacturer:</b> GlobeTek INC.  <b>Model:</b> GTM46101-1005-USB  <b>Dimensions without blade or cable (LxWxH):</b> 41mm x 71mm x 31.5mm  <b>Mass:</b> 50g  <b>Classification:</b> Protection class II, </p>

<b>Technical Information</b>	<b>VisCalor® Dispenser Gun</b>
<b>Handpiece</b>	<p> <b>Battery:</b> 3.6V nominal, 3000mAh NiMH Li-ion, 10.8 Wh. IEC 62133 certified  <b>Battery Pack Manufacturer:</b> LG CHEM LTD  <b>Battery Pack Model:</b> LG-HG2-18650-1NR  <b>Dimensions (LxWxH):</b> 150mm x 120mm x 25mm  <b>Mass:</b> 153g  <b>Classification:</b> Type BF,   <b>Intermittent Operation:</b> The device has been designed solely for short-term operation.  <b>Operating Time:</b>            Approximately 2 hrs. fully charged.            Approximately 15 heat cycles on high setting.         </p>

<b>Technical Information</b>	<b>VisCalor® Dispenser Gun</b>
<b>Operating Conditions</b>	<p> <b>Time to Charge Empty Battery Pack:</b> Approximately 3 hrs  <b>Nominal Consumption:</b> 6W max  <b>Temperature:</b> 10°C - 40°C (59°F - 104°F)         </p> <div style="text-align: center;">  </div> <p> <b>Atmospheric Pressure:</b> 697hPa - 1013hPa         </p> <div style="text-align: center;">  </div>

Technical Information	VisCalor® Dispenser Gun
<p><b>Transport and Storage Conditions</b></p>	<p><b>Temperature:</b> -20°C - 40°C (-4°F - 104°F)</p>  <p><b>Relative Humidity:</b> 30% - 90% (non condensing)</p>  <p><b>Atmospheric Pressure:</b> 500hPa - 1400hPa</p> 

## Symbol Identification

Description for additional symbols.

	Serial Number		Consult instructions for use
	Manufacturer		Temperature Limitation
	Manufacturer Date		Humidity Limitation
	Class II Medical Electrical Equipment		Pressure Limitation
	Type BF Patient Applied Part		Batch Code / Lot Number
	Keep Dry		European Representative
	Part Number		Do not use if seal or packaging is compromised

## Symbol Identification

	<p>This symbol is a mandatory marking for devices entering the European market to indicate conformity with the essential health and safety requirements set out in European Directives</p>	<p><b>Rx ONLY</b></p>	<p><b>CAUTION:</b> U.S. federal law restricts this device to sale by or on the order of a dental professional.</p>
	<p>Do not reuse</p>		<p>This symbol refers to the special disposal of electrical and electronic devices in EU countries. Please do not discard this device in household garbage. Check the proper means of disposal in your country at your community recycling, waste center or at your dealer. Take care to dispose of properly.</p>
	<p>Warning/Caution</p>		<p>Autoclavable up to the temperature specified</p>

## Battery Disposal

Batteries contain toxic material and should not be disposed of in landfills or incinerators. Dispose of depleted batteries as directed by your local solid waste handling regulations.

### **Safety Notes, Warnings and Precautions**

Read all instructions before operating this unit. The VisCalor® Dispenser heated composite delivery system emits high intensity light waves and must only be used as indicated in this manual.

The VisCalor® Dispenser is a medical device which is subject to IEC 60601-1 (EN 60601-1) and EMC directives IEC 60601-1-2 (EN 60601-1-2) Edition 4.0, as well as the 93/42/EEC Medical Device Directive. The VisCalor® Dispenser complies with the relevant EU regulations.

The device has been shipped from the manufacturer in a safe and technically sound condition. In order to maintain this condition and to ensure risk-free operation, the notes and regulations in these Instructions for Use have to be observed. To prevent damage to equipment and risk for patients, users and third parties, the following safety instructions have to be observed.

#### **WARNING**

- The user should test the product before use to ensure proper functionality.
- As with any heavily used medical device, the user needs to ensure a functional backup is readily available.

- **DO NOT** keep or position near flammable materials, or materials that could combust.
- **DO NOT** look directly into the IR light output when the device is on.
- **DO NOT** insert fingers, instruments, or other objects into the handpiece when the battery pack is removed.
- **DO NOT** autoclave the handpiece, battery, charging cord, or stand.
- **ONLY** heat composite starting at room temperature.
- **DO NOT** heat composite if it is already at an elevated temperature from use in the VisCalor® Dispenser or other warming devices.
- **ONLY** recommended batteries can be used with the VisCalor® Dispenser. Reference table on page 20 for battery specifications.
- **ONLY DISPENSE** composite after warming. Not intended to dispense composite at room temperature or below.

### CAUTION

- U.S. Federal law restricts the sale of this device by or on the order of a healthcare professional. Use of the device is restricted to qualified and trained personnel only in accordance with the operation instructions. The manufacturer assumes no liability for any damage arising from any other or improper use of this device.
- Only use the charger which is provided with the device. The use of any other charger can result in damage to the device.
- Condensation resulting from the device being transferred from a cold to a warm environment may be a potential risk. Never begin operating the device until it has reached the ambient temperature.
- Use only components and accessories listed in the instructions associated with the device. Failure to do so will void the warranty, may decrease the performance, and may lead to unsafe operation.
- In order to avoid electric shock, do not introduce any objects into the device or remove the device enclosure.
- Should you have any reason to suspect the safety of the device to be compromised, the device must be taken out of operation and labeled accordingly to prevent third parties from inadvertently using a possibly defective device. Safety may be compromised,

e.g., if the device malfunctions or is noticeably damaged.

- Keep solvents, flammable liquids, and sources of intense heat away from the device as they may damage the plastic housing of the device, the seals, or the operating buttons.

## **Operating Conditions and Safety Consideration**

### ***Heat Generation***

The VisCalor® Dispenser heated composite delivery system has been designed not to overheat to the point of discomfort or injury during standard operating durations. However, care should be taken to allow the product to completely cool between uses (approximately two-three minutes), to ensure overheating does not occur.

### ***Cool-Down***

VisCalor® Dispenser features a built in fan which will turn on automatically during operation. The fan will remain on after operation until the device reaches a lower temperature.

### ***Adverse Reactions***

There are no known adverse reactions.

### ***Contraindications***

VisCalor® Dispenser should not be used for any application outside of composite restorations.

- **DO NOT** heat composite if it is already at an elevated temperature from in the VisCalor® Dispenser or other warming devices.
- **ONLY DISPENSE** composite after warming. Not intended to dispense composite at room temperature or below.

### *Warranty*

VOCO gives a 12-month warranty period if the device cannot be used as intended due to quality defects. This warranty begins on the day of purchase. VOCO reserves the right to either repair the device or replace individual parts or the entire device at its discretion. This warranty does not cover fuses, wearing parts or mains adapters.

### *General information*

- We cannot accept any liability for malfunctions or damage resulting from improper use, the removal of components, modifications or inadequate maintenance.
- We cannot accept any liability for malfunctions or damage resulting from operation in unsuitable environments.
- We cannot accept any liability for malfunctions or damage resulting from acts of God.
- The warranty shall also be void in the case of product defects which can be definitively attributed to natural disasters or

excessive temperatures that do not conform to the specifications in this user manual.

- We cannot accept any liability for damages caused by the use of unsuitable composite caps or the use of the false program.

## **Appendix – Electromagnetic Compatibility and Electrical Safety Information**

The VisCalor® Dispenser heated composite delivery system is tested according to IEC 60601-1-2, Edition 4.0. Medical electrical devices are subject to particular preventive action according to EMC rules and must be installed and operated according to the EMC guidelines in the accompanying documents.

### **Guidance and Manufacturer's Declaration – Electromagnetic Emission**

The following tables are guidelines according to the 4th edition of the medical standard IEC 60601-1-2.

VisCalor® Dispenser heated composite delivery system is intended for use in the electromagnetic environment specified below. The customer or the user of the VisCalor® Dispenser heated composite delivery system should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	VisCalor® Dispenser uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	VisCalor® Dispenser is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Complies	N/A
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	N/A

Table: According to IEC 60601-1-2, Edition 4.0

## Guidance and Manufacturer's Declaration – Electromagnetic Emission

The VisCalor® Dispenser heated composite delivery system is intended for use in the electromagnetic environment specified below. The customer or the user of the VisCalor® Dispenser heated composite delivery system should

assure that it is used in such an environment.

Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 15 kV air	+/- 6 kV contact +/- 15 kV air	Floors should be concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electric fast transient / burst IEC 61000-4-4	+/- 2 kV for power supply lines	+/- 2 kV for power supply lines	Mains power quality should be that of typical commercial or dental environment.
Surge IEC 61000-4-4	+/- 1 kV - line +/- 2 kV - line - earth	+/- 1 kV - line No prot. earth	Mains power quality should be that of typical commercial or dental environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles	Mains power quality should be that of typical commercial or dental environment. If the user of the VisCalor® Dispenser requires continued operation during power mains interruptions, it is recommended that the VisCalor® Dispenser be powered from an uninterruptible power supply or battery

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
	70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or dental environment.

Portable and mobile RF communications equipment should not be used closer to any part of VisCalor® Dispenser, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

### Recommended separation distance

Conducted RF IEC 61000-4-6	3 Vrms 150kHz to 80 MHz	3V	VisCalor® Dispenser $d = 0.35 \sqrt{P}$
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5 GHz	10V/m	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.70 \sqrt{P}$ 80 MHz to 2.5 MHz

<p>Radiated RF IEC 61000-4-3</p>	<p>3 V/m 80MHz to 2.5 GHz</p>	<p>10V/m</p>	<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b). Interference may occur in the vicinity of equipment marked with the following symbol.</p> 
--------------------------------------	-----------------------------------	--------------	--

Table: According to IEC 60601-1-2, Edition 4.0

**NOTE:** UT is the a.c. mains voltage prior to application of the test level.

**NOTE:** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**a** - Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered.

If the measured field strength in the location in which the VisCalor®

Dispenser heated composite delivery system is used exceeds the applicable RF compliance level above, the VisCalor® Dispenser should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the VisCalor® Dispenser heated composite delivery system.

**b** - Over the frequency range 150 kHz to 80 MHz, field strength should be less than 10 V/m.

### **Recommended Separation Distances Between Portable and Mobile RF communications equipment and the “VisCalor® Dispenser Heated Composite Delivery System”**

The VisCalor® Dispenser heated composite delivery system is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the VisCalor® Dispenser heated composite delivery system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the VisCalor® Dispenser heated composite delivery system as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150kHz to 80MHz $d = 0.35 \sqrt{P}$	80 MHz to 800MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	0.035	0.12	0.23
0.1	0.11	0.38	0.73
1	0.35	1.2	2.3
10	1.1	3.8	7.3
100	3.5	12	23

Table: According to IEC 60601-1-2, Edition 4.0

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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**VOCO GmbH**
**Inter-Med, Inc.**
Racine, WI 53404-2500 - U.S.A.