BIG-ETCH

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Anwendung

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Instructions for use

Notes on application

Phosphoric acid etching gel for the conditioning of dental enamel and dentine in the dental field. M+W Big Etch is a 36% phosphoric acid etching gel with optimal thixotropic flow characteristics. It is used for etching tooth substance within the framework of the enamel etching technique and total-etch technique. According to current knowledge, the product has no side effects if the safety and application instructions are observed. Underfilling materials can be dissolved by the etching gel during longer exposure times.

- For dental use only!
- Do not use beyond the shelf life date
- Close refill syringe immediately after use!

Instructions of use

The sealing cap is twisted off counter clockwise and an applicator cannula is attached.

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Preparation

The areas of the tooth substance to be etched must be clearly delimited.

- 1.) The use of a rubber dam is recommended! Any contact with soft tissue (skin, oral mucosa, eyes etc.) is to be avoided. If etching gel ends up on soft tissue, immediately rinse the area with plenty of water and, if necessary, contact a doctor.
- 2.) Before preparation, clean the tooth using a rubber cup and fluoride- and fat-free polishing paste.
- 3.) Use water to properly remove cleaning residues, clean, and dry the surface of the tooth with oil-free compressed air.
- 4.) Prepare the cavity accordingly; bevel the melted edges to achieve permanent enamel junctions. Rinse the cavity thoroughly with water and dry with oil-free compressed air. Avoid over-drying!

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5.) Protect dental pulp in deep cavities (e.g. cover with a thin layer of suitable underfilling material). The cavity should be dry and free of contamination, whereby the over-drying of the dentine is to be prevented.

Etching process

6.) Apply etching gel to the enamel margins of the cavity and, if necessary, to the dentin.

Recommended application times:

- 15 - 30 sec. on tooth enamel
- 10 - 15 sec. on dentin
- 50 - 80 sec. on milk teeth and teeth with a high content of fluoride

7.) Rinse for at least 20 sec. with water.

8.) Dry with oil- and water-free compressed air. The etched area must not be touched or contaminated by saliva until the treatment is

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completed. However, should contamination occur after the etching and drying, the area must be etched again briefly (approx. 10 sec.). The etched enamel should look chalky white! After etching, you can continue with the application of adhesives etc. in accordance with the manufacturer's instructions for use. Do not reuse application cannulas

Safety instructions

H 314 Causes severe skin burns and eye damage. Wear protective gloves / eye protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. IF SWALLOWED: Rinse mouth. DO NOT induce vomiting. IN ON SKIN: Wash with plenty of water/soap.

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PRISMAN GmbH · Otto-Hahn-Ring 6-18 · D-64653 Lorsch · Tel.: +49 (0) 6251-866980-0
Müller & Weygandt GmbH · Reichardsweide 40 · D-63654 Büdingen · Tel.: +49 (0) 60 42 - 88 00 88