Instructions

SUGGESTED PIN PLACEMENT PROCEDURE

- Determine pin position at least 1mm from the dentino-enamel junction and indent channel sites lightly with round bur (**Fig. 1**).
- Drill pin channel in dentin with KODEX Depth Limiting Drill. Drill in two passes; first to half depth, remove drill and clean debris from drill and pin channel, then drill to full depth. During drilling and withdrawal, the drill must be kept rotating continuously (**Fig. 2**). Dry pin channel and apply cavity varnish with a Hygenic Paper Point. Use a latch-type slow speed handpiece for proper channel sizing and pin placement.
- Insert LINK PLUS pin into handpiece or Universal Hand Driver. Position pin over channel before activating for insertion.
- Using very low speed, apply light downward pressure. The LINK PLUS pin will automatically align and thread into the drilled channel and shear when seated at the shoulder stop (**Fig. 3**).
- TMS Bending Tool B-61-R or B-61-M may be used to bend pins to match the contours of the final restoration. Build up core or restoration with material of choice (**Fig. 4**).

NOTES AND PRECAUTIONS

Rotate drill in a clockwise direction. The drill must be kept in a continuous clockwise rotation during the entire drilling procedure until the drill has been fully withdrawn from the channel.

Dull drills may cause overizing of channels. Always use sharp drills. (Recommend maximum of 20 drilled channels per drill.)

Do not bend pins more than 45-degrees.

CLEANING, DISINFECTION AND STERILIZATION

Drills, TMS Bending Tool and Hand Wrench: Prior to each use, remove debris ultrasonically and sterilize by autoclave ONLY between patient use. (To minimize corrosion, place drills in separate bags from other instruments). Immediately air blow dry. A disinfectant may be used prior to use.

Pins: Prior to each use, disinfect or sterilize pins (autoclave ONLY).

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