FotoDent® denture

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FotoDent® denture

Important note regarding use

Please read these instructions carefully and in their entirety before using the product for the first time. The product must only be used as described in these instructions. Do not use the product for any purpose other than that specified. Please follow all safety precautions and instructions contained herein. The product must only be used by persons who have read and understood these instructions.

We assume no liability for damage resulting from improper use, mishandling, or failure to follow these instructions. This does not affect mandatory liability under the Product Liability Act, liability for injury to life, limb, or health, or liability arising from willful misconduct or gross negligence.

Intended use

Light curing material for production of dental denture base.

Indication for use: Light curing material for production of full and partial re-movable dentures using 3D printing process.

Indication: Removable denture basis

Contraindication: An allergic reaction might occur in case of any allergy to (meth)acrylates.

Side effects: There are no known side effects when used and processed properly. However, immune reactions (e.g. allergies) or local disorders of sensation (e.g. irritation of the oral mucosa) cannot be ruled out in principle. If in doubt, an allergy test should be performed before use. Information about potential allergic ingredients can be obtained from Dreve Dentamid GmbH.

Patient target group: Patients with a partially toothed or completely edentulous upper and/or lower jaw.

Intended users: The dental technicians or specialist staff in the dental laboratory is responsible for the production of denture bases in 3D printing. The dentist uses the finished denture base on the patient.

1. Requirements

Printer: Rapid Shape: D20iI, D30iI, D90iI, D20+, D30+ (385nm), Asiga: MAX UV, PRO 4K (385nm), Dreve: Phrozen Sonic XL 4K (qualified by Dreve) / Phrozen Sonic XL 4K 2022 (qualified by Dreve)

Software: Autodesk Netfabb®, Composer, Dreve ElementS **Post curing uni:** PCU LED N, FotoDent® flash, Ofotflash G171, RS cure Please contact us if you are using any other type of equipment.

Basic Material Properties

See table 1. Please refer to Technical Data Sheet (TDS). Available on request.

2. Specific Manufacturing Considerations

2.1 Digital denture base file

2.1.1 File format: STL file

2.1.2 Digital design: denture base or baseplate

2.2 Asiga printer and software

2.2.1 Hardware: Please refer to the printer manufacturer's manual for this information.

2.2.2 Asiga Printer Software (Composer): Please refer to the printer manufacturer's manual for this information

2.2.3 Printing parameters

- Slice thickness: 50µm
- Optimal Orientation: 0–15 degree tilted orientation
- Support point size: varies based on support type chosen

d) Support density: perimeter of denture base and labial region. The labial side faces the build platform.

2.2.4 Environmental Conditions

- a) MAX UV/PRO 4K print room temperature: 35 °C ± 3 °C
- b) Humidity: 30–70 %

2.3 Rapid Shape printer and software

2.3.1 Hardware: Please refer to the printer manufacturer's manual for this information.

2.3.2 Rapid Shape Printer Software (Autodesk netfabb®): Please refer to the printer manufacturer's manual for this information.

2.3.3 Printing Parameters

- a) Slice thickness: 50µm
- b) Optimal Orientation: 0–15 degree tilted orientation
- c) Support point size: varies based on support type chosen
- d) Support density: perimeter of denture base and labial region. The labial side faces the build platform.

2.3.4 Environmental conditions

- a) Print temperature: Please refer to the printer manufacturer's manual for this information.
- b) Humidity: 30–70 %

2.4 Dreve Phrozen printer and software

2.4.1 Hardware: Please refer to the printer manufacturer's manual for this information.

2.4.2 Dreve Phrozen Printer Software (Dreve ElementS): Please refer to the printer manufacturer's manual for this information.

2.4.3 Printing Parameters

- a) Slice thickness: 100µm
- b) Optimal Orientation: 10 degrees ± 2.5 degrees inclined orientation
- c) Support point size: varies based on support type chosen
- d) Support density: perimeter of denture base and labial region. The labial side faces the build platform.

2.4.4 Environmental conditions

- a) Print temperature: Please refer to the printer manufacturer's manual for this information.
- b) Humidity: 30–70 %

2.5 Recommended curing light equipment (Post curing units)

2.5.1 UV curing Flood Type System

See overview post-curing 2.5.1

2.6 Notification

The specifications and biocompatibility have been validated using printer software (Autodesk Netfabb®, Composer, Dreve ElementS) and aforementioned printer and process parameters indicated.

3. Process description

3.1 Preparation

- Optimal orientation of denture base for printing: 0–15 degree tilted orientation. The labial side faces the build platform, for for example as shown in the figures under 3.1

- During the application of FotoDent® denture we recommend wearing personal protective equipment i.e. suitable gloves, safety goggles etc.
- FotoDent® denture should be well shaken for a minimum of 30 seconds, 1 hour before use.

- Carefully pour FotoDent® denture into the pre-determined container of the production unit.
- Please see further steps within printer manual.
- Remove all bubbles with a cleaned object.

Processing temperature	
Rapid Shape	23 ± 3 °C
Asiga	35 ± 3 °C
Phrozen	23 ± 3 °C

3.2 Printing process

- Select the corresponding profile for FotoDent® denture in the printer software. Make sure that all software is up to date.

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• Contattare Dreve Dentamid GmbH per un elenco dei software e dell' hardware di processo convalidati.