Pontiform automix/ Pontiform direct mix

Ⓓ Produktbeschreibung

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Pontiform automix/ Pontiform direct mix

ⒼB Product description

M+W Pontiform automix/ Pontiform direct mix is a bis-acrylate based temporary crown and bridge material that is free from peroxide and methylmethacrylate. Stress to the pulp is avoided by the tissue-compatible low polymerization temperature. M+W Pontiform automix/ Pontiform direct mix offers special "Snap-Set" setting behavior and consequently has a long plastic working stage, during which the temporary restoration can be easily removed from the patient's mouth. Thanks to its excellent mechanical values, M+W Pontiform automix/ Pontiform direct mix is extremely fracture-resistant and color stable, and thus most suitable for fabricating long-term provisional restorations. M+W Pontiform automix/ Pontiform direct mix guarantees a highly accurate fit of the finished temporary restoration with optimal marginal integrity. M+W Pontiform automix/ Pontiform direct mix is mixed automatically.

Indication

- Extensive tooth defects: Tooth decay, severe wear, traumatic injuries of the dental enamel
- Shape and positional anomalies of the teeth
- Tooth loss

Contraindication

Do not use the material where there are contact allergies or known allergies to any of the components.

Intended purpose

Fabrication of temporary crowns, bridges, inlays, onlays, partial crowns, veneers and long-term provisional restorations

Limitations of use

Use the material to create temporary restorations only

Intended user

Dentist, dental assistant staff and dental technician

Patient target group

People treated in the course of a dental procedure.

Recommended use

A dental impression should be taken using alginate or silicone before preparing a stump for a crown or bridge or before a planned extraction.

The interdental flash should be cut away to improve the stability of the temporary restoration to be fabricated at a later stage. In the case of alginate impressions, dispense the M+W Pontiform automix/ Pontiform direct mix as soon as possible after taking the impression.

Elastic phase

Dry the prepared teeth and lightly coat the stumps, surrounding tissue and any existing synthetic material with a separating agent (e.g., Vaseline). M+W Pontiform automix/ Pontiform direct mix should be placed in the deepest part of the impression and then filled towards the gingival areas. Always keep the mixing tip immersed in the material to avoid air bubbles. M+W Pontiform automix/ Pontiform direct mix has a firm, elastic consistency approximately 2 to 3 min after mixing and can be easily removed from the patient's mouth. The curing process must be monitored intra-orally (e.g., with a probe), as the temporary restoration can only be easily removed during the elastic phase.

Final curing and processing

The final curing of M+W Pontiform automix/ Pontiform direct mix takes place in the next 3 min outside of the mouth. To this end, the temporary restoration should be placed back into the impression. After the removal of the oxygen-inhibited layer with a solvent (e.g., alcohol), the temporary restoration can be refined and polished. Commercially available temporary cements can be used for luting.

Repair of temporary restorations

A feature of temporary restorations fabricated using M+W Pontiform automix/ Pontiform direct mix is their high mechanical strength. If, however, a temporary restoration should fracture, we recommend the following procedure:

Fracture of the temporary restoration shortly after fabrication: Bond the fractures with freshly mixed

M+W Pontiform automix/ Pontiform direct mix.

Fracture of an older or in situ temporary restoration: Roughen the fractured sections and prepare undercuts. Bond the prepared fracture with freshly mixed Bis-acrylic temporary crown & bridge material. We recommend using a composite bonding agent to improve the adhesion. Press the parts together for approx. 3 min.

Timing

Injection into the impression	0 s
Insert the impression into the mouth	0 – 45 s
Removal from the mouth	120 – 180 s
Removal of the oxygen-inhibited layer with solvent.	360 – 420 s

Notes for use

If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

Interactions

- Residues of methacrylate-based materials may alter the setting behavior of silicone impression and bite registration materials.
- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.

Using the double-barreled syringe

The double-barreled syringes and cartridges enable the automatic dispensing and mixing, which allows a time-saving direct application. When using a new double-barreled syringe, remove the transport cap first. The mixing tip can then be fitted. After use and until the next treatment, leave the mixing tip on the double-barreled syringe as a seal. The used mixing tip should be replaced with a new tip when the device is used for the next treatment.

Note: When using a double-barreled syringe or cartridge for the first time, discard a small amount of the initially ejected material (e.g., a pea-sized amount).

Cartridge handling

The double-barreled syringes and cartridges enable the automatic dispensing and mixing, which allows a time-saving direct application. When using a new cartridge, remove the transport cap first. The mixing tip can then be fitted. After use and until the next treatment, leave the mixing tip on the cartridge as a seal. The used mixing tip should be replaced with a new tip when the device is used for the next treatment. **Note:** When using a double-barreled syringe or cartridge for the first time, discard a small amount of the initially ejected material (e.g., a pea-sized amount).

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

Warnings/precautions

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

Compressive strength: ≥ 200 MPa
Flexural strength REF000 290: ≥ 50 MPa
Flexural strength REF000 216: ≥ 60 MPa
Flexural strength REF000 214: ≥ 60 MPa

Storage/disposal

- Store at room temperature (15 – 25 °C/59 – 77 °F).
- Do not use after the expiration date.
- Leave the used mixing tip on the cartridge or syringe as a cap after application!
- Disposal must comply with national regulations.

Composition

Dental glass, EBPADMA, aliph. Polyester diurethanes, unsaturated polyester resin, SiO₂, additives.

Packaging

REF 000 290 1 Double-barreled syringes (A2) 5 ml, 5 Mixing tips
REF 000 214 A3, 1 Cartridge 50 ml, 15 Mixing tips
REF 000 216 A2, 1 Cartridge 50 ml, 15 Mixing tips

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