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ENGLISH

General description and purpose of the medical device:

Arti-BruX® is a red indicator paint which is applied with a brush (single use product) to the occlusal surfaces of restorative and prosthetic appliances as well as hard splints. After drying, Arti-BruX® forms a thin and even lacquer coating. The lacquer coating is intended for short-term use (up to max. 24 hours). It is saliva resistant and is abraded by tooth contacts (occlusal surfaces).

In addition to the known occlusion test materials (articulating papers, foils, sprays etc.), which only present a momentary view of the occlusal contacts, Arti-BruX® can be used as an Occlusogram that shows contacts and movement traces of teeth. Contact surfaces as well as traces of movement can be clearly identified by means of the removed lacquer. A short-term occlusion test (max. 24 hours) is recommended for patients with parafunction. **For professional use only.**

Indication

Arti-BruX® BK 89 is a red liquid contact color:
• To control active dynamic occlusion on restorative and prosthetic restorations as well as on bite splints.
• For short-term control of occlusion (max. 24 h) - Occlusogram.
• For visualization of enamel fissures as well as margin fissures of fillings and crowns.

Contraindication

• Arti-BruX® consists mainly of concentrated alcohol (ethanol 96%). Arti-BruX® should not be applied directly invasively to patients who are addicted to alcohol, pregnant women or children. Also contact with injured mucosa should be prevented.
• Arti-BruX® must not be used in case of iodine allergy.
• Arti-BruX® must not be applied to open dentin or prepared teeth.

Side Effects

• With insufficiently polymerized acrylics or vacuum splints, made of different thermoplastics, discoloration of the plastic may occur.
• If the surfaces to be tested are not properly polished or sand blasted, the lacquer can penetrate into micro cracks of the plastic or ceramic surfaces and cause discoloration.
• When Arti-BruX® is applied, the color can penetrate into enamel cracks or margins of fillings and cause temporary or even permanent discoloration.
• After longer use of Arti-BruX® (up to max. 24 hours) discoloration of the saliva may also occur; the lacquer consists of a natural resin and red food coloring. After use, the mouth should be rinsed with water.

Notes

• The product must be stored in the sealed bottle at temperatures between 0°C and 40°C.
• The product contains a solvent (ethanol 96%) which can evaporate quickly. After use, the bottle must be closed immediately.
• The solvent contained in the product is highly flammable. Do not use the product near ignition sources such as Bunsen burners (see also point Labeling / Symbols).
• Always apply the color to the surface to be tested with a single use brush (for example, Master-Brush™) to prevent cross-contamination. Arti-BruX® can also be applied invasively. The occlusal surfaces need to be dried.
• Arti-BruX® can be diluted with concentrated, medical alcohol (ethanol 96%) to change the color and viscosity.
• If the product is has been diluted too much, leave the residue to thicken in the provided single-use dappen dish.
• The set contains components (Master-Brush™ application brush, dappen dish) that are for single use only and must not be reused.
• Do not reuse application brush and dappen dish after use to avoid transfer of germs.
• The brush holder (Multi-Hold™) can be reprocessed according to the instructions (Refer to point „Reconditioning“).

Package Content

• 1 bottle Arti-BruX®, content 15ml
• 10 plastic dappen dishes (single use product)
• 1 holder for brush (Multi-Hold™ for Master-Brush™)
• 40 brushes Master-Brush™ (single use product)
• 1 goat hair polishing brush for dental laboratory handpiece
• Instructions for Use

Preparation

• Before each use clean, disinfect and sterilize the brush holder (Multi-Hold™). This applies in particular for the first-time use, because the supplied product set is delivered non-sterile. Please also note the detailed requirements for reprocessing under point „Reconditioning“.
• If necessary fix the brush in the brush holder.
• Before use, shake the lacquer bottle (mixing beads are audible), so that the color is evenly mixed. If the lacquer is too thick, it can be diluted at any time with medical alcohol (ethanol 96%, undenatured).

Application

• Arti-BruX® is used, taking into account usual practice hygiene requirements.
Lacquer application (and start of application)
• Clean and dry the surfaces to be tested (no saliva should be on the occlusal surfaces). If the occlusal contacts or tracks are not sufficiently visible, repeat the cleaning and drying, as well as the color application.
• Angle the application brush (e.g. Master-Brush™) at the tip as required to facilitate invasive application.
• Apply the lacquer with the application brush to the occlusal surfaces to be tested.
• Do not reuse the application brush and the dappen dish to prevent transmission of germs!
• To control bruxism, especially sleep bruxism, the patient should insert the splint shortly before going to bed.

Effect/Diagnosis

Check the surfaces to be tested after the end of the observation period. The coating is rubbed off by tooth contacts (occlusal surfaces). These contact patterns are used to check the occlusal contacts in static and dynamic. A uniform wear, e.g. on bite splints, is an indication that the individual track is evenly loaded. Strong one-sided wear indicates parafunctions (e.g. bruxism) or an incorrectly adjusted splint. The same applies to prosthetic restorations.

Cleaning (of the surfaces to be tested)

After successful inspection, remove the coating from the surface. Appropriate tools are: Tooth brush and appropriate tooth paste, enclosed polishing brush with a laboratory hand piece; use only with low speed and low pressure on splints, to avoid damage to the acrylic, alcohol bath or wadding soaked with alcohol. The lacquer dissolves well with concentrated alcohol (Ethanol 96%, undenatured). Splints can also be polished with commercially available polishing agents for the dental branch (e.g. www.hatho.de). Prosthetic restorations can be cleaned with a toothbrush or with standard equipment for professional dental cleaning.

Reprocessing

In principle, the following should be observed: „Requirements for hygiene in the processing of medical devices Recommendation of the Commission for hospital hygiene and infection prevention (KRINKO) at the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Devices (BfArM)“, Bundesgesundheitsblatt 2012-55: 1244-1310.

- The instrument must be reprocessed immediately after each treatment (within a maximum of 2 h).
- Thoroughly clean the brush holder by hand with a small brush with firm bristles.
- Note: The cleaning should be done in a water bowl without further additives below the water surface in order to achieve adequate cleaning of the instrument parts whilst avoiding protein fixation as well as to protect the environment from contamination due to splash water.
- Intermediately rinse the brush holder with water (minimum drinking water quality).
- Place the instrument in a standard cleaning and disinfection bath. Examples: Becht Premium Konzentrat, Dürr Dental ID 213 Instrument disinfection, Pluradent instrument bath PLULINE, Schülke & Mayr gigasept™ instru AF (see also "List of disinfectants and disinfection methods tested and approved by the Robert Koch Institute" or disinfectant list VAH).
- Note: Formaldehyde containing cleaning and disinfecting agents may only be used after sufficient cleaning in order to avoid protein fixation.
- Note: The instructions for use of the manufacturer of the cleaning and disinfecting agent must be followed carefully. In particular, the concentrations to be used and the reaction times must be observed!
- Note: For machine cleaning, the manufacturer's instructions for use must be followed carefully!
- Visual inspection for corrosion, damaged surfaces, chipping, mold damage and impurities. Damaged instruments must be discarded (For a limited restriction of reuse, see the chapter „Reusability“). If residual contamination still exists, the entire cleaning process must be repeated with all steps (Cleaning, intermediate rinsing, disinfection, final rinsing and drying).
- The instrument must be free of any residue and dry before further processing.
- Maintenance of the instruments is not required.
- The instrument is packaged and sealed in a sufficiently large single sterilization bag according to EN 868-2ff ISO 11607 (suitable for steam sterilization). Observe the instructions of the manufacturers of the sterilization bags and the sealing machines and the current normative requirements.
- Sterilization must be carried out in a validated process with moist heat in an autoclave according to DIN EN 13060 Type B or DIN EN 285 resp. ANSI AAMI S179. Observe the manufacturers operating instructions of the autoclave.
- Sterilize the instrument with humid heat (saturated water vapor) using a pre-vacuum procedure for 5 minutes at 134°C.
- The recommended storage life for sterile medical devices is described in standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling.

Reusability

• Frequent reprocessing restricts the life of the plastic instrument. The instrument must be disposed of in case of cracks in the plastic.
• It is not permitted to reuse single use items (applicator (Master-Brush™) as well as single use dappen dish). Safe use cannot be guaranteed if these products are reused, since there is an infection risk and / or the safety of the products is not liable.
• The use of damaged or soiled instruments is the responsibility of the user.
• In case of non-compliance, all liability is excluded.

Disposal Instructions

Arti-BruX® consists of a red food color and a natural resin, which is also approved for food and pharmaceutical applications. All components of the kit can be disposed of with the usual practice waste.

Notification of Incidents

Serious adverse events occurring in connection with the product must be reported to the manufacturer and the competent authority of the Member State.

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H225 Highly flammable liquid and vapor
P210 Keep away from heat, hot surfaces, sparks, open flames and other sources of ignition.
No smoking.

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H319 Causes severe eye irritation.
P305 IF IN EYES.
P351 Remove cautiously with water for several minutes.
P338 Rinse contact lenses, if present and easy to do. Continue rinsing.

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