Permaplast LH One Viscous Flow

Stand: 03/2026
50022/260115

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GB Instructions for use

Intended purpose

Composite-based dental filling materials are used to build up and maintain tooth substance.

Product description and User

Product Description

M+W Permaplast LH One Viscous Flow is a light-curing, flowable, highly radiopaque (210% Al) composite with high viscosity. The guidelines and specifications of DIN EN ISO 4049 apply.

M+W Permaplast LH One Viscous Flow composite is available in screw syringes and cavilfiles. The cavilfiles are for single use. Please do not reuse them since contamination and germ formation otherwise cannot be excluded.

Patient target group

All patients requiring tooth substance to be replaced or built up in order to restore or maintain dental function.

Users

For use in the dental practice by dental professionals

Composition

Glass powder, diurethane dimethacrylate, silicon dioxide, tetramethylene dimethacrylate

Filler content: 77% by weight (57% by volume) inorganic fillers (0.005 - 40 µm)

Indications

- Fissure sealing
- Extended fissure sealing on molars and premolars
- Fillings of class V according to Black (cervical caries, root erosions, wedge-shaped defects)
- Minimally invasive fillings of classes I, II and III according to Black
- Correction of dental enamel defects
- Blocking out of undercuts
- Minor shape and colour corrections on the enamel

Contraindications

If a patient has allergies or hypersensitivities to a component of this product, it should not be used or used only under the strict supervision of the attending physician/dentist.

Hazard and safety information

Contains: tetramethylene dimethacrylate, di-urethane dimethacrylate.

Warning: May cause an allergic skin reaction. Avoid breathing vapours/spray. Wear protective gloves. If skin irritation or rash occurs: Get medical advice/medical attention. Avoid release to the Environment.

Interactions with other materials

Phenolic substances (such as eugenol) inhibit polymerisation. Therefore, do not use any lining materials (such as zinc oxide-eugenol cements) which contain substances of this type. The dentist should consider known interactions and cross-reactions of the medical device with other materials already in the patient's mouth before using the product.

Application

Prior to the treatment, clean the hard tooth tissue with a fluoride-free polishing paste.

Cavity preparation

Preparation of the cavity with protection of the hard tooth tissue according to the general rules of adhesive technique. In the region of the anterior teeth, all enamel edges must be bevelled. In the region of the posterior teeth, by contrast, do not bevel the edges and avoid feather edges. Then clean the cavity with water spray, remove all residues, and dry it. Isolation is necessary. The use of a rubber dam is recommended.

Pulp protection/lining

Lining is not necessary if an enamel-dentin adhesive is used. If there are very deep cavities near the pulp, cover corresponding areas with a calcium hydroxide preparation

Approximal contact areas

For cavities with approximal sections, a transparent matrix must be created and fixed in place.

Adhesive system

Etching and bonding according to the manufacturer's instructions.

Application of composite

Application from syringes

Rotate the delivery tip clockwise until it locks in place on the syringe and dispense M+W Permaplast LH One Viscous Flow in thin layers (max. 2 mm) directly into the cavity. Use the curved delivery tips included with the kit for application. For hygiene reasons, tips are intended only for single use. Avoid introducing air bubbles during application with the tip. Ensure the prepared hard tooth tissue is thoroughly wetted.

Application from cavilfiles

Insert the cavilfil in the dispenser. Remove the cap. Secure the cavilfil such that the opening is aligned at the correct angle for application into the cavity. Introduce the material into the cavity using slow, even pressure. Do not use excessive force! The layer thickness must not exceed 2 mm. When finished, pull back the spindle to remove the cavilfil from the dispenser. Then remove the cavilfil. Note: For hygiene reasons, cavilfiles are intended only for single use.

Curing

Cure each layer after coating for 40 seconds using a commercially available polymerisation unit. The light guide must be held as close as possible to the surface of the filling. A dispersion layer forms on the surface during polymerisation which should not be touched or removed. It is used to bond the following composite layers.

Finishing

M+W Permaplast LH One Viscous Flow can be finished and polished immediately after polymerisation. Finishing diamonds, flexible discs, silicone polishers, and polishing brushes can be used for finishing. Check the occlusion and articulation and grind so that no early contacts or undesirable paths of articulation remain on the surface of the restoration.

Special information

- In case of time-consuming restorations, the surgical light should be temporarily moved further away from the work area to prevent premature curing of the composite or the material should be covered with an opaque foil.
- The curved application needles supplied with the product are for single use only, as contamination of the material and possible transmission of germs cannot be ruled out..
- A light polymerisation unit with an emission spectrum in the range of 350 – 500 nm is to be used for polymerisation. The physical properties required are achieved only with properly working lamps. For this reason, it is necessary to regularly check the light intensity according to the manufacturer's instructions.

Light intensity for curing $\geq 650 \text{ mW/cm}^2$

Wavelength for curing 350 – 500 nm

Curing time 40 sec.

Information on storage and handling

Store at 10 - 25°C (50 - 77°F). Close syringes tightly immediately after use. The material should be at room temperature before use. Retract the syringe spindle slightly after use to prevent the discharge opening from becoming clogged.

Shelf life

The maximum shelf life is printed on the label of each syringe or directly on the cavilfil. Do not use after the expiry date.

Side effects

With proper preparation and use of this medical device, adverse effects are extremely rare. However, immune reactions (such as allergies) or local discomfort cannot in principle be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case.

Disposal

Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

Troubleshooting

Problem	Cause	Remedy
Composite does not cure	Light output of the light-curing lamp is inadequate	Check the light output Clean light guide if soiled. Replace the light source, if necessary
	Emitted wavelength range of the light-curing lamp is inadequate	Consult manufacturer of the polymerisation lamp. Recommended wavelength range: 350 - 500 nm
Composite appears to be too hard and solid in the syringe	Material stored for a longer period of time at temperatures < 10°C	Allow composite to warm up to room temperature before use
	Syringe not closed correctly, composite polymerised	Close the syringe correctly with the cap each time after removing composite
Composite does not correctly cure (dark or opaque shades)	Composite layer too thick per curing cycle	Maintain max. layer thickness of 2.0 mm per layer
Restoration appears too yellow in comparison to the colour reference	Inadequate polymerisation of the composite layering	Repeat exposure cycle several times, min. 40 sec.

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M+W Dental Müller & Weygandt GmbH · Reichardsweide 40 · 63654 Büdingen/Germany · Tel.: +49 (0) 60 42 - 88 00 88
MANI MEDICAL GERMANY GmbH · Hertha-Sponer-Straße 2 · 61191 Rosbach v.d. Höhe/Germany