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Intended use

Bite registration material based on vinylpolysiloxane.

Indication for use

StoneBite® is an addition-curing silicone for bite registrations.

StoneBite® is indicated for:

- bite registrations
- key for intraoral registration

Patient target group

Patients who need dental care.

Warning

StoneBite® should be used solely on the occlusal surface of the teeth. Do not use StoneBite® for conventional impressions. Due to the high end hardness the removal is difficult and may lead to inadvertent extraction of teeth.

Description of material

StoneBite® is an addition-curing silicone for bite registration with high final hardness. The high final hardness enables a problem-free milling or grinding. In case of breaks, clearly defined pieces will be obtained. The material is easily injectable, even through mixing canulas with a thin diameter.


Technical Data according to DIN EN ISO 4823

Consistency	type 2
Working time	30 ± 5 sec. (23 °C / 73.4 °F)
Intraoral working time (incl. setting time)	≥ 80 sec. (37 °C / 98.6 °F)
Elastic recovery	> 99 %
Linear dimensional change	< 0.2 %
Final hardness	48 ± 2 Shore D

The processing time refers to a room temperature of 23 °C ± 1 °C (73.4 °F ± 1.8 °F) and a relative air humidity of 50 %. Lower temperatures prolong, higher temperatures accelerate the processing and setting times. Technical data according to DIN EN ISO 4823.

Processing

StoneBite® is exclusively available in the S50-double cartridge system. By means of an Injector and mixing tips an ideal 1:1 mixing of the components will be reached without losing time. With the aid of application tips StoneBite® can be applied properly and without loss on the teeth rows. Due to the thixotropy, undesired flow of StoneBite® into dental isterices is prevented. This working step is to be carried out without interruption, as the mouth temperature leads to an accelerated curing. Ensure that the bite is brought into the maximum intercuspidation within the given processing time of 30 sec. and keep it closed for further 50 sec. The fast polymerisation evokes a slight, physiologically harmless temperature increase. After curing the registrate may be removed from the patient’s mouth easily and without deformation.

Disinfection

Impression – After removal from the mouth rinse the bite registration under running water. Use a 2 % glutaraldehyde for 15 min. or another appropriate solution per manufacturer’s instructions. Disinfection procedure will not alter the impression.

Mixing canula and application tip – After use leave the mixing canula on the double cartridge, disinfect it by wiping the mixing canula with a disinfectant and do not exchange it before the next application. It serves as a sealing cap. The mixing canula and application tip are for single use only and must be replaced before the next application.

DS 50 Injektor – DS 50 Injektor can be used several times if properly disinfected. Disinfection can be carried out using a commercially available dipping disinfectant (e.g glutaraldehyde) according to the manufacturer’s instructions.

Cautions

- Due to technical filling reasons there might be an unequal transportation of the two components of new cartridges. For this reason it is recommended to dispose the initially mixed pea-sized amount.
- Remove residues of retraction agent with water spray before taking the bite registration.
- The curing may be disturbed by the use of latex gloves, hand creams, cleaning agents etc.
- Remove all residues of registration material from the patient’s mouth.
- The product is exclusively designed for the field of application described above and only to be used by qualified and trained personnel.
- After use leave the mixing canula as a cap on the double cartridge.
- Mixed material is chemically persistent. Spots onto clothing cannot be removed.

Until today, skin irritations or allergies are not known in connection with addition-curing silicone materials.

Keep away from children!

Only for dental use!

Contraindication

In case of allergies to any of the ingredients do not use StoneBite®. Not desired side effects of this medical device are not to be expected in case of proper processing and application. However, immune reactions (e. g. allergies) or irritations principally cannot

be excluded. Should you get to know non-desired side effects – also in doubtful cases – we appreciate in any case to be informed with a detailed description of the attendant circumstances and symptoms. We will follow each hint.

Serious incidents

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

Single-use advice

The hardened bite register can not be used for a further registration due to its physical features.

Disposal

Disposal of contents / container in accordance with local / regional / national / international regulations and safety data sheet.

Storage conditions

At room temperature, i. e. 18 °C / 64 °F–28 °C / 82.4 °F

Durability: see date of expiry.

StoneBite® is not to be used when expiry date has passed.

Form of delivery

D504872

StoneBite®	type 2
2 x 50 ml double cartridges	12 mixing tips
	12 application tips
StoneBite®	4 x 50 ml double cartridges

Mixing tips	40 pieces
Mixing tips	100 pieces
Application tips	40 pieces
DS 50 Injektor 1:1	

This product was developed for dental use and has to be applied according to the instructions. In the case of damages resulting from off-use or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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