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en Instructions for use

Patient target group
All persons for whom individual trays are required for treatment.

Intended user
Dentist

Composition Qualitative data:
Oligomeric urethane acrylate, vinyl ester resin, 2,2-dimethoxy-1, 2-diphenylethane-1-one or by 2,2-dimethoxy-2-phenylacetophenone, (+)-2,3-dioxo-1,7,7-trimethyl-bicyclo-(2,2,1)-heptane di-bornane-2,3 dione, diphenyl(2,4,6-trimethylbenzoyl) phosphinoxide, 2-hydroxy-2-methyl-1-phenyl-propane-1-one, 2-ethylhexyl-4-(dimethylamino)benzoate, 2,2-diethoxyacetophenone, soda-lime glass, amorphous silica, dyes (Pigment Red 146 = CI 12485 / Pigment Blue 29 = CI 77007 / D&C Red 30 = CI 73360 / Pigment Blue CI 74160 / Pigment Violet 19 = C.I. 73900), flavourings (Optamint).

Quantitative data:
In 100 g there are 19 g acrylates/methacrylates and 80 g glass.

Intended purpose of the product
M+W LC-Löffelplatten is a light-curing tray material in the form of prefabricated, baseplates to make rigid, custom

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impression trays and bite-blocks, as well as for other indications. M+W LC-Löffelplatten can be cured in all ultraviolet and LED light-curing units. The baseplates are supplied in a re-usable light-proof box that can be used for transporting or storing models, etc.

Indications
- Custom impression tray

Contraindications
The material contains methacrylates. Do not use in patients with an allergy to methacrylates or other substances contained in the product.

Performance characteristics of the product
M+W LC-Löffelplatten polymerises with light to form a solid and compact polymer.

Specific precautions
Avoid any contact between the eyes and skin and the uncured material, i.e. also with the inhibited layer. Always use a dust extractor when trimming the tray.

Preparation/Manufacture
- Mark the outline of the future tray on the model.
- Block out the undercuts on the model.
- To create a spacer for the impression material, block out

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the entire model (e.g. with a sheet of wax) and cut out stops.
- If the plaster is very dry, soak the model briefly in water or apply dye isolation or an alginate isolation.
- Remove the baseplate from the box and close it.
- Once the model has been blocked out, adapt the baseplate to it and cut to shape using, e.g. a scalpel.
- Use the remaining material to form a handle and attach it to the tray.
- Place the model and tray in a light-curing unit and cure without vacuum.
- Remove the tray from the model and polymerize once again from the underside.
- If the tray has been cured without a vacuum, use alcohol to remove the smear layer.

Times:
Working at daylight 20 minutes

Polymerisation in an ultraviolet or LED light-curing unit, depending on the unit in use:
The top of the tray, on the model 3-5 minutes.
The underside of the tray, off the model 3-5 minutes.

Trimming:
- Trim the custom tray with tungsten carbide burs and smooth it with sandpaper.

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Hygiene:
Immerse impressions taken with M+W LC-Löffelplatten customs trays in a disinfectant solution based on glutar aldehyde. The time depends on the disinfectant manufacturer's directions.

The finished medical product is used by the dentist, e.g. for individual trays with impression material.

Caution
The product has been developed for use in the dental sector and must be used according to the directions. The manufacturer accepts no liability for damage that results from use other than the above. Furthermore the user is obliged under his own responsibility to examine the material for its suitability and its usage possibilities for the intended applications, especially if the latter are not detailed in the directions for use.

Side effects
Adverse side effects are not known.

Interactions
Interaction of the product with other medical devices, medicinal products or other substances is not known.

Storage
The Product should be stored between 0° C - 22°C.

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Müller & Weygandt GmbH · Reichardsweide 40 · 63654 Büdingen/Germany · Tel.: +49 (0) 60 42 - 88 00 88