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REF 000 711 Inhalt/Content/Contenu/Contenuto/Obsah:

e 10 L

- Ⓒ GB SPRAY DISINFECTANT ALCOHOL-FREE
- Ⓕ F VAPORISATEUR DÉSINFECTANT SANS ALCOOL
- Ⓘ I DISINFETTANTE SPRAY SENZA ALCOL
- Ⓒ CZ SPREJOVÁ DEZINFEKCE BEZ OBSAHU ALKOHOLU

Stand:300120

LOT

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40°C

0124

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☞ SPRAY DESINFECTANT ALCOHOL-FREE

Type of application:

As per German Medical Product Law (MPG): Aqueous ready-to-use surface disinfection for prophylactic surface disinfection of non-invasive medical devices. Particularly suited for use on alcohol-sensitive surfaces such as acrylic glass.

According to the Biocide Directive: Aqueous ready-to-use surface disinfection for the disinfection of surfaces. Particularly suited for use on alcohol-sensitive surfaces such as acrylic glass.

Application: The ready-for-use solution is applied undiluted from a distance of approx. 20-30 cm as foam (using a suitable foam spray head) onto the surface to be disinfected and evenly distributed with a lint-free, undyed disposable wipe and

allowed to react. Ensure complete wetting. If necessary, remove excess quantity of the preparation with a disposable undyed wipe. Before using on large surfaces with sensitive materials, for example, artificial leather, test first for material compatibility on a small hidden section. To avoid streaking, the treated surfaces should be wiped at regular intervals with a moist (water) disposable wipe. The disinfectant solution may only be used for soaking wipes that are compatible with the solution according to EN 16615.

Exposure times of the impregnating liquid:

Bacteria according to EN 13727, EN 1276 and EN 13697 (high load): Conc. 1 min

Candida albicans according to EN 13624, EN 1650 and EN 13697 (high load): Conc. 1 min

Enveloped viruses (i.e. HBV, HIV, HCV) according to RKI/DVV with high load: Conc. 1 min

Rotaviruses according to DVV without load: Conc. 1 min

Noroviruses according to EN 14476 (low load): Conc. 30 min

Exposure times for cloth wipe systems:
Polyoma-SV 40 Viruses (RKI/DVV with contamination): Conc. 30min

Bacteria and Candida albicans according to VAH with mechanics and EN 16615 (high load): Conc. 2 min

Bacteria and Candida albicans:
Entered in VAH list: Conc. 5 minutes
(shortest exposure time for surface disinfectants in VAH list)

100 g contain:

0,24g benzyl-C12-16-alkyldimethyl ammonium chloride
0,24g N,N-didecyl-N,N-dimethyl ammonium chloride
0,24g C12-14-alkyl[(ethylphenyl)methyl] dimethyl ammonium chloride
Only pass completely emptied containers to recycling. Dispose of residual material acc. to EAK 07 06 99
Directive 648/2004/EC: Contains: < 5% cationic surfactants, fragrances.

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