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Intended use

Bite registration material on vinyl polysiloxane basis.

Indications for use

StoneBite® scan is an addition-vulcanising silicone for bite registrations with high final hardness.

Description of material

StoneBite® scan is an addition vulcanising silicone (VPS) for bite registration with an extremely high final hardness. The high final hardness enables a problem-free milling resp. grinding. In case of breaks, clearly defined pieces will be obtained. The material is easily injectable, even through mixing canulas with a thin diameter. **StoneBite® scan** is reflection-active and can be used for optical collection of data (registration of antagonists) in CAD / CAM / CIM Technology.

Indication

- bite registrations
- registration key for the intraoral bearing device registration
- optical collection of data in CAD / CAM / CIM systems

Technical Data according to DIN EN ISO 4823


Consistency	type 2
Processing time	30 ± 5 sec. (23 °C)
Demoulding (incl. processing time)	≥ 80 sec. (37 °C)
Recovery from deformation	> 99,9 %
Linear change of measures	< 0,2 %
Final hardness	32 ± 2 Shore D

Processing

StoneBite® scan is exclusively available in the S50-double cartridge system. By means of an Injector and mixing tips an ideal 1:1 mixing of the components will be reached without loosing time.

With the aid of application tips **StoneBite® scan** can

be applied cleanly onto the teeth rows and without loss. Due to the high stability the material does not flow away. This working step is to be carried out without interruption, as the mouth temperature leads to an accelerated curing. Ensure that the bite is brought into the maximum intercuspidation within the given processing time of 30 sec. and keep it closed for further 50 sec. The fast polymerisation evokes a slight, physiologically harmless temperature increase. After curing the registrate may be removed from the patient's mouth easily and without deformation.

Disinfection

After removal from the mouth rinse the bite registration under running water. Disinfection for 15 min. in a 2 % glutaraldehyde or another appropriate solution.

Important advice

- Due to technical filling reasons there might be an unequal transportation of the two components of our cartridges. For this reason it is recommended to dispose the initially mixed pea-sized amount.
- The vulcanisation may be disturbed by the use of latex gloves, hand creams, cleaning agents etc.
- The product is exclusively designed for the field of application described above and only to be used by qualified and instructed personnel.
- Remove all residues of registration material from the patient's mouth.
- After use leave the mixing canula as a cap on the double cartridge

Keep away from children!

Only for dental use!

Storage conditions

At room temperature, i.e. 18 °C–28 °C.

Durability: see date of expiry.

StoneBite® scan is not to be used when expiry date has passed.

Form of delivery

REF D504870

StoneBite® scan

2 x 50 ml double cartridges
12 mixing tips
12 application tips

Accessories

REF D4970X1

REF D49701X1

REF D49603

REF D5007

Mixing canulas 40 pieces
Mixing canulas 100 pieces
Application tips 40 pieces
DS 50 Injector 1:1

This product was developed for dental use and has to be used according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according to suitability and possible use in self dependence before using it subject to its purpose, especially in case the purpose is not mentioned in the instructions.

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