DENTA-O-RESIN

Stand 290818

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© M+W Dental-O-Resin was developed specifically for fabricating orthodontic appliances.

If there is anything in this patient information leaflet that you do not understand, please contact our customer service department before using the product.

Directions of use

Indications

Orthodontic appliances produced in the sprinkling and kneading techniques.

Properties

- Available in clear
- outstandingly brilliant colors
- sprinkles very well soaks up the monomer well (when sprinkled on)
- acrylic will not slump whilst being built-up = rapid application
- high final hardness, not brittle! easy finishing, does not smear

Instructions:

The sprinkling method

Carefully apply a separating agent to the model. A high-quality alginate separator, should be used. Should the model be very dry, briefly soak it in water prior to applying the separator. Sprinkle the powder onto the model with the plastic spray bottle. Moisten the powder with KFO-monomer. Once the entire liquid has been soaked up and a homogeneous dough has been formed, apply the next layer. Build-up the desired form layer by layer. Once the entire appliance has been formed, soak the surface once more. The upper layer will quickly lose its shine and the model can be placed in the polymerization vessel.

Polymerization under pressure

The water in the pressure vessel should be at 37 °C (100 °F). The acrylic must be fully immersed in water. A pressure of 2.2 atmospheres must be applied and maintained for 20 minutes in order to polymerize the acrylic. When M+W Denta-O-Resin auto polymerizing acrylic is processed in this manner its density is so high that its mechanical properties correspond to those of heat curing acrylics. M+W Denta-O-Resin does not tend to form bubbles or porosities and exhibits a high precision of fit after being cured in a pressure vessel.

Instructions for the color concentrates

Obtain an individual color with help of the M+W Denta-O-Resin color concentrates (red, green, blue, yellow, 10 ml each). Either add the concentrates to the liquid or to the mixed product (3 - 5 drops, as desired).

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e.g. allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

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GDF GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany
Vertrieb durch Müller & Weygandt GmbH · Reichardsweide 40