Befestigungszement Perma

Stand:30/03/26
500266 1260327

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Instructions for Use

Intended purpose

Fixing composites/root pin cements are used for the adhesive, permanent bonding of restorative materials to restorative materials, to hard tooth structure or to a denture base, or to secure root pins in the root canal.

Description of product and users

Description of product

M+W Befestigungszement Perma is a dual-curing attachment composite for a natural tooth. The material is fluorescent and radiopaque and is used to permanently attach indirect restorations to the tooth. This prevents the restoration from becoming dislodged whilst at the same time filling the minimum gap between dental preparation (or implant abutment) and restoration.

Patient target group

M+W Befestigungszement Perma is used in patients requiring attachment of an indirect restoration to dentine, construction of a stump superstructure or securing of an endodontic post.

Users

For use in the dental practice by dental professionals.

Composition

Composition

Glass powder, diurethane dimethacrylate, prepolymerisable, tetramethylene dimethacrylate, silicon dioxide, dibenzoyl peroxide

Filler content

75% by weight, inorganic fillers (0.005 - 5,0 µm)

Indications

- Adhesive bonding of indirect constructions such as veneers, inlays, onlays, crowns and bridges. The indirect constructions may be made of vitreous ceramic, zirconia, composite, zirconia-based composite, stainless steel, non-stainless steel and titanium.
- Adhesive bonding of endodontic posts
- Stump superstructures (as core build-up material)

Contraindications

- If a patient has allergies or hypersensitivities to a component of this product, it should not be used or used only under the strict supervision of the attending physician/dentist.
- Application of the material is contraindicated if leaving to dry, or the specified application technique is not possible.

Warnings

Contains: di-urethane dimethacrylate, tetramethylene dimethacrylate, dibenzoyl peroxide.

Important: May cause allergic skin reactions. Harmful to aquatic life with long-lasting effects.

Safety instructions

Avoid breathing vapours/ spray. Wear protective gloves. If skin irritation or rash occurs: Get medical advice/attention. Avoid release to the environment.

Interactions with other medicinal products

- Phenolic substances (such as eugenol) inhibit polymerisation. Therefore, do not use any lining materials (such as zinc oxide-eugenol cements) which contain substances of this type.
- Oxidative disinfectants such as hydrogen peroxide may interact with the curing system and impair the curing process.

Application / Preparation

Notes: The material should be at room temperature during processing as it may be difficult to squeeze out and mix at refrigerator temperatures.

M+W Befestigungszement Perma should be processed quickly after removal from the Automix syringe and the structure inserted rapidly.

Preparatory work

Removal of temporary

Remove temporary, clean cavity or stump, remove residues of temp. cement, rinse thoroughly with water and dry with oil-free compressed air.

Trial of framework structure

Try on structure and check fit and contact points. Clean the objects to be veneered thoroughly with alcohol (e.g. isopropanol) and then dry with oil-free compressed air. Avoid repeated ingress of moisture or contamination. Exercise extreme caution when checking occlusion of brittle ceramic workpieces, as otherwise there is a risk of cracking when in an unluted state.

Adjustments may be made using fine diamond-tip drills at moderate speed and light pressure if required. Milled surfaces should be repolished. Relative drying of the work area is essential when composites are used for adhesive fixing, preferably using wadding and dry angles. Absolute drying, e.g. by using a rubber dam, is recommended.

Note: Alcohol must NOT be used to clean PMMA materials. They should be cleaned, for example, with disinfectants based on quaternary ammonium compounds.

Pre-treatment of natural tooth

A coating of 37% phosphoric acid etching gel is applied to the preparation surface in accordance with the manufacturer's instructions. Next, rinse off the etching gel thoroughly with water. Wet the surface with a suitable dual-curing dental adhesive system in accordance with the manufacturer's instructions. We recommend using a micro-brush to apply the adhesive. Remove water and solvent residue with a gentle air stream.

Pre-treatment of structure

The framework manufacturer's instructions will generally apply.

Pre-treatment of	Procedure
Metals	<ul style="list-style-type: none"> Sand blast the inside surface of the restoration with 30 - 50 µm aluminium oxide at max. 2 bar Cleaning with alcohol and drying of objects Apply an even coating of a suitable luting agent with a thin brush in accordance with the manufacturer's instructions. <p>Note: In order to achieve optimum bonding, do not clean with phosphoric acid.</p>
Lithium disilicate glass-reinforced ceramic	<ul style="list-style-type: none"> Etch with hydrofluoric acid gel according to manufacturer's instructions or in accordance with restoration manufacturer's instructions. Rinse off the restoration thoroughly with water and dry with oil-free air Apply an even coating of a suitable luting agent with a thin brush in accordance with the manufacturer's instructions.
Composite	<ul style="list-style-type: none"> Sand blast the inside surface of the restoration with 30 - 50 µm aluminium oxide at max. 2 bar Cleaning with alcohol and drying of objects Apply an even coating of a suitable adhesive system with a thin brush in accordance with the manufacturer's instructions.
PMMA materials	<ul style="list-style-type: none"> Sand blast the inside surface of the restoration with 30 - 50 µm aluminium oxide at max. 2 bar Clean with suitable agent, e.g. disinfectant based on quaternary aluminium compounds and drying of objects Apply an even coating of a suitable luting agent with a thin brush in accordance. <p>Note: Avoid contact with alcohol as this can cause fracturing (craquelure).</p>

Attachment of indirect restorations

Application of M+W Befestigungszement Perma

Squeeze out a small amount of paste and discard to balance the catalyst/base paste ratio in the Automix syringe. Attach a new disposable Automix tip to the dual chamber syringe. The initial quantities of mix from the syringe should not be used for luting. Squeeze out the M+W Befestigungszement Perma from the Automix syringe and apply the desired quantity directly into the restoration and spread evenly. Since the attachment material hardens in the mixing tip, the tip can be used to seal the contents of the syringe until next used.

Insertion of restoration and removal of surplus

Apply light pressure to bring restoration in situ and fix/hold. Remove surplus attachment composite immediately using a disposable brush, foam pellet, dental floss or an implant scaler. Pay particular attention to ensuring that the surplus is properly removed from difficult to reach areas (approximal, gingival margins). Light-cure composite or vitreous ceramic constructions for 20 sec on each side after surplus is removed. Chemical curing is complete after 6 min. Please light-cure all margins again for 20 sec immediately after surplus is removed. Formation of an oxygen inhibition layer can be prevented by coating restoration margins immediately afterwards with glycerine gel / airblocker and rinsing off with water after curing is complete.

Preparing the finished restoration

Remove the retraction cord if applied. Check occlusion and functional movement and adjust as required. Polish restoration margins and structure.

Attachment of endodontic posts

Prepare the root canal and cavity as usual. Clean the root canal with a 3% solution of hydrogen peroxide (H₂O₂) or 2.5 - 5.25% solution of sodium hypochlorite (NaOCl). Next, immediately rinse with water and dry with paper points.

Important: Since residues of oxidative disinfectants such as hydrogen peroxide interact with the material and may prevent it from functioning properly, ensure that it is completely removed.

Select a suitable post and try it on. Clean the post with alcohol. Wet the post with a suitable dual-curing dental adhesive system in accordance with the manufacturer's instructions. A coating of 37% phosphoric acid etching gel is applied to the preparation surface of the cavity (coronal and radicular portion) in accordance with the manufacturer's instructions. Next, rinse off the etching gel thoroughly with water and remove all acid residues. Aspirate the water and dry the cavity with paper points. To keep the dentine moist, do not blow dry in air. A suitable dual-curing dental adhesive system is applied to the root canal on the occlusal surface of the stump with a thin micro-brush in accordance with the manufacturer's instructions. Surplus is removed with paper points.

Important: Ensure that the micro-brush reaches all areas of the canal and the bonding agent is evenly spread throughout. The handle of the micro-brush must not come into contact with the tooth surface and the brush must not be caught on the surface.

Remove water and solvent residue with air. Position the post to check canal depth and press adhesive further down into the dentine tubules. Then remove the post. Attach a new disposable Automix tip to the M+W Befestigungszement Perma dual chamber syringe. The initial quantities from the syringe should not be used for luting. Fill the canal from the bottom of the canal to the coronal part with M+W Befestigungszement Perma; do not use Lentulo spirals or other rotating instruments. Wet the surface of the post with sufficient M+W Befestigungszement Perma and position the post fully into the canal, letting surplus attachment composite drain out. Fix post and carefully remove surplus cement. Once introduced, light-cure for 40 sec. Restoration can now begin.

Stump superstructure

Etch the coronal surface with 37% phosphoric acid etching gel in accordance with the manufacturer's instructions. Rinse off the etching gel thoroughly with water and remove all acid residues. Aspirate water and dry the surface carefully with a paper point. To keep the dentine moist, do not blow off in air.

We recommend using a micro-brush to apply the adhesive. Remove water and solvent residue with a gentle air stream. Coat with the light-curing dental bonding and light-cure in accordance with the manufacturer's instructions. Place a new mixing tip on the M+W Befestigungszement Perma dual cartridge and apply a sufficient quantity around the post. Fill the stump shaper, if used, and place over the post. Remove surplus and then light-cure for 40 sec. Chemical curing is complete after 6 min. Finish and polish, the impression can then be taken.

Notes on processing

Processing time

Processing and setting times depend on the ambient temperature. The following times apply once M+W Befestigungszement Perma is removed from the Automix syringe:

	Room temperature approx. 21°C (approx. 69°F)	Intraoral 37 °C (98.6 °F)
Processing time	3-4 min	approx. 2 min
Curing time, incl. processing time	approx. 7-9 min	approx. 3-4 min

Light polymerisation

Light intensity	Lightwave range	Exposure time per unit of area
> 500 mW/cm ²	350 - 500 nm	20 sec.

Note: Please light-cure all margins again for 20 sec after surplus material is removed.

Mixing ratio

Use of the Automix syringe including mixing tip ensures that M+W Befestigungszement Perma is always mixed at the optimum ratio of 1:1.

Troubleshooting / FAQ List

Problem	Cause	Remedy
Attachment does not grip after attachment composite is applied	Bonding surface contaminated	Make sure that the working environment is free of contamination after cleaning with alcohol and drying
	Components of the attachment composite are not evenly mixed (air bubbles)	Use a mixing tip and discard the first 0.25 g or so of the material squeezed out of the tip
	Curing of attachment composite takes longer than processing time	Keep within processing time and only mix small quantities
	Objects move during curing time for attachment composite	Keep within curing time and do not touch objects, fix in place if necessary
Fixing composite does not cure	Contamination of bonding surfaces with products containing Eugenol	Use products that are free of Eugenol
	Attachment composite stored incorrectly (too warm)	Purchase new attachment composite and store in refrigerator
	Unequal quantities of the two components mixed in	Use original mixing tips for the syringe and discard the initial material squeezed out

Information on storage and handling

M+W Befestigungszement Perma must be stored in a cool place. The storage temperature is 3 - 9 °C (37 °F - 48 °F). To close the automix syringe after use, leave the used mixing tip plugged in.

Shelf life

The maximum shelf life is printed on the label of each syringe. Do not use after the expiry date.

Warnings on side effects

With proper preparation and use of this medical product, adverse effects are extremely rare. However, immune reactions (such as allergies) or local discomfort cannot in principle be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case.

Instructions for the consumer

Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

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