

3M | Astringent Retraction Paste

ESPE





- (de)** Adstringierende Retraktionspaste
- (fr)** Pâte astringente pour l'évasement sulculaire
- (it)** Pasta astringente per retrazione
- (es)** Pasta de retracción astringente
- (pt)** Pasta de Retracção Adstringente
- (nl)** Adstringerende retractiepaste
- (el)** Στυπτική πάστα απώθησης
- (sv)** Astringerande retraktionspasta
- (fi)** Adstringoiva retraktiopasta
- (da)** Adstringerende retraktionspasta
- (no)** Adstringerende retraksjonspasta

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en ENGLISH

en
ENGLISH**Product Description**

Astringent Retraction Paste is an astringent retraction paste for any situation in which a temporary displacement of the marginal gingiva and/or the assurance of a dry and clean sulcus is required. The acidic paste containing aluminum chloride is applied directly from the retraction capsule into the sulcus. A hemostatic effect results from the compression on the gums by the paste containing aluminum chloride as an astringent. This hemostatic effect quickly closes any small blood vessels which might be unintentionally opened. Depending on the clinical situation and working technique of the person performing the treatment, the paste can be used as an alternative to or in combination with retraction cords or other retraction methods. Any suitable, commonly available dispensers, e.g., the Filtek™ Restoratives Dispenser can be used for application of the retraction paste. A retraction capsule contains sufficient material for up to 3 teeth, but may not be used on more than one patient.

 These Instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly readable. For details on all additionally mentioned products please refer to the corresponding Instructions for Use.

Intended Purpose

Intended purpose: astringent paste for temporary displacement of the marginal gingiva and for hemostasis and moisture control.

Intended users: educated dental professionals, i.e. general dentists, dental assistants, dental hygienists, who have theoretical and practical knowledge on usage of dental products.

Clinical benefit: temporary displacement of the marginal gingiva and creation of a dry and clean sulcus as required for dental procedures.

Indications

- All indications for temporary displacement of the marginal gingiva and to provide for a dry sulcus when the periodontium is healthy, such as:
 - Taking impressions, either with impression material or as digital impression
 - Cementation of temporary and permanent restorations
 - Preparation of Class II and V fillings

Contraindications

The astringent retraction paste must not be used in patients suffering from a diseased periodontium, open furcations, or exposed bone.

Precautionary Measures

The user must wear protective goggles during treatment with the astringent retraction paste. If eye contact occurs, wash immediately with plenty of water and consult a physician if necessary.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

3M Safety Data Sheets/Safety Information Sheets can be obtained from www.3M.com or contact your local subsidiary.

Application

- ▶ Remove a retraction capsule from the blister (immediately before use only), place it in a suitable, commonly available dispenser, and check to see that the retraction capsule is firmly seated in the dispenser.
- ▶ Remove the sealing cap from the capsule tip (immediately before application only) and discard.
 - The capsule tip is rounded and must not be shortened because this would result in sharp edges.
- ▶ Dispense a small quantity of paste onto a mixing block and discard (Fig. 1).
- ▶ Rinse the sulcus thoroughly with water and dry lightly with air before application of the paste.
- ▶ Introduce the capsule tip into the sulcus, opening the sulcus (Fig. 2).
- ▶ Move the capsule tip in the sulcus slowly and evenly around the tooth, pressing out the paste and filling the sulcus all around with enough retraction paste so that an excess appears (Fig. 3).
 - If required by the clinical situation, the retraction paste can be used in combination with a cord to open the sulcus further (Fig. 4).
- ▶ Allow the astringent retraction paste to react for at least 2 minutes in the sulcus (Fig. 5). Keep all moisture away from the site during this period.
 - An overly long residence time may cause damage to tissue.
- ▶ Remove the retraction paste completely from the sulcus after the residence time, using a mixture of air and water and with the aid of a suction device (Fig. 6).
 - The patient must not swallow the retraction paste.

Concluding the Treatment

- ▶ Carefully check the sulcus of the treated teeth and the surrounding areas. Remove any astringent retraction paste from the mouth.

Notes

Astringent Retraction Paste is a single-use device.

After cleaning, disinfection and sterilization the paste would not work any longer as intended.

Disposal

Dispose of the contents or container in accordance with the applicable regulations. Please pay special attention to the disposal of contaminated waste to avoid health risks because of improper handling. Dispose of empty retraction capsules in the normal waste of the medical practice.

Storage and Stability

Store the product at 15–25° C/59–77° F.

Do not use after the expiration date.

Customer Information

No person is authorized to provide any information that deviates from the information provided in this instruction sheet.

Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Symbol Glossary

Reference Number and Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.1.1 Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
ISO 15223-1 5.1.3 Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.4 Use-by date		Indicates the date after which the medical device is not to be used.
ISO 15223-1 5.1.5 Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6 Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.3.7 Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed.

Reference Number and Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.4.2 Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure.
ISO 15223-1 5.4.4 Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
CE Mark		Indicates conformity to European Union Medical Device Regulation or Directive with notified body involvement.
Medical Device		Indicates the item is a medical device.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a dental professional.
Recycle		Indicates packaging is recyclable with mixed paper.)
Low-density polyethylene		Indicates plastic component made of low-density polyethylene is recyclable.
High-density polyethylene		Indicates plastic component made of high-density polyethylene is recyclable.
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

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3M, ESPE and Filtek are trademarks of 3M or 3M Deutschland GmbH.
Used under license in Canada. © 2019, 3M. All rights reserved.
U.S.A. 1-800-634-2249.
3M Deutschland GmbH. Utilisées sous licence au Canada. © 2019, 3M.