# POWER WHITENING XTRA 38% (MIXED 30% HP)

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**Instructions for POWER WHITENING XTRA 38% (MIXED 30% HP)**

**Introduction:** POWER WHITENING XTRA 38% (MIXED 30% HP) is a dualsyringe system to whiten the teeth in the dental office quickly and easily – normally in a single session. POWER WHITENING XTRA 38% (MIXED 30% HP) contains 38 % Hydrogen Peroxide ( $H_2O_2$ ) - mixed 30%. The gel is highly viscous to prevent dripping. It can be applied easily to the teeth using the mixing tips. The mixed gel is light red and small air bubbles show the reaction during whitening. The intensity of the red colour has no influence on the bleaching result. The used syringe can be re-closed allowing the remaining gel to be used for the next patient, with a new mixing tip.

**Indications:** POWER WHITENING XTRA 38% (MIXED 30% HP) serves for quick, medically indicated brightening of strongly discoloured vital and devitalised teeth, whose colour cannot be returned to a natural or normal tooth colour through the professional cleaning of teeth. Medically indicated brightening of teeth applies to discolouration which leads to a deterioration (health-related) in the patient's quality of life due to the abnormal deviation in the colour of teeth. If necessary, professional bleaching of teeth can also be employed as an alternative to veneering (significantly more invasive) of the teeth concerned with veneers or crowns. Tooth colours which are common in the population in terms of brightness, translucency and colour shade are regarded as being natural or normal. Tools for establishing the diagnosis and indication, for example, include the Vita A1-D4 classical colour scale (for the most common natural tooth colour shades) and the separately available questionnaire (OHIP-G19) to determine oral health related quality of life. Examples for the causes of such (intrinsic) discolourations include: tetracycline / minocycline discolouration, moderate to severe dental fluorosis (TF Index > 2), devitalized teeth (pulp necrosis) - i.e. due to trauma, genetic or age-related (degenerative) discolouration, porphyria or icterus, pronounced food deposits with penetration into the tooth substance (intrinsic discolouration through extrinsic dyes). It is used to brighten one or more teeth or parts of a tooth. It is also suited for bleaching individual, endodontically treated devitalized teeth including intracoronal bleaching.

To validate the indication and to determine the success of therapy a special questionnaire to examine the oral health related quality of life (MLQ) of your patients is available. For this purpose by the patient weighted questions ("never"=0, "hardly"=1, "from and to"=2, "often"=3 and "very often"=4), are added for a total value. Theoretically any value between 0 result here(=no impaired quality of life) to 56 (=maximum impairment of MLQ) The questionnaire should be filled by the patient

GB

alone, possibly in the waiting room- in reach case before and 2 weeks after treatment. This questionnaire is only for assistance for your diagnosis; a need for treatment because of a certain value cannot be derived. Wherein a sum of a value of > 10 and a tooth colour of B4-C4 or darker or a difference of more than 8 levels of brightness to the adjacent teeth is generally indicated a change in tooth colour. The questionnaire can be obtained free of charge from M+W Dental.

**Warning:** Because of the high concentration of Hydrogen Peroxide POWER WHITENING XTRA 38% (MIXED 30% HP) is highly reactive. To avoid contact with skin or eyes we recommend the use of gloves and safety goggles. If the gel does come into contact with skin, remove it and rinse immediately with water. White spots may appear in the affected area and you may experience a burning sensation. The symptoms should disappear after approx. 30 – 60 minutes. If the gel comes into contact with the eyes, rinse immediately and consult a doctor. Any dental work (crowns, fillings etc.) should be done 1 – 2 weeks after the treatment in order to match the colour of the newly whitened teeth. Please note that the adhesive strength of any dental work may be lessened in the 2 – 3 days following bleaching treatment. Do not perform whitening treatments under anaesthetic. The manufacturer does not recommend the use of Bleaching gel on patients who are pregnant, lactating or under the age of 18.


### 1.) Prior to treatment

A diagnosis for the indication of medically indicated tooth brightening must be made by the dentist prior to bleaching. After examining and advising the patient, you should determine and document the existing tooth colour shade. Maybe you should perform professional tooth cleaning, best 1-2 weeks prior to brightening treatment. Use a lip retractor. If required, also use cotton wads and apply lip care to the lips. **Important:** as restorative materials cannot be brightened, we recommend medical tooth bleaching to be performed prior to placing aesthetical restorations. To provide correct colour matching, you should perform possible restorations no earlier than 1-2 weeks after bleaching treatment.


## 2.) Application of Gingiva Protector

Attach the dispensing tip to the top of the syringe of the Gingiva Protector. Dry the gingiva carefully to ensure that the Gingiva Protector will hold and seal during the complete treatment. Place the Gingiva Protector directly on the gingiva and cover 3 – 4 mm of the gingiva and 0.5 – 1 mm of the teeth. Now light-cure the Gingiva Protector with a UV-lamp for 30 – 40 seconds in an arching movement over the material. Instead of the Gingiva Protector a rubber dam can be used.

## 3.) Whitening procedure

Attach the mixing tip to the top of the dual-syringe. Apply the gel directly to the surface of the teeth (to whiten single teeth, Gingiva Protector and the POWER WHITENING XTRA 38% (MIXED 30% HP) Gel can also be applied lingual). The gel should be applied about 1 – 2 mm thick (avoid contact with other parts of the body and clothes). Allow the gel to work about 12 – 15 minutes. Reduce the application time and cycles if patient experiences sensitivity. Remove the gel from the teeth using suction – do not rinse – then apply a fresh layer of gel to the teeth. This process can be repeated three times. After the last treatment remove the gel first using suction and then by rinsing. Now the Gingiva Protector, lip holder and cotton wool balls can be removed. Finally the patient should rinse their mouth with warm water. The tooth colour after treatment should also be documented for your records. After the treatment apply a fluoride varnish to the teeth.

**Internal bleaching of endodontically treated teeth:**

Wait at least one week following the endodontic procedure before performing non-vital tooth whitening and ensure the tooth is endodontically sound.

1. Check the X-rays, to ensure that an adequate apical endodontic sealer is in place and to check where the alveolar bone margin runs.
2. Remove the filling material from the access hole until the depth that you have defined on basis of the X-ray image. Continue to remove the filling material and / or temporary endodontic materials in order to allow the tooth whitening agent contact with the internal structure of the tooth. Only this procedure allows best results.
3. Seal the root with conventional glass ionomer cement. Apply a 1-2 mm thick layer to ensure a secure seal.
4. The internal whitening is performed according to the instructions for vital teeth, but additionally intracoronally. Cover the gingiva around the tooth with Gingival Protector. Apply POWER WHITENING XTRA 38% (MIXED 30% HP) gel into the cavity and make sure about the tooth structure is covered with 1-2mm gel. Apply the gel also on the outside surface of the tooth. Wait 15-20 minutes. Suck of the gel. Procedure might be repeated for 1 or 2 times. After last application suck of gel and rinse remaining gel from cavity.
5. Use eugenol temporary restorative material to seal the opening.
6. Before starting work on adjacent teeth, it is recommended to wait at least two weeks after the whitening treatment - in order to stabilise the tooth colour

**Side effects:** If the patient feels a burning sensation in the mouth, treatment should be stopped at once and the gel should be removed. Sensitivity can be reduced by using a fluoride gel. For safety reasons, the patient should also wear safety goggles. For further information please refer to the general instructions for dentists.

**Storage and shelf life:** POWER WHITENING XTRA 38% (MIXED 30% HP) should be stored in a cool place or refrigerator (4 – 12 °C). The expiry date on the package or syringe refers only to products that have been stored correctly. Do not freeze the gel and protect it from heat and sunlight.

**After whitening treatment:** For at least 4 hours after whitening session, the patient should avoid consuming tea, coffee, red wine, cigarettes etc. as this will reduce effectiveness of treatment.

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WHITEmile GmbH Weinheimer Str. 6 · D- 69488 Birkenau Tel. +49(0)6201-8432190
**M+W Dental**® Vertrieb durch Müller & Weygandt GmbH · Reichardsweide 40