CE 0123

Dental Scaler Instruction Manual

(Model: PT X)

Please read this manual before operating the device

Guilin Refine Medical Instrument Co., Ltd.
SE1-PTX-M001 Version: 1.0 20260112

Safety Precautions

⚠ Warning: Neglecting these precautions may result in personal injury, such as electric shock, fire, or damage to the product.

- 1) This product is not a home appliance. It is intended solely for use in hospitals and dental clinics. The product must comply with the relevant operating specifications, as well as applicable laws and regulations in the medical field. The user must be a professionally trained and qualified dentist or technician. Adequate protection (e.g., goggles, mask) should be worn during operation to prevent cross-infection.
- 2) Please use an independent power outlet with protective grounding. Do not unplug the power cord with wet hands.
- 3) Please fully insert the power plug into the power socket. Do not use any power source other than the specified voltage.
- 4) Do not damage, modify, pull, bend, or twist the power cord excessively. Do not place heavy objects on the power cord.
- 5) After using the machine, turn off the power switch and remove the power plug.
- 6) When the equipment recovers after being powered off, please wait for the equipment to stabilize before stepping on the pedal. The operator should release the pedal promptly when the power is off.
- 7) Do not place the product on an unstable workbench, such as a swaying table, inclined surface, or a location subject to vibration.
- 8) This equipment is a reusable product. Keep it clean before and after use. Before each use, the handpieces, tips, torque wrench, etc., must be disinfected and sterilized. It is recommended to disinfect and sterilize according to the method described in Chapter 5 of this manual.
- 9) Do not knock or scrape the handpiece. Do not pull the cable forcefully during use to avoid damage.
- 10) Do not bend or polish the tips. The tips must be screwed into the handpiece with a torque wrench, and water spray must be generated during operation.
- 11) Before each operation, work outside the patient's mouth for more than 10 seconds with water running to drain the water in the pipeline at the back of the handpiece.
- 12) When the tip is damaged or worn out, the vibration intensity may decrease. The operator should replace the tip promptly according to the clinical situation. It is not recommended to use the tip if the abrasion exceeds 2mm. Use the tip comparison indicator to check. If the length is beyond the first line, the power remains the same; if it is between the first and second line, the power decreases to 80%; if it reaches or exceeds the second line, the power drops significantly, and tip replacement is recommended.
- 13) Do not screw or unscrew the tip while stepping on the foot switch or when the machine is operating.
- 14) The inner thread of tips manufactured by certain manufacturers may be rough, rusty, and could damage the thread system, resulting in irreparable harm to the scaler. Please use corresponding tips from the Refine brand.
- 15) Choose the appropriate power based on the type of tip. Refer to the Annex Tip Power Chart for the recommended power and water volume for different types of tips in the system. You can double-click the current mode to select the tip.
- 16) Do not use unclean water.
- 17) Improper cleaning and treatment of titanium implants, porcelain restorations, etc., can cause loosening of the adhesive, cracking of the porcelain restorations, or even damage to the porcelain. The cleaning or treatment of the oral cavity should be done with caution.
- 18) Noise will be generated during the scaling process. Those sensitive to noise can wear earplugs.
- 19) We are responsible for safety only when maintenance, repair, and modification of the machine are performed by the company or its authorized distributors, replacement parts are from the company, and operation follows the instructions in the manual.
- 20) This equipment may cause electromagnetic interference. Do not use it near electronic surgery equipment, and exercise caution in environments with strong electromagnetic interference.
- 21) This equipment does not contain toxic or hazardous substances and can be disposed of in accordance with relevant laws and regulations governing the disposal of medical devices.
- 22) Do not step on the foot pedal when the air polish handpiece cable is disconnected from the main unit.
- 23) Before using the ultrasonic handpiece, ensure the air polish handpiece cable is correctly placed on the handle bracket. Similarly, before using the air polish handpiece, ensure the ultrasonic handpiece cable is correctly placed on the ultrasonic handle bracket.
- 24) The tip must be securely tightened.
- 25) The air polishing handpiece nozzle should not be pointed at people.
- 26) Injury may occur if the powder accidentally sprays into the eyes. We strongly recommend that all personnel (doctors, nurses, patients) wear goggles during air polishing treatment.
- 27) Do not remove the powder tank while stepping on the pedal or when the machine is operating.
- 28) Before replacing the air polishing handpiece or nozzle, use a syringe to blow moisture from the joints at both ends (especially the gas interface) to prevent moisture from entering the gas path and avoid clogging of the powder in the pipeline.
- 29) The device must not be used in an MRI environment, as it is easily affected by electromagnetic emissions and may not function properly.
- 30) If the main unit expires, or if misuse leads to a short circuit of the circuit board or accidental dropping causes damage to components, the device should no longer be reused.

Contraindications

- 1) This equipment cannot be used on patients with hemophilia.
- 2) This equipment cannot be used on patients or doctors with a heart pacemaker.
- 3) This equipment should be used with caution on patients with heart disease, pregnant women, and children.
- 4) This equipment cannot be used on patients with respiratory diseases such as asthma and chronic bronchitis

Symbol instruction

Symbol	Implication	Symbol	Implication
	Trademark		Foot switch
	Caution		Humidity limitation: 10%-93%
	Recovery		Waste electrical and electronic equipment
	Date of manufacture		Clean mode
	ON (power)		Serial number
	Functional earth		Unique device identifier
	Atmospheric pressure limitation: 70kPa-106kPa		Protective earth (ground)
	Fragile, handle with care		Sterilizable in a steam sterilizer (autoclave) at 154°C
	Anhydrous mode		For indoor use only
	Authorized representative in the European Community		Alternating current
	Batch number		Power switch
	Consult instructions for use		Temperature limit: -20°C- +55°C
	Anti-drip equipment		Keep dry
	Manufacturer		Boost mode
	Type B applied part		CE marking with identification number of the Notified Body
	OFF (power)		Medical device

Forward

Guilin Refine Medical Instrument Co., Ltd. is a high-tech enterprise integrating R&D, production and sales of dental products, with a complete quality assurance system.

1. Product introduction

1.1. Brief introduction

The PT X dental scaler combines an ultrasonic scaler and an air polishing scaler. Its main features are:

- 1) Automatically switches the working mode according to the selected handpiece.
- 2) Touchscreen LCD display.
- 3) Special chemical solutions, such as hydrogen peroxide, sodium hypochlorite, and chlorhexidine, can be used to enhance the clinical treatment effect.
- 4) The LED handpiece enables convenient clinical operation.
- 5) The automatic frequency tracking system ensures optimal working conditions, resulting in more stable device performance.
- 6) The detachable ultrasonic handpiece and air polishing handpiece can be sterilized at a high temperature of 134°C and a pressure of 0.22 MPa.

1.2. Model: PT X

1.3. Intended use

This product includes an ultrasonic system and an air polish system. The ultrasonic system is used for periodontal treatment in oral clinical care to remove tartar and plaque, as well as for root canal cleaning. The air polish system is used to remove plaque and pigment and can also be used for implant maintenance.

1.4. Intended Patient Population

Adults and pediatric patients

1.5. Device Configuration

Please refer to the packing list for device configuration details.

1.6. Structure and Components

The machine consists of the main unit, handpieces (ultrasonic scaler handpiece, air polishing handpiece), cables, bottle, dental air polishing powder, torque wrench, foot pedal, and power cord.

1.7. Replacement instruction for main accessories

Name	Model	Replacement instruction
Air polishing handpiece	XP-1 / XP-2	Reusable and durable; replace the handpiece when it is blocked and cannot be cleared, when sealing performance is reduced, or when the ventilation pipe is significantly worn out.
ultrasonic scaler handpiece	HY-2L	Reusable and durable; replace the handpiece when the output power drops significantly or it no longer functions.
Power cord	/	Reusable and durable; replace the power cord if it is cracked or non-conductive.
Water pipe	/	Reusable and durable; replace the water pipe when it is cracked or leaking.
Torque wrench	/	Reusable and durable; replace it as per the instruction manual when it is broken.
Foot pedal	RFS05	Reusable and durable; replace it as per the instruction manual when it is broken.
Water bottle	/	Reusable and durable; replace it as per the instruction manual when it is broken.
Powder tank	/	Reusable and durable; replace it as per the instruction manual when it is broken.

1.8. Device safety classification

- 1) Classified by operation mode: Continuous operating device
- 2) Type of protection against electric shock: Class I
- 3) Degree of protection against electric shock: B type
- 4) Applied part of the equipment: tip
- 5) Degree of protection against harmful ingress of water: Ordinary equipment (IPX0). The foot pedal is an anti-drip device (IPX1).
- 6) Degree of safety application in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide: This equipment cannot be used in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide.

1.9. Main technical specification

Specification	PT X
Tip thread specifications for compatibility	M3*0.5
Main Unit Size (mm) Length*Width*Height	320mm*330mm*120mm
Main unit weight	4.0kg
Handpiece model	HY-2L/XP-1/XP-2
Control mode	Touch control
Water supply mode	Water tank
Power adapter input	100VAC - 240VAC 50Hz/60Hz
Input power	200VA
Main unit fuse	T1.6AH 250V
Foot pedal shell protection Level	IPX1
Tip output features	Frequency: 30kHz±5kHz; Output main vibration offset of tip 1µm - 90µm; Semi-off-set force 0.1-5N; Output power of tip 3W-20W
Water inlet pressure	0.1bar-5bar (0.01MPa-0.5MPa)
Air inlet pressure	5.5bar-7.5bar (0.55MPa-0.75MPa)
Water temperature	≤ 45°C

1.10. Operation environment

Environmental temperature: +5°C - +40°C
Atmospheric pressure: 70kPa - 106kPa

Relative humidity: 30% - 75%
Cooling water temperature: +5°C - +25°C

2. Installation

2.1. Main unit diagram


Figure 1 Front View


Figure 2 Back View

2.2. Touch panel


Figure 3 Touch panel diagram

2.3. Symbol Implication

Symbol	Implication	Symbol	Implication
Tip P1	Tip currently in use	⚙️	Setting
E	Endo mode	🔊	Turn up water vol-ume / power / air pressure
💾	Data storage function	P	Periodontal mode
🔊	Turn down water vol-ume/power/ air pressure	C	Cleaning mode
💧	water volume	🔌	Water heating switch
G	Scaling mode	📊	Power adjustment
A	Air polishing mode		

2.4. Schematic diagram of handpiece


Figure 4 Ultrasonic handpiece

Figure 5 Air polishing handpiece (supragingival)

Figure 6 Air polishing handpiece (subgingival)

2.5. Schematic diagram of tip installation


Figure 7
Schematic diagram of
tip installation


Figure 8
Schematic diagram
of external air supply
installation

2.6. Installation procedures

- 2.6.1. Open the package, check whether the equipment is complete according to the packing list, and place the main unit on a stable surface, facing the operator.
- 2.6.2. Plug the external air pipe connector into the air intake connector on the back of the main unit (as shown in Figure 8).
- 2.6.3. Fill the water bottle with an appropriate amount of water, and place the water bottle into the sink seat directly above the main unit. (It is recommended to apply a layer of petroleum jelly on the O-ring at the bottom of the water bottle to facilitate plugging and unplugging.)
- 2.6.4. Installation of the Foot Pedal
 - 1) Wired mode: Connect one end to the foot pedal port on the foot pedal, and the other end to the foot pedal port on the main unit.
 - 2) Wireless mode: Remove the pedal battery cover and install 2 AA batteries. Switch on the main unit, and the foot pedal will be connected and controlled. (When leaving the factory, the Bluetooth foot pedal and the main unit are pre-paired. If you replace the foot pedal, please enter "System Settings" to pair the new Bluetooth foot pedal.)
- 2.6.5. Connect the ultrasonic handpiece and the air polishing handpiece to the corresponding cables, and place the handpieces on the brackets on both sides of the main unit. The ultrasonic handpiece should be placed on the left, and the air polishing handpiece on the right.
- 2.6.6. Turn off the power switch, connect the power cord to the power socket on the back of the main unit, and then plug the power cord into the power supply, as shown in Figure 2.

⚠️ Warning! When the power cord is connected to the network power supply, it must be connected to a protective ground. Do not place or install the product in a location where it is difficult to disconnect the power supply.

3. Function and operation

3.1. Description of the Foot Pedal Function

Different machines are equipped with various types of foot pedals; please refer to the actual packing list.

3.1.1. Multifunctional Foot Pedal

- a) Choose the foot pedal mode according to section 2.6.4, and place it on a flat surface.
- b) The pedal and its functions are shown below:

Button	Operating mode	Function	
		Ultrasonic system	Air polishing system
A	Standard	Vibration + Water	AirPowder + water
B	No-water mode	Vibration	Air only
C(A+)	Boost mode	5-Level increase in Power	5-Level increase in air pressure
D	Cleaning mode	Water only	Water + Air

Note: Under Boost mode, Power or air pressure would be increased by 5 grades, the maximum is the 12th grade. After releasing button C, it will automatically return to the previously set grade.

Figure 9
Schematic diagram of multi-function foot pedal

- 3.1.2. Standard foot pedal (round)
See description of foot pedal key A in the table above.


Figure 10
Schematic diagram of standard foot pedal (round)

3.2. Ultrasonic system

3.2.1. Scaling

- 1) Turn on the power switch and pick up the ultrasonic handpiece. The panel will automatically switch to Ultrasonic mode.
- 2) Press the "G" button to enter Scaling mode.
- 3) Select the tip and tighten it onto the handpiece using a torque wrench.
- 4) When stepping on button A, the tip vibrates with the LED on and water emits. Release the button, and the vibration and water stop. The LED will turn off 10 seconds later.
- 5) The handpiece is generally held in a pen-holding position.
- 6) The frequency is extremely high when the machine is working. It can remove tartar as the tip softly strokes the tooth surface in an erasing motion. DO NOT apply excessive pressure or stay in one spot for too long.
- 7) Vibrating intensity: Adjust the vibrating intensity according to your needs. The default power setting is level 3. Please adjust the intensity based on tooth sensitivity and the hardness of calculus during clinical application.
- 8) Water volume: Adjust the water volume on the panel.
- 9) During clinical scaling, keep the side of the tip in contact with the tooth surface horizontally, applying light lateral pressure to allow the tip to vibrate freely.
- 10) After operation, keep the device running for 30 seconds with water to flush the handpiece and tips. Remove the tip and sterilize it.

3.2.2. Periodontal Treatment

- 1) Use a torque wrench to tighten the periodontal treatment tip to the ultrasonic handpiece. Click the "P" button on the panel to enter Periodontal Treatment mode.
- 2) The operation and adjustment methods are similar to Scaling mode.

3.2.3. Endodontic Irrigation

- 1) Attach the endo file to the handpiece using the endo wrench.
- 2) Click the "E" button to enter Endo mode.
- 3) The default power is level 1 in Endo mode. The operator can adjust the power as needed for clinical treatment.
- 4) Select the appropriate Endo file and gently place it into the root canal of the patient's tooth. Step on the foot pedal to begin ultrasonic endodontic irrigation.
- 5) During clinical cleansing, the file should not be pressed too firmly in the root canal.
- 6) Do not step on the pedal before the file is positioned in the root canal.
- 7) The recommended power range for endodontic irrigation is level 1 to level 5.

3.3. Air Polishing System

- 1) Add an appropriate amount of powder to the powder tank (the amount should be between the "Max" and "Min" markings on the tank), then tighten the upper cover of the powder tank and secure the tank on the main unit.
- 2) Pick up the air polishing handpiece, and the panel will automatically switch to Air Polishing mode.
- 3) Adjust the water volume and air pressure, aim the nozzle at a pool, press the foot button A, and verify that the nozzle sprays gas, powder, and water properly before use.
- 4) Before air polishing treatment, ensure the patient wears goggles and a veil. Users should also wear goggles or a protective mask.
- 5) The handpiece is generally held in a pen-holding position.
- 6) Adjust the water volume and air pressure to the appropriate levels. The recommended water volume starts at level 5, and air pressure starts at level 3. Depending on the sensitivity of the patient's teeth and the hardness of the dental plaque, adjust the water volume and air pressure during the clinical process. Increasing the air pressure enhances the cleaning effect but weakens the polishing effect, while increasing the water volume enhances the polishing effect but weakens the cleaning effect.
- 7) The nozzle should be aligned with the tooth surface during cleaning, but avoid direct contact. Keep the nozzle 3-5mm from the tooth surface at a 30-60 degree angle. The smaller the angle, the larger the cleaning area. During cleaning, perform a small circular motion on the tooth surface.
- 8) The air/powder mixture reflected from the tooth surface should be evacuated using a strong suction device on the dental unit during treatment.
- 9) After treatment, adjust the water volume to the maximum level and polish the tooth surface.

3.4. Cleaning Mode

- It is recommended to flush and disinfect the device's pipeline daily. "Cleaning" mode enables cleaning and disinfection of the pipeline to reduce crystal buildup and bacterial accumulation.
- 1) Fill the water tank with distilled or demineralized water.
 - 2) Pick up the ultrasonic handpiece, point it at the sink, click the "Cleaning" button on the screen, and step on the pedal to begin cleaning.
 - 3) After 30 seconds of cleaning, release the pedal to stop.
 - 4) After cleaning, place the ultrasonic handpiece back into its bracket. Then, pick up the air polishing handpiece, point the nozzle at the pool, and click the "Cleaning" button again. The device will automatically blow out residual powder from the pipeline and release high-pressure gas from the powder tank.
 - 5) After cleaning for 15 seconds, the device will automatically exit Cleaning mode. You can also click "Purge" on the screen to stop cleaning.

3.5. Settings

- 1) Click "ON"/"OFF" to turn the heating function on or off. (Water temperature will vary depending on the actual water output.)
- 2) Click the "Settings" button at the bottom right corner of the screen to access the "System Settings" interface.
- 3) Click the "Calibration" button, and the system will automatically calibrate the air pressure to the value set in the program.
- 4) The Bluetooth pedal is pre-paired with the machine. Any subsequent replacement of the foot pedal or main unit requires pairing again.

4. Troubleshooting

4.1. Troubleshooting

Fault	Possible cause	Solutions
The tip does not vibrate, and no water flows out after powering on and stepping on the foot switch.	Poor power plug connection	Ensure the power plug is securely connected.
	Poor foot pedal switch connection	Ensure the foot pedal plug is securely connected.
	Blown fuse	Replace the T1.6AL 250V fuse.
	Battery depleted in wireless pedal mode	Replace the battery.
	Poor connection in wired pedal mode	Securely connect and fasten the wire.
The tip does not vibrate, but water sprays after powering on and stepping on the foot pedal.	Main unit malfunction	Contact your local distributor or the manufacturer.
	The tip is loose.	Securely tighten the tip.
	Loose connection between the tail wire and the circuit board.	Contact your local distributor or the manufacturer.
	Handpiece malfunction.	Contact your local distributor or the manufacturer.
The tip vibrates, but no water sprays after powering on and stepping on the foot pedal.	Cable malfunction.	Contact your local distributor or the manufacturer.
	Increased resistance in the diaphragm pump causing insufficient torque for water retention.	Increase the water volume.
	Contaminants in the solenoid valve.	Contact your local distributor or the manufacturer.
Water continues to spray after power is turned off.	Clogged water line.	Contact your local distributor or the manufacturer.
	Contaminants in the solenoid valve.	Contact your local distributor or the manufacturer.

Handpiece overheating.	Insufficient water flow.	Increase the water volume.
	Device malfunction.	Contact your local distributor or the manufacturer.
Water spray is insufficient.	Insufficient water volume.	Increase the water volume.
	Clogged water line.	Contact your local distributor or the manufacturer.
Weak tip vibration.	The tip is not securely tightened.	Securely tighten the tip.
	The tip has loosened.	Securely tighten the tip.
The endo file does not vibrate.	Damaged tip (See Note 1).	Replace the tip.
	Thermocouple malfunction.	Contact your local distributor or the manufacturer.
No air or water spray after powering on and stepping on the foot pedal.	The nut is loose.	Securely tighten the nut.
	Loose power supply connection.	Ensure the power supply plug is securely connected.
	Poor connection in wired pedal mode.	Securely connect and fasten the wire.
The nozzle does not spray air, but water sprays after powering on and stepping on the foot pedal.	Depleted battery in wireless pedal mode.	Replace the battery.
	Clogged nozzle.	Clear the nozzle blockage.
	Clogged handpiece.	Clear the handpiece blockage.
Air flows but no water sprays after powering on and stepping on the foot pedal.	Clogged handpiece tail cord.	Remove the tail cord from the main unit, clear the blockage, or replace it.
	Solenoid valve malfunction.	Contact your local distributor or the manufacturer.
	Increased resistance in the diaphragm pump causing insufficient torque for water retention.	Increase the water volume
Air leakage in the powder tank.	Contaminants in the solenoid valve.	Contact your local distributor or the manufacturer.
	Clogged water line.	Contact your local distributor or the manufacturer.
	Damaged O-ring at the base of the powder tank.	Unscrew the powder tank base, remove the damaged O-ring, replace it with an O-ring of the same specification, and securely fasten the powder tank base.
	Defective or worn O-ring.	Install a new O-ring.
	Powder residue in the threading prevents proper sealing.	Clean any residual powder from the threaded area.
Handpiece water leakage	The upper cover of the powder tank is damaged.	Replace the upper cover of the powder tank.
	The powder tank thread is damaged, preventing the screw from securing properly.	Replace the powder tank.
	Damaged O-ring in the handpiece	Install a new O-ring
Reduced air powder scaling efficiency	Insufficient powder in the tank	Refill the powder tank.
	Powder residue in the pipe, handpiece, or nozzle passage	Clear the passage with a fine needle and use compressed air to remove residue.
Touch panel displays a prompt message	Powder tank is missing	Verify the powder tank and reinstall if necessary.
	Air pressure is outside the acceptable range	Adjust the external air pressure to the required range: 5.5 bar - 7.5 bar (0.55 MPa - 0.75 MPa).
	Two handpieces were picked up simultaneously	Ensure only one handpiece is picked up before stepping on the pedal.
	No handpiece has been detected	Ensure only one handpiece is picked up before stepping on the pedal.
	The selected handpiece does not match the configured mode	Select the appropriate handpiece and corresponding mode.

Note: If the problem still can't be solved, please contact your local dealer or the manufacturer.

4.2. Notes

Note 1: If the tip is securely tightened and a water mist is still sprayed, the tip is considered damaged if the following phenomena occur:

- 1) The vibration intensity of the tip and the degree of water atomization are significantly reduced.
- 2) The tip produces a harsh clicking sound during operation.

5. Cleaning, Disinfecting, and Sterilizing

- 5.1. Follow the recommendations outlined in the manual "Reprocessing Instructions for Cleaning, Disinfecting, and Sterilizing" provided with your product for the proper procedures of cleaning, disinfecting, sterilizing, and packing the components. Adhere to the current reprocessing regulations enforced in your country. Avoid using scouring powder or abrasive sponges, as they may damage the surface.

- 5.2. The handpiece, endo wrench, torque wrench, LED lamp, and light conductor (for handpieces with an LED lamp) can be sterilized.

- Recommended sterilization cycles for the handpiece: 300 cycles
- Recommended sterilization cycles for the torque wrench and endo wrench: 300 cycles

Notice:

- a) For cleaning and disinfection, use ethanol or isopropanol disinfectant without additives (such as quaternary ammonium salts). Using other disinfectants may cause discoloration or cracking.

* For more details, contact the disinfectant manufacturer.

* In this manual, ethanol or isopropanol disinfectant without additives (such as quaternary ammonium salts) is referred to as "disinfectant alcohol."

- b) Clean the handpiece with compressed air before sterilization.
- c) Ensure that the scaling tip is unscrewed from the handpiece, as it must not be sterilized with other components.
- d) Check for any damage to the exterior of the handpiece during treatment and sterilization. Do not apply any protective oil to the surface of the handpiece.
- e) The handpiece has two waterproof "O" rings at its end. Lubricate them frequently with dental lube, as sterilization and repeated insertion and removal will shorten their lifespan. Replace them if they become damaged or excessively worn.
- f) The following sterilization methods are prohibited:
 - Boiling in water.
 - Immersion in iodine, alcohol, or glutaraldehyde.
 - Baking in an oven or microwave.

⚠ Notice: We are not responsible for any damage caused by the above methods.

6. Transportation, storage and maintenance

6.1. Transportation

- 6.1.1. Prevent excessive impact and shaking during transportation. Handle it carefully and do not invert it.
- 6.1.2. Do not transport it with dangerous goods.
- 6.1.3. Avoid exposure to direct sunlight and protect it from rain or snow during transportation.

6.2. Storage

- 6.2.1. The device should be handled carefully and lightly. Ensure that it is kept away from vibration and installed or stored in a cool, dry, and well-ventilated place.
- 6.2.2. Do not store the machine with articles that are poisonous, combustible, caustic, or explosive.
- 6.2.3. This machine should be stored in a room where the relative humidity is 10% to 93%, atmospheric pressure is 70kPa to 106kPa, and the temperature is between -20°C and +55°C.
- 6.2.4. When the device is not in use, turn off the power supply and unplug the power plug. If it is not used for an extended period, it should be powered on and connected to water and air for five minutes once a month.

6.3. Maintenance

- 6.3.1. Air Filter
 - 1) When water accumulates in the filter, turn the knob at the bottom of the filter counterclockwise to drain the water, then tighten the knob clockwise.
 - 2) Replacement of the filter element: Use a filter wrench to unscrew the transparent cover of the air filter, then unscrew the black nut at the lower end of the filter element. Remove the white filter element and discard it in the trash, replace it with a new filter element, and reinstall the black nut and transparent shell. It is recommended to replace the filter element every 24 months, and a spare filter element is included in the accessories.


Figure 11 Schematic diagram of filter element replacement

7. Environmental Protection, Disposal, and Scrapping

The instrument contains no hazardous components and can be disposed of and destroyed in accordance with the relevant local regulations.

Part	Toxic or harmful substances or elements					
	Pb	Hg	Cd	Cr6+	PBB	PBDE
Main unit	○	○	○	○	○	○
handpiece	○	○	○	○	○	○
Tip	○	○	○	○	○	○
Nozzle Foot pedal	○	○	○	○	○	○
Mechanical elements, including bolts, nuts, washers, etc.	○	○	○	○	○	○
○ : Indicates that the content of the toxic substance in all homogeneous materials of the part is below the limit requirement stipulated in SJ/T-11363-2006 Limit Requirements for Toxic and Hazardous Substances in Electronic Information Products. X : Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceed the limit requirement specified in SJ/T-11363-2006. Please dispose according to the local laws or consult with dealer from whom you purchased it about waste disposal.						

8. After-Sales Service

- 8.1. If the machine fails to function properly due to quality issues after purchase, our company will be responsible for repairs with a valid warranty card. For specific details, please refer to the warranty instructions provided in the warranty card.
- 8.2. This machine does not contain user-serviceable parts. All maintenance, adjustments, calibration, and modifications of technical parameters must be performed exclusively by the company's technicians or authorized repair centers. If the customer chooses to perform repairs independently, our company may, upon request, provide circuit diagrams, component lists, legends, calibration guidelines, or other necessary information to assist qualified technicians in repairing designated parts. However, our company will not be held responsible for any resulting consequences.
- 8.3. Users must use only original accessories. Please contact your local dealer or our company to purchase them. The use of accessories from other brands is strictly prohibited to prevent equipment damage or other hazards.
- 8.4. Users should not attempt to repair damaged handpieces, scaling tips, or other accessories. Instead, they must purchase new parts and replace them before use. For further information, please contact the manufacturer.

9. EMC - Declaration of conformity

9.1. Instructions for use:

The ME EQUIPMENT or ME SYSTEM is suitable for use in hospitals or dental clinics.

- ⚠ Warning: Do not operate near active HF surgical equipment or the RF-shielded room of an ME system for magnetic resonance imaging (MRI), where electromagnetic (EM) disturbances are intense.
- ⚠ Warning: Avoid placing this equipment adjacent to or stacked with other devices, as it may cause improper operation. If such placement is necessary, both this equipment and the other devices should be monitored to ensure normal operation.
- ⚠ Warning: Using accessories, transducers, or cables other than those specified or provided by the manufacturer may increase electromagnetic emissions or reduce electromagnetic immunity, leading to improper operation.
- ⚠ Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should not be used closer than 30 cm (12 inches) to any part of this equipment, including manufacturer-specified cables. Otherwise, this may lead to performance degradation.
- ⚠ Note: The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If used in a residential environment (where CISPR 11 Class B is typically required), this equipment may not provide adequate protection for radio-frequency communication services. Users may need to implement mitigation measures, such as relocating or reorienting the equipment.

List of all cables:

Number	Name	Length	Shielded or not	Detachable or not	Note
1	Power cord	2.0m	No	Yes	/
2	Foot pedal cord	2.8m	No	Yes	/
3	Ultrasonic handpiece cord	2.0m	No	No	/
4	Air polishing handpiece cord	2.0m	No	No	/

Replaceable accessories please refer to 1.7 replacement instruction for main accessories.

Performance of the ME Equipment:

- The Dental Scaler PT X achieves ultrasonic cleaning through vibration of the ultrasonic handpiece. Depending on the application, the "G", "P", and "E" modes can be selected. The "A" mode adjusts the air polishing function by controlling the air pressure to modify the polishing power; "C" mode is the cleaning mode for either the ultrasonic or air polishing handpiece.
- After the operation, the ultrasonic handpiece pipeline and air polishing handpiece pipeline must be cleaned.
- When the equipment's essential performance is lost or degraded due to electromagnetic disturbances, the doctor should immediately stop using it to avoid treatment errors.
- Remove the source of disturbances or adjust the direction or position of the ME equipment to ensure it can be used in normal operating condition.

9.2. Technical description

9.2.1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the excepted service life:

- 1) Portable and mobile RF communications equipment may affect the performance of this device. Avoid using the equipment in areas with strong electromagnetic interference, and keep it away from mobile phones, microwave ovens, and similar devices.
- 2) Avoid placing this equipment adjacent to or stacked with other devices, as it may cause improper operation. If such placement is necessary, both this equipment and the other devices should be monitored to ensure they are operating normally.
- 3) Except for cables sold by the manufacturer as spare parts for internal components, the use of accessories and cables other than those specified or provided by the manufacturer may increase electromagnetic emissions or decrease electromagnetic immunity, leading to improper operation.
- 4) Using accessories, transducers, and cables other than those specified or provided by the manufacturer in conjunction with this equipment may increase electromagnetic emissions or decrease electromagnetic immunity, leading to improper operation.

9.2.2. Guidance and Manufacturer's Declaration - Electromagnetic Emissions and Immunity

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class A
Harmonic emissions IEC 61000-3-2	Not Applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable

Table 2

Guidance and manufacturer's declaration - electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/ burst IEC 61000-4-4	±2 kV power supply lines ±1 kV signal input/output 100 kHz repetition frequency	±2 kV power supply lines Not applicable 100 kHz repetition frequency
Surge IEC 61000-4-5	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode	±0.5 kV, ±1 kV differential mode Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle	0 % UT; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF IEC61000-4-3	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
Note: UT is the a.c. mains voltage prior to application of the test level.		

Table 3

Guidance and Declaration - Electromagnetic Immunity						
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
	385	380–390	TETRA 400	Pulse modulation 18 Hz	27	27
	450	430–470	GMRS 460,FRS 460	FM ± 5 kHz deviation 1 kHz sine	28	28
	710	704–787	LTE Band 13,17	Pulse modulation 217 Hz	9	9
	745					
	780					
	810					
	870	800–960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28	28
	930					
	1720					
	1845	1700–1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	28	28
	1970					
	2450	2400–2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28	28
	5240	5100–5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9	9
	5500					
	5785					

Table 4

Guidance and manufacturer's declaration - electromagnetic Immunity				
Radiated RF IEC61000-4-39 (Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields)	Test Frequency	Modulation	IEC 60601-1-2 Test Level (A/m)	Compliance level (A/m)
	30 kHz	CW	8	8
	134,2 kHz	Pulse modulation 2.1 kHz	65	65
	13,56 MHz	Pulse modulation 50 kHz	7,5	7,5

Attachment: Reprocessing Instructions for Cleaning, Disinfecting, and Sterilizing

1. Getting Started

- 1.1. Read these operating instructions carefully, as they explain essential details and procedures. Pay special attention to the safety precautions. Always keep this manual readily accessible.
- 1.2. To prevent injury and property damage, follow the specified directives.
- 1.3. The instructions in this manual apply only to the product with which it was delivered.

2. Introduction

- 2.1. These reprocessing instructions provide guidelines for cleaning, disinfection, sterilization, and packaging of the manufacturer's reusable products intended for reprocessing in medical facilities.
- 2.2. The goal of reprocessing reusable products is to reduce bioburden and achieve sterility, thereby eliminating the risk of infection related to product reuse. Decisions regarding cleaning, disinfecting, or sterilizing the manufacturer's medical and dental instruments are based on the potential risk of infection associated with their use.
- 2.3. Steam sterilization is recommended.
- 2.4. Sterilization or high-level disinfection cannot be achieved unless all components of the assembly are thoroughly cleaned first.
- 2.5. If the manufacturer's reprocessing instructions appear to be inadequate, please inform the manufacturer.
- 2.6. We encourage you to report any adverse events related to device reprocessing directly to the manufacturer.

3. Reprocessing Instructions for Reusable Products

- 3.1. These instructions apply to the reprocessing of all reusable products (hereafter referred to as "products") from the manufacturer. Additional product-specific instructions may be included to provide further details as needed.
 - ⚠ **Important: Before use, carefully read the operating instructions for the manufacturer's instruments and devices with which the product will be used.**
- 3.2. Reusable products must be cleaned, disinfected, and sterilized before first use. Reprocessing procedures have limited impact on this device. The number of reprocessing cycles is primarily determined by the device's function and wear. There is no maximum number of allowable reprocessing cycles from a processing perspective. The device should no longer be used if there are signs of material degradation.
 - ⚠ **If the product is damaged, it must be reprocessed before being sent back to the manufacturer for repair.**

4. Preparation – Basic Principles

- 4.1. Effective sterilization can only be achieved after thorough cleaning and disinfection. As part of your responsibility for maintaining product sterility, ensure that only validated equipment and product-specific procedures are used for cleaning, disinfection, and sterilization. Additionally, validated parameters must be followed during every cycle.
- 4.2. Follow the applicable legal requirements in your country, as well as the hygiene regulations of the hospital or clinic. Special attention should be given to additional requirements for prion inactivation.

5. Preparation at the Point of Use

Disconnect the product. Immediately after use, remove heavy soiling from the products using cold water (<40°C). Do not use a fixating detergent or hot water (>40°C), as this can cause the fixation of residues, potentially affecting the outcome of the reprocessing process. Store the products in a humid environment.

6. Transportation

Ensure safe storage and transportation to the reprocessing area to prevent damage and contamination.

7. Preparation for Decontamination

Whenever possible, the products must be reprocessed in a disassembled state.

8. Pre-cleaning

Perform manual pre-cleaning until the products are visibly clean. Submerge the products in a cleaning solution and flush the lumens using a water jet pistol with cold tap water for at least 10 seconds. Clean all surfaces with a soft-bristle brush.

9. Cleaning

In terms of cleaning, disinfection, rinsing, and drying, a distinction should be made between manual and automated reprocessing methods. Automated reprocessing methods are preferable due to their better standardization potential and increased industrial safety.

Automated Cleaning:

Use a washer-disinfector (WD) meeting the requirements of the ISO 15883 series.

Put the products into the machine on a tray. Connect the products with the WD by using suitable adapter and start the program:

4 min pre-washing with cold water (<40°C)

Emptying

5 min washing with a mild alkaline cleaner at 55°C

Emptying

3 min neutralising with warm water (>40°C);

Emptying

5 min intermediate rinsing with warm water (>40°C)

Emptying

The automated cleaning processes have been validated by using 0.5% neodisher MediClean forte (Dr. Weigert).

- ⚠ **Acc. to EN ISO 17664 no manual reprocessing methods are required for these devices. If a manual reprocessing method has to be used, please validate it prior to use.**

- ⚠ **Purified water shall be used for cleaning.**

10. Disinfection

Automated Thermal Disinfection in washer/disinfector under consideration of national requirements in regards to A0 value (see EN 15883).

A disinfection cycle of 5 min disinfection at 93°C has been validated for the product to achieve an A0 value of 3000.

Purified water shall be used for disinfection.

11. Drying

Automated Drying:

Dry the exterior of the products at 40°C for 5 minutes using the drying cycle of a washer/disinfector. If necessary, additional manual drying can be performed using a lint-free towel. Insufflate the cavities of the products with sterile compressed air.

12. Functional Testing and Maintenance

Perform a visual inspection to ensure the products are clean, and reassemble them if required. Conduct functional testing according to the instructions for use. If necessary, repeat the reprocessing process until the products are visibly clean.

Before packaging and autoclaving, ensure that the products have been maintained according to the manufacturer's instructions.

13. Packaging

Package the products using appropriate sterilization packaging materials. The packaging material and system must comply with EN ISO 11607.

14. Sterilization

Sterilize the products using a fractionated pre-vacuum steam sterilization process (in accordance with EN 285, EN 13060, and EN ISO 17665), taking into account the respective country-specific requirements. Minimum requirements: 3 min at 134 °C (in EU: 5 min at 134 °C)

Maximum sterilization temperature: 138°C

Drying time:

For steam sterilization, we recommend a drying time of 15 to 40 minutes. Choose a suitable drying time, depending on the autoclave and load. Refer to the autoclave's instructions for use.

After sterilization:

- a) remove the product from the autoclave.
- b) Let the product cool down at room temperature for at least 30 minutes.
Do not use additional cooling.

Check that the sterilization wraps or pouches are not damaged.

- ⚠ **Flash sterilization is not allowed on lumen products.**

- ⚠ **The manufacturer assumes no responsibility for the use of other sterilization procedures (e.g. ethylene oxide, formaldehyde and low temperature plasma sterilization). In such cases, please observe the respective valid standards (EN ISO 14937/ANSI AAMI ISO 14937 or the procedure-specific standard) and verify the suitability and effectiveness in principle of the procedure (if necessary, including investigations on sterilizing agent residue), taking into account the specific product geometry as part of the validation.**

•Maximum sterilization temperature 138°C

15. Storage

Store sterilized products in a dry, clean, and dust-free environment with a relative humidity of 10% to 93%, an atmospheric pressure of 70 kPa to 106 kPa, and a temperature range of -20 °C to +55 °C. Refer to the label and instructions for use.

After sterilization, the product should be packaged in a medical sterilization bag or a clean, sealed container and stored in a designated storage cabinet. The storage period should not exceed 7 days. If exceeded, the product must be reprocessed before use.

16. Service Life

The products are designed for multiple sterilization cycles, with materials specifically selected for durability. However, repeated sterilization exposes the devices to thermal and chemical stresses, leading to material aging. If there are restrictions on the number of permissible re-sterilization cycles, they will be specified in the product-specific instructions.

- ⚠ **The use of ultrasound baths and strong cleaning or disinfection fluids (alkaline pH >9 or acidic pH <5) may shorten the device's lifespan. The manufacturer accepts no liability in such cases.**
- ⚠ **Devices must not be exposed to temperatures above 138 °C.**

The user is responsible for ensuring that reprocessing procedures, including resources, materials, and personnel, meet the required standards. Best practices and, in many cases, national regulations require these processes and resources to be properly validated and maintained.

Disposal and Scrapping

- The product does not contain a battery or toxic substances. Additionally, there are no components that require special removal from the main unit for disposal and scrapping.
- When the device reaches the end of its service life, do not dispose of it with household waste. Comply with the Waste Electrical and Electronic Equipment (WEEE) directives and the medical waste disposal regulations in your country.
- The torque wrench, endo wrench, and handpiece, which may come into contact with biological material and pose biological hazards, must be detached from the main unit and reprocessed according to the reprocessing guidelines in Chapter 5 before disposal and scrapping.

Statement

Note: Refine reserves the right to modify the product technology, accessories, instruction manual and product packaging content at any time without notice. The product is subject to the actual product, and the pictures are for reference only. The final interpretation right belongs to Guilin Refine Medical Instrument Co., Ltd.

Refer to the product label for the date of manufacture.

MedNet EC-REP C llb GmbH, Borkstrasse 10, 48163 Münster, Germany

Guilin Refine Medical Instrument Co., Ltd.
No.8-3, Information Industrial Park, High-Tech Zone, Qixing District,
541004 Guilin, Guangxi, PEOPLE'S REPUBLIC OF CHINA
Tel: +86-773-7796686 | Email: refine@refine-med.com
Website: http://www.refine-med.com

Scan and log in to the website for more information.

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MedNet EC-REP C.Ilb GmbH, Borkstrasse 10, 48163 Münster, Germany