solventum

Impregum™ F

(en) **Polyether Impression Material**
(de) **Polyether Abformmasse**
(fr) **Matériau pour empreintes à base de polyéther**
(it) **Materiale per impronte a base di poliestere**
(nl) **Materiaal de impressie a base de polyëter**
(es) **Material de polióter para impresiones**
(pt) **Polyether afdrukmateriaal**
(sv) **Avtrycksmaterial av polyeter**
(da) **Polyetherbaseret aftrykskompositmateriale**
(no) **Polyeter avtrykksmateriale**
(fi) **Polyeetteri-jäljennösmassa**
(el) **Πολυαιθέρας υλικό από πολυαιθέρα**

Instructions for Use
Gebrauchsinformation
Notice d'utilisation
Informazioni per l'uso
Instrucciones de empleo
Instruções de uso

Gebruiksinformatie
Πληροφορίες χρήσεως
Bruksanvisning
Käyttöohjeet
Brugsanvisning
Bruksveiledning

Solventum GmbH
Edisonstrasse 6
82194 Gröbenzell, Germany
Solventum.com

© Solventum 2025. Solventum, the S logo, Garant and Impregum are trademarks of Solventum or its affiliates.

D3283210
2025-01

(en) ENGLISH

Intended Use
Dental impression material.

Product Description
Solventum™ Impregum™ F is a medium-bodied consistency polyether impression material for manual mixing. The mixing ratio is 7 volumes base paste : 1 volume catalyst paste.

Impregum F is a multiple patient multiple use device.

Solventum Syringe Tips are single-use devices.

These Instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly readable. For details on all additional components/products please refer to the corresponding Instructions for Use. The Instructions for Use are available at: eIFU.Solventum.com.

Indications for Use
- Impression of inlay, onlay, crown, bridge, and veneer preparations
- Functional impressions
- Fixation impressions
- Implant impressions

Intended Patient Population
All patients requiring dental treatment, unless the patient's condition limits the use.

Contraindications
Polyether impression materials are contraindicated for use with patients who have a known allergy to polyether impression materials or contained allergens, as listed in section "Precautions and Warnings".

Intended Users
Educated dental professionals who have theoretical and practical knowledge and use of dental products.

Clinical Benefit
Create a precise reproduction of the clinical situation within dental workflows.

Precautions and Warnings

For Patients
The base paste and catalyst paste of this product contain substances that may cause an allergic reaction by skin contact in certain individuals. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If an adverse reaction is observed in a patient discontinue the use of the product in this patient and seek medical attention as needed.
Known allergens are 1-Dodecylimidazole, Carthamus Tinctorius (Safflower flowers) and Sulfonium salt.

For Dental Personnel
The base paste and catalyst paste of this product contain substances that may cause an allergic reaction by skin contact in individuals. In order to reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a face protection shield is recommended.

The safety and performance of the product has not been clinically evaluated in pregnant or lactating women, or children under 18 years old.

Potential Side-Effects
Potential side effects are local irritations, or local or systemic allergic reactions.

Residual Risks
In rare cases impression material can be swallowed or aspired. Swallowing of impression material can cause harm to the gastrointestinal tract. As the impression material cures in the mouth of the patient, it is possible that loose teeth or restorations are removed with the impression material.
In situations with improper in-office infection control or if the IFU was not followed, the following risks may occur: Injury, allergy, cytotoxic reaction or distorted impression which could lead to unspecified patient effect.

Instructions for Use
- Block-out undercuts or areas where gingival recession exists to prevent the impression material from "locking" into the tooth structure. Failure to do so out may make tray removal difficult, or cause extraction of natural teeth or prosthesis. This block-out is especially important when using impression materials with high Shore hardness.

Impression Tray
Any impression trays generally used for precision impressions are suitable.
- For sufficient adhesion, apply a thin layer of Solventum™ Polyether Adhesive/Solventum™ Polyether Contact Tray Adhesive to the tray and allow to dry completely according to the Instructions for Use.

Retraction
Solutions based on aluminium hydroxide chloride or aluminum sulfate are suitable retraction materials. Use non-impregnated threads or rings containing epinephrine (adrenaline), 8-hydroxyquinoline sulfate, and iron (III) sulfate for retraction may impair the setting of polyether impression materials.
- Keep the area free of blood and contaminants during biting.
- In subgingival preparations, hemostatic threads or rings may be used.
- Prior to taking the impression, completely remove any residual retraction agent by rinsing and drying.

Dosing and Mixing
- Dispense equal lengths of base paste and catalyst to a mixing block. Overand underdosing of the catalyst does not affect processing time, but impairs the quality of the impression.
- Mix pastes for 45 sec with a spatula. Touch all of the edge of the block unmixed quantities of paste adhering to the spatula, and repeatedly take up and streak out until the colour of the mixture is uniform. Do not mix by stirring in order to avoid the formation of bubbles!

Times

Mixing	Processing time from start of mixing	Setting time from start of mixing	Residence time in the mouth
min:sec	min:sec	min:sec	min:sec
00:45	02:45	06:00	03:15

These are the processing times at 23° C. Higher temperatures shorten, and lower temperatures prolong total processing times.

Impression Taking
Mono-phase technique
- By means of the filling device, place the mixture of pastes in the Elastomer Syringe.
- Then load the tray previously prepared with adhesive.
- Apply the mixture into the sulcus or cavity. Shortly before up, always keeping the tip of the application nozzle immersed in the material and in contact with the surface of the tooth.
- Immediately after application, position the loaded tray in the mouth, and hold steady until setting is complete without applying pressure.
- To remove the impression, first break the seal between the gingiva and the tooth surrounding dentition. Remove any residual cured impression material from the mouth.

Model Preparation
- Prepare a cast from the impression with a commercial specialised stone plaster no earlier than 30 min and no later than 14 days after taking the impression.
- To avoid introducing bubbles in the model, briefly pre-rinse with water and dry with air before pouring the stone slurry. Make sure to rinse impair the quality of polyether impressions and are not necessary!
- Polyether impressions can be silver-coated, whereas copper-coating is not feasible.

Storage Conditions/Shelf Life
Store the product at 5 - 25 °C/59 - 77° F. Do not refrigerate!
Do not use after the expiration date.
Shelf life (unopened) at 30° C/86° F in the dark.

Handling Before and/or After Use
- Elastomer Syringe Tips are single-use devices and would be damaged when removing and reinserting them.
- Please keep the single-use accessory in the original bag to guarantee the product's sterility.
- **Use the product at 18–25 °C/64–77 °F.**
- At temperatures below 18 °C/64 °F the viscosity of the pastes may increase sufficiently to make mixing in the device difficult. After keeping the pastes at 18–25 °C/64–77 °F for about 24 hours, the processability is re-established without compromising quality.
- Direct exposure to sunlight and damp storage conditions may damage the impression.
- Polyether impressions should not be exposed to solvent-containing liquids under any circumstances. Swelling and imprecise modeling may result.
- Solventum polyether impressions are not suitable for electrolytic copy polyethers. The working and setting times of the individual products must be strictly observed.

Cleaning and Disinfection

Multiple Use and Cleaning of Medical Device
If the tube is not visibly dirty and has not been used, contrary to expectations, there is contamination with blood or saliva during the treatment, for example, discard the tube immediately.

Specifications for Product Use
Before use of the tube, put on new, unused disposable gloves during all clinical applications of the tube to avoid the risk of contamination.
The material should be applied in advance to prevent contamination during the treatment.

Cleaning
The tube must not be placed into a disinfection bath or a cleaning-disinfection device (washer-disinfector). A superficial cleaning of the tube can be carried out as part of hygiene management, i.e., in the event of soiling from the environment, e.g., dirt, tooling does not mean contamination with blood, tissue, secretion, etc.). CaviWipes™ have been tested for use as a cleaning agent. Always observe all applicable legal and hygiene regulations for dental offices and/or hospitals. Cleaning is done using ready-to-use disinfectant wipes (e.g., CaviWipes™ for surfaces in accordance with the active substance compounds; alternatively, use equivalent materials). Follow the cleaning agent manufacturer's directions for use. Clean the medical device by wiping it for one min with a new cloth.

Inspection, Maintenance, Testing
Examine the medical device for damage, discoloration, and contamination before every use. Do not use damaged and contaminated devices under any circumstances.

Cleaning the Elastomer Syringe
Please refer to the pertinent Instructions for Use for information about cleaning and disinfecting the Elastomer Syringe.

Cleaning
Rinse parts that not set may be removed with ethanol or by rinsing with water and soap.
- The adhesive can be removed from metal trays with acetone.

SDS and Disposal
See the Solventum Safety Data Sheet (SDS)/Safety Information Sheet (SIS), available at www.Solventum.com or contact your local subsidiary, for disposal and ingredient information.
Dispose of the contents or content in accordance with the applicable regulations. Please pay special attention to the disposal of contaminated waste, waste with sharps to avoid health risks because of improper handling.

Customer Information
No person is authorized to provide any information which deviates from the information presented in this Instructions for Use.
Please report an undesirable outcome (e.g. adverse event, serious incident (SI)) occuring in relation to the device to Solventum and the local competent authority (EU) or local regulatory authority.

Symbol Glossary

Reference Number and Symbol Text	Symbol	Description of Symbol
ISO 15223-1 5.7.10 Unique device identifier	UDI	Indicates a carrier that contains unique device identifier information.
Non-corrugated fiberboard	21 PAP	Indicates product packaging is made of non-corrugated fiberboard. Official Journal of the EC, Commission Decision (97/129/EC)
Aluminum	41 ALU	Indicates product packaging is made of aluminum. Official Journal of the EC, Commission Decision (97/129/EC)

For more information, see eIFU.Solventum.com
Information valid as of December 2025

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