



Trade name: BioBleach soft

Substance number: S0018SOFT

Version: 1 / GB

Date revised: 06.02.2024

Replaces Version: - / GB

Print date: 06.02.2024

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

BioBleach soft

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/preparation

Dental splinting material short-term

1.3. Details of the supplier of the safety data sheet

Address/Manufacturer

Dreve Dentamid GmbH

Max-Planck-Straße 31

59423 Unna

Telephone no.

+49 2303 8807-0

Fax no.

+49 2303 8807-29

Information provided

Department Research & Development: Fax: +49 2303 8807-562

by / telephone

E-mail address of

sicherheitsdatenblatt@dreve.com

person responsible

for this SDS

1.4. Emergency telephone number

Henkel Fire Department / 24h-Emergency-Contact-No.: +49 211 797-3350

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

This product is not classified hazardous in accordance with Regulation (EC) No 1272/2008.

2.2. Label elements

Labelling according to regulation (EC) No 1272/2008

Supplemental information

EUH210

Safety data sheet available on request.

2.3. Other hazards

No special hazards have to be mentioned.

The Substance does not meet PBT-criteria. This substance does not meet the vPvB-criteria. This substance does not have endocrine disrupting properties with respect to humans. This substance does not have endocrine disrupting properties with respect to non-target organisms.

SECTION 3: Composition/information on ingredients

3.2. Mixtures



Trade name: BioBleach soft

Substance number: S0018SOFT

Version: 1 / GB

Date revised: 06.02.2024

Replaces Version: - / GB

Print date: 06.02.2024

Chemical characterization

Ethylene vinyl acetate copolymer (EVA)

Hazardous ingredients

vinyl acetate

CAS No. 108-05-4

EINECS no. 203-545-4

Concentration \geq 0,1 < 1 %

Classification (Regulation (EC) No. 1272/2008)

Flam. Liq. 2 H225

Acute Tox. 4 H332

Carc. 2 H351

STOT SE 3 H335

Additional remarks:

CLP

Regulation (EC) No 1272/2008, Annex VI, Note D

SECTION 4: First aid measures

4.1. Description of first aid measures

General information

No special measures required

After inhalation

Ensure supply of fresh air.

After skin contact

After contact with molten product cool quickly with cold water. Consult a doctor if skin irritation persists.

After eye contact

Separate eyelids, wash the eyes thoroughly with water (15 min.). In case of irritation consult an oculist.

After ingestion

Rinse out mouth and give plenty of water to drink.

Adhere to personal protective measures when giving first aid

First aider: Pay attention to self-protection!

4.2. Most important symptoms and effects, both acute and delayed

Until now no symptoms known so far.

4.3. Indication of any immediate medical attention and special treatment needed

Hints for the physician / hazards

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media

Recommended: alcohol resistant foam, CO₂, powders, water spray/mist, Extinguishing measures to suit surroundings



Trade name: BioBleach soft

Substance number: S0018SOFT

Version: 1 / GB

Date revised: 06.02.2024

Replaces Version: - / GB

Print date: 06.02.2024

5.3. Advice for firefighters

Special protective equipment for fire-fighting

In case of combustion use a suitable breathing apparatus.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Avoid contact with skin, eyes and clothing.

6.2. Environmental precautions

Do not allow to enter drains or waterways.

6.3. Methods and material for containment and cleaning up

Pick up mechanically. Clean contaminated floors and objects thoroughly, observing environmental regulations.

6.4. Reference to other sections

Refer to protective measures listed in Sections 7 and 8.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling

Observe the usual precautions for handling chemicals.

Advice on protection against fire and explosion

No special measures required.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Store product in closed containers.

Hints on storage assembly

Do not store together with foodstuffs.

Further information on storage conditions

Keep container tightly closed and dry.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Other information

Contains no substances with occupational exposure limit values.

8.2. Exposure controls

General protective and hygiene measures

Observe the usual precautions for handling chemicals.

Respiratory protection

Do not inhale vapours



Trade name: BioBleach soft

Substance number: S0018SOFT

Version: 1 / GB

Date revised: 06.02.2024

Replaces Version: - / GB

Print date: 06.02.2024

Hand protection

Heat-resistant gloves when using hot material

Eye protection

Safety glasses

Body protection

Clothing as usual in the chemical industry.

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties****Physical state**

Thermoforming blank

Colour

clear, transparent

Odour

characteristic

Melting point

Value 60 to 110 °C

Freezing point

Remarks not determined

Boiling point or initial boiling point and boiling range

Remarks not determined

Flammability

evaluation not determined

Upper and lower explosive limits

Remarks not determined

Flash point

Remarks Not applicable

Ignition temperature

Remarks not determined

Decomposition temperature

Remarks not determined

pH value

Remarks not determined

Viscosity

Remarks not determined

Solubility(ies)

Remarks not determined

Partition coefficient n-octanol/water (log value)

Remarks not determined

Vapour pressure

Remarks not determined

Density and/or relative densityValue 0,95 g/cm³

Temperature 20 °C

Relative vapour density

Remarks not determined



Trade name: BioBleach soft

Substance number: S0018SOFT

Version: 1 / GB

Date revised: 06.02.2024

Replaces Version: - / GB

Print date: 06.02.2024

9.2. Other information**Odour threshold**

Remarks not determined

Evaporation rate (ether = 1) :

Remarks not determined

Solubility in water

Remarks virtually insoluble

Explosive properties

evaluation not determined

Oxidising properties

Remarks not determined

Other information

None known

SECTION 10: Stability and reactivity**10.1. Reactivity**

No hazardous reactions when stored and handled according to prescribed instructions.

10.2. Chemical stability

No hazardous reactions known.

10.3. Possibility of hazardous reactions

No hazardous reactions known.

10.4. Conditions to avoid

Protect from heat and direct sunlight

10.5. Incompatible materials

None known

SECTION 11: Toxicological information**11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008****Acute oral toxicity**

Remarks Based on available data, the classification criteria are not met.

Acute oral toxicity (Components)**vinyl acetate**Species rat (male)
LD50 3470 mg/kg**Acute dermal toxicity**

Remarks Based on available data, the classification criteria are not met.

Acute dermal toxicity (Components)**vinyl acetate**Species rabbit
LD50 7440 mg/kg



Trade name: BioBleach soft

Substance number: S0018SOFT

Version: 1 / GB

Date revised: 06.02.2024

Replaces Version: - / GB

Print date: 06.02.2024

Acute inhalational toxicity

Remarks Based on available data, the classification criteria are not met.

Acute inhalative toxicity (Components)**vinyl acetate**

Species	rat		
LC50	4000		ppm(V)
Duration of exposure	4	h	
Administration/Form	Vapors		

Skin corrosion/irritation

Remarks Based on available data, the classification criteria are not met.

Serious eye damage/irritation

Remarks Based on available data, the classification criteria are not met.

Sensitization

Remarks Based on available data, the classification criteria are not met.

Subacute, subchronic, chronic toxicity

Remarks not determined

Mutagenicity

Remarks Based on available data, the classification criteria are not met.

Reproductive toxicity

Remarks Based on available data, the classification criteria are not met.

Carcinogenicity

Remarks Based on available data, the classification criteria are not met.

Carcinogenicity (Components)**vinyl acetate**

evaluation Suspected of causing cancer.

Specific Target Organ Toxicity (STOT)**Single exposure**

Remarks Based on available data, the classification criteria are not met.

Repeated exposure

Remarks Based on available data, the classification criteria are not met.

Specific Target Organ Toxicity (STOT) (Components)**vinyl acetate****Single exposure**

evaluation May cause respiratory irritation.

Aspiration hazard

Based on available data, the classification criteria are not met.

11.2 Information on other hazards**Endocrine disrupting properties with respect to humans**

This substance does not have endocrine disrupting properties with respect to humans.

Other information

No toxicological data are available.

SECTION 12: Ecological information

Trade name: BioBleach soft

Substance number: S0018SOFT

Version: 1 / GB

Date revised: 06.02.2024

Replaces Version: - / GB

Print date: 06.02.2024

12.1. Toxicity

General information

not determined

Fish toxicity (Components)

vinyl acetate

Species	Fathead minnow (<i>Pimephales promelas</i>)		
NOEC	0,551		mg/l
Duration of exposure	34	d	
Method	OECD 210		

Daphnia toxicity (Components)

vinyl acetate

Species	Daphnia magna		
EC50	12,6		mg/l
Duration of exposure	48	h	
Method	OECD 202		

Algae toxicity (Components)

vinyl acetate

Species	Pseudokirchneriella subcapitata		
ErC50	12,7		mg/l
Duration of exposure	72	h	
Method	OECD 201		

12.2. Persistence and degradability

General information

not determined

Biodegradability (Components)

vinyl acetate

Value	82	to	98	%
Duration of test evaluation	14	d		
	Readily biodegradable (according to OECD criteria)			

12.3. Bioaccumulative potential

General information

not determined

Partition coefficient n-octanol/water (log value)

Remarks not determined

Octanol/water partition coefficient (log Pow) (Components)

vinyl acetate

log Pow	0,73		
Temperature	25	°C	

12.4. Mobility in soil

General information

not determined

12.5. Results of PBT and vPvB assessment

General information

not determined

Results of PBT and vPvB assessment



Trade name: BioBleach soft

Substance number: S0018SOFT

Version: 1 / GB

Date revised: 06.02.2024

Replaces Version: - / GB

Print date: 06.02.2024

The Substance does not meet PBT-criteria.
This substance does not meet the vPvB-criteria.

12.6 Endocrine disrupting properties**Endocrine disrupting properties with respect to the environment**

This substance does not have endocrine disrupting properties with respect to non-target organisms.

12.7. Other adverse effects**General information**

not determined

General information / ecology

Do not discharge product unmonitored into the environment.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations for the product**

Must not be disposed together with household garbage.
Dispose of waste according to applicable legislation.

Disposal recommendations for packaging

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

SECTION 14: Transport information

	Land transport ADR/RID	Marine transport IMDG/GGVSee	Air transport ICAO/IATA
14.1. UN number or ID number	The product does not constitute a hazardous substance in land transport.	The product does not constitute a hazardous substance in sea transport.	The product does not constitute a hazardous substance in air transport.
14.2. UN proper shipping name	-	-	-
14.3. Transport hazard class(es)		-	-
Label			
14.4. Packing group		-	-

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****Other information**



Trade name: BioBleach soft

Substance number: S0018SOFT

Version: 1 / GB

Date revised: 06.02.2024

Replaces Version: - / GB

Print date: 06.02.2024

All components are contained in the TSCA inventory or exempted.

15.2. Chemical safety assessment

For this preparation a chemical safety assessment has not been carried out.

SECTION 16: Other information

Hazard statements listed in Chapter 2/3

H225	Highly flammable liquid and vapour.
H332	Harmful if inhaled.
H335	May cause respiratory irritation.
H351	Suspected of causing cancer.

CLP categories listed in Chapter 2/3

Acute Tox. 4	Acute toxicity, Category 4
Carc. 2	Carcinogenicity, Category 2
Flam. Liq. 2	Flammable liquid, Category 2
STOT SE 3	Specific target organ toxicity - single exposure, Category 3

Supplemental information

Relevant changes compared with the previous version of the safety data sheet are marked with: ***
This information is based on our present state of knowledge. However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.