



Trade name: GumQuick Comp.A

Substance number: 36601

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 26.09.2025

Print date: 26.09.2025

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

GumQuick Comp.A

This substance/mixture contains components in nanoform

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/preparation

Gingiva mask material

### 1.3. Details of the supplier of the safety data sheet

#### Address/Manufacturer

Dreve Dentamid GmbH

Max-Planck-Straße 31

DE-59423 Unna

Telephone no. +49 2303 8807-0

Fax no. +49 2303 8807-29

Information provided by / telephone Department Research &amp; Development: Fax: +49 2303 8807-562

E-mail address of person responsible for this SDS [sicherheitsdatenblatt@dreve.com](mailto:sicherheitsdatenblatt@dreve.com)

### 1.4. Emergency telephone number

Henkel Fire Department / 24h-Emergency-Contact-No.: +49 211 797-3350

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

This product is not classified hazardous in accordance with Regulation (EC) No 1272/2008.

### 2.2. Label elements

#### Labelling according to regulation (EC) No 1272/2008

#### Supplemental information

EUH210

Safety data sheet available on request.

### 2.3. Other hazards

No special hazards have to be mentioned.

The product contains no PBT substances. The product contains no vPvB substances. This product does not contain a substance that has endocrine disrupting properties with respect to human. The product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms.

## SECTION 3: Composition/information on ingredients

Trade name: GumQuick Comp.A

Substance number: 36601

Version: 2 / GB

Date revised: 26.09.2025

Replaces Version: 1 / GB

Print date: 26.09.2025

### 3.2. Mixtures

#### Chemical characterization

Addition-vulcanising 2-component silicone

#### Hazardous ingredients

##### Cristobalite

CAS No.	14464-46-1				
EINECS no.	238-455-4				
Concentration	>= 25	<	50		%
Classification (Regulation (EC) No. 1272/2008)	STOT RE 1		H372		

##### White mineral oil

CAS No.	8042-47-5				
EINECS no.	232-455-8				
Registration no.	01-2119487078-27				
Concentration	>= 1	<	10		%
Classification (Regulation (EC) No. 1272/2008)	Asp. Tox. 1		H304		

#### Further ingredients

##### Silanameine, 1,1,1-trimethyl-N-(trimethylsilyl)

CAS No.	68909-20-6				
EINECS no.	272-697-1				
Registration no.	01-2119379499-16				
Concentration	>= 10	<	25		%
Advice: [7]					
Classification (Regulation (EC) No. 1272/2008)	STOT RE 2		H373		

Name of set of nanoform	Nanoagglomerate				
Particle size distribution	d50	60-300	nm		
	Method	Transmission Electron Microscopy (TEM)			
Shape and aspect ratio of particles	amorphous				
	Nano particles				
Particle size distribution	d50	6-50	nm		
	Method	Transmission Electron Microscopy (TEM)			
Shape and aspect ratio of particles	amorphous				
40 to 350	m <sup>2</sup> /g				

#### Note

[7] Nanoforms

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

#### General information

No special measures required

#### After inhalation

Ensure supply of fresh air. In the event of symptoms take medical treatment.

Trade name: GumQuick Comp.A

Substance number: 36601

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 26.09.2025

Print date: 26.09.2025

**After skin contact**

In case of contact with skin wash off with warm water. Consult a doctor if skin irritation persists.

**After eye contact**

Separate eyelids, wash the eyes thoroughly with water (15 min.). In case of irritation consult an oculist.

**After ingestion**

Do not induce vomiting - aspiration hazard. Summon a doctor immediately.

**Adhere to personal protective measures when giving first aid**

First aider: Pay attention to self-protection!

**4.2. Most important symptoms and effects, both acute and delayed**

Until now no symptoms known so far.

**4.3. Indication of any immediate medical attention and special treatment needed****Hints for the physician / hazards**

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

## SECTION 5: Firefighting measures

**5.1. Extinguishing media****Suitable extinguishing media**

Recommended: alcohol resistant foam, CO<sub>2</sub>, powders, water spray/mist, Extinguishing measures to suit surroundings

**Non suitable extinguishing media**

Full water jet

**5.2. Special hazards arising from the substance or mixture**

In case of combustion evolution of dangerous gases possible.

**5.3. Advice for firefighters****Special protective equipment for fire-fighting**

In case of combustion use a suitable breathing apparatus.

**Other information**

Fire residues and contaminated fire-fighting water must be disposed of in accordance with the local regulations.

## SECTION 6: Accidental release measures

**6.1. Personal precautions, protective equipment and emergency procedures**

Avoid contact with skin, eyes and clothing.

**6.2. Environmental precautions**

Do not allow to enter drains or waterways.

**6.3. Methods and material for containment and cleaning up**

Pick up mechanically. Clean contaminated floors and objects thoroughly, observing environmental regulations.

**6.4. Reference to other sections**

Trade name: GumQuick Comp.A

Substance number: 36601

Version: 2 / GB

Date revised: 26.09.2025

Replaces Version: 1 / GB

Print date: 26.09.2025

Refer to protective measures listed in Sections 7 and 8.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

#### Advice on safe handling

Observe the usual precautions for handling chemicals. For personal protection see Section 8.

#### Advice on protection against fire and explosion

No special measures required.

### 7.2. Conditions for safe storage, including any incompatibilities

#### Requirements for storage rooms and vessels

Store product in closed containers.

#### Hints on storage assembly

Do not store together with foodstuffs.

#### Further information on storage conditions

Keep container tightly closed and dry.

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

#### Other information

Contains no substances with occupational exposure limit values.

#### Derived No/Minimal Effect Levels (DNEL/DMEL)

##### White mineral oil

Type of value	Derived No Effect Level (DNEL)	
Reference group	Worker	
Duration of exposure	Repeated exposure	
Route of exposure	inhalative	
Mode of action	Systemic effects	
Concentration	164,6	mg/m <sup>3</sup>

Type of value	Derived No Effect Level (DNEL)	
Reference group	Worker	
Duration of exposure	Repeated exposure	
Route of exposure	dermal	
Mode of action	Systemic effects	
Concentration	217,1	mg/kg/d

Type of value	Derived No Effect Level (DNEL)	
Reference group	Consumer	
Duration of exposure	Repeated exposure	
Route of exposure	inhalative	
Mode of action	Systemic effects	
Concentration	34,78	mg/m <sup>3</sup>

Type of value	Derived No Effect Level (DNEL)	
Reference group	Consumer	
Duration of exposure	Repeated exposure	

Trade name: GumQuick Comp.A

Substance number: 36601

Version: 2 / GB

Date revised: 26.09.2025

Replaces Version: 1 / GB

Print date: 26.09.2025

Route of exposure	dermal	
Mode of action	Systemic effects	
Concentration	93,02	mg/kg/d
Type of value	Derived No Effect Level (DNEL)	
Reference group	Consumer	
Duration of exposure	Repeated exposure	
Route of exposure	oral	
Mode of action	Systemic effects	
Concentration	25	mg/kg/d

## 8.2. Exposure controls

### General protective and hygiene measures

Observe the usual precautions for handling chemicals. Do not eat, drink or smoke during work time.  
Wash hands before breaks and after work. Avoid contact with eyes.

### Respiratory protection

Not necessary.

### Hand protection

Not necessary.

### Eye protection

Not necessary.

### Body protection

Clothing as usual in the chemical industry.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

<b>Physical state</b>	liquid, viscous	
<b>Colour</b>	white	
<b>Odour</b>	characteristic	
<b>Melting point</b>		
Remarks	not determined	
<b>Freezing point</b>		
Remarks	not determined	
<b>Boiling point or initial boiling point and boiling range</b>		
Value	> 300	°C
<b>Flammability</b>		
evaluation	not determined	
<b>Upper and lower explosive limits</b>		
Remarks	not determined	
<b>Flash point</b>		
Value	> 130	°C
Method	closed cup	
<b>Auto-ignition temperature</b>		
Remarks	not determined	
<b>Decomposition temperature</b>		

Trade name: GumQuick Comp.A

Substance number: 36601

Version: 2 / GB

Date revised: 26.09.2025

Replaces Version: 1 / GB

Print date: 26.09.2025

Remarks not determined

**pH value**

Remarks not determined

**Viscosity**

Remarks not determined

**Solubility(ies)**

Remarks not determined

**Partition coefficient n-octanol/water (log value)**

Remarks not determined

**Vapour pressure**

Remarks not determined

**Density and/or relative density**

Value	1,23		g/cm <sup>3</sup>
Temperature	20	°C	

**Relative vapour density**

Remarks not determined

**9.2. Other information****Odour threshold**

Remarks not determined

**Evaporation rate (ether = 1) :**

Remarks not determined

**Solubility in water**

Remarks virtually insoluble

**Explosive properties**

evaluation not determined

**Oxidising properties**

Remarks not determined

**Other information**

None known

## SECTION 10: Stability and reactivity

**10.1. Reactivity**

No hazardous reactions when stored and handled according to prescribed instructions.

**10.2. Chemical stability**

No hazardous reactions known.

**10.3. Possibility of hazardous reactions**

No hazardous reactions known.

**10.4. Conditions to avoid**

No hazardous reactions known.

**10.5. Incompatible materials**

None known

**10.6. Hazardous decomposition products**

Trade name: GumQuick Comp.A

Substance number: 36601

Version: 2 / GB

Date revised: 26.09.2025

Replaces Version: 1 / GB

Print date: 26.09.2025

No hazardous decomposition products known.

## SECTION 11: Toxicological information

### 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

#### Acute oral toxicity

Remarks Based on available data, the classification criteria are not met.

#### Acute oral toxicity (Components)

##### White mineral oil

Species	rat	
LD50	> 5000	mg/kg
Method	OECD 401	

#### Acute dermal toxicity

Remarks Based on available data, the classification criteria are not met.

#### Acute dermal toxicity (Components)

##### White mineral oil

Species	rabbit	
LD50	> 2000	mg/kg
Method	OECD 402	

#### Acute inhalational toxicity

Remarks Based on available data, the classification criteria are not met.

#### Acute inhalative toxicity (Components)

##### White mineral oil

Species	rat	
LC50	> 5,2	mg/l
Duration of exposure	4	h
Administration/Form	Dust/Mist	
Method	OECD 403	

#### Skin corrosion/irritation

Remarks Based on available data, the classification criteria are not met.

#### Serious eye damage/irritation

Remarks Based on available data, the classification criteria are not met.

#### Sensitization

Remarks Based on available data, the classification criteria are not met.

#### Subacute, subchronic, chronic toxicity

Remarks not determined

#### Mutagenicity

Remarks Based on available data, the classification criteria are not met.

#### Reproductive toxicity

Remarks Based on available data, the classification criteria are not met.

#### Carcinogenicity

Remarks Based on available data, the classification criteria are not met.

#### Specific Target Organ Toxicity (STOT)

##### Single exposure

Remarks Based on available data, the classification criteria are not met.

Trade name: GumQuick Comp.A

Substance number: 36601

Version: 2 / GB

Date revised: 26.09.2025

Replaces Version: 1 / GB

Print date: 26.09.2025

**Repeated exposure**

Remarks Based on available data, the classification criteria are not met.

**Specific Target Organ Toxicity (STOT) (Components)****Cristobalite****Repeated exposure**evaluation Causes damage to organs through prolonged or repeated exposure  
Route of exposure inhalative**Aspiration hazard**

Based on available data, the classification criteria are not met.

**Aspiration hazard (Components)****White mineral oil**

Harmful: may cause lung damage if swallowed.

**11.2. Information on other hazards****Endocrine disrupting properties with respect to humans**

The product does not contain a substance that has endocrine disrupting properties with respect to humans.

**Other information**

No toxicological data are available.

## SECTION 12: Ecological information

**12.1. Toxicity****General information**

not determined

**Fish toxicity (Components)****White mineral oil**

Species	golden orfe ( <i>Leuciscus idus</i> )			
LC50	>	10000		mg/l
Duration of exposure		96	h	
Method	OECD 203			

**Daphnia toxicity (Components)****White mineral oil**

Species	Daphnia magna			
LL50	>	100		mg/l
Duration of exposure		48	h	
Method	OECD 202			

**White mineral oil**

Species	Daphnia magna			
NOEL	>=	1000		mg/l
Duration of exposure		21	d	
Method	OECD 211			

**Algae toxicity (Components)****White mineral oil**

Species	Pseudokirchneriella subcapitata			
LOEC	>=	100		mg/l
Duration of exposure		72	h	
Method	OECD 201			



Trade name: GumQuick Comp.A

Substance number: 36601

Version: 2 / GB

Date revised: 26.09.2025

Replaces Version: 1 / GB

Print date: 26.09.2025

## 12.2. Persistence and degradability

### General information

not determined

### Biodegradability (Components)

#### White mineral oil

Value	31	%
Duration of test evaluation	28 d	
Method	Moderately/partially biodegradable OECD 301 F	

## 12.3. Bioaccumulative potential

### General information

not determined

### Partition coefficient n-octanol/water (log value)

Remarks not determined

### Octanol/water partition coefficient (log Pow) (Components)

#### White mineral oil

log Pow	4,3	to	18,02
Temperature	20	°C	
Source	ECHA		

## 12.4. Mobility in soil

### General information

not determined

## 12.5. Results of PBT and vPvB assessment

### General information

not determined

### Results of PBT and vPvB assessment

The product contains no PBT substances  
The product contains no vPvB substances.

## 12.6 Endocrine disrupting properties

### Endocrine disrupting properties with respect to the environment

The product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms.

## 12.7. Other adverse effects

### General information

not determined

### General information / ecology

Do not allow to enter soil, waterways or waste water canal. Avoid release into the atmosphere.

## SECTION 13: Disposal considerations

## 13.1. Waste treatment methods

### Disposal recommendations for the product

Must not be disposed together with household garbage.

Trade name: GumQuick Comp.A

Substance number: 36601

Version: 2 / GB

Date revised: 26.09.2025

Replaces Version: 1 / GB

Print date: 26.09.2025

Dispose of waste according to applicable legislation.

**Disposal recommendations for packaging**

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

**SECTION 14: Transport information \*\*\***

	<b>Land transport ADR/RID ***</b>	<b>Marine transport IMDG/GGVSee ***</b>	<b>Air transport ICAO/IATA ***</b>
<b>14.1. UN number or ID number</b>	The product does not constitute a hazardous substance in land transport.	The product does not constitute a hazardous substance in sea transport.	The product does not constitute a hazardous substance in air transport.
<b>14.2. UN proper shipping name</b>	-	-	-
<b>14.3. Transport hazard class(es)</b>		-	-
Label			
<b>14.4. Packing group</b>		-	-

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****Other information**

All components are contained in the TSCA inventory or exempted.

**15.2. Chemical safety assessment**

For this preparation a chemical safety assessment has not been carried out.

**SECTION 16: Other information****Hazard statements listed in Chapter 2/3**

H304 May be fatal if swallowed and enters airways.  
H372 Causes damage to organs through prolonged or repeated exposure.

**CLP categories listed in Chapter 2/3**

Asp. Tox. 1 Aspiration hazard, Category 1  
STOT RE 1 Specific target organ toxicity - repeated exposure, Category 1

**Supplemental information**

Relevant changes compared with the previous version of the safety data sheet are marked with: \*\*\*  
This information is based on our present state of knowledge. However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.